

## ZooMe Auto-Flex Folding Travel Scooter



HCPCS:



FIGURE A



FIGURE B



FIGURE C



FIGURE D




FIGURE E



FIGURE F




FIGURE G



FIGURE H

- Fully electronic folding and unfolding in 15 seconds or less with a wireless key fob
- Aluminum alloy frame is lightweight for superior portability
- Airline safe lithium battery
- High-gloss shrouds for a protective finish
- Easy throttle control
- Delta tiller for users with limited dexterity
- Auto-folding seat with padded, adjustable armrests
- Flat-free, non-marking tires for worry-free travel
- Easy-to-adjust tiller places the scooter's controls in the ideal position
- Secondary manual-release levers and toggle switch for folding and unfolding
- Front and rear reflectors
- 2 ignition keys and 1 key fob included
- Anti-tip wheels for additional safety

### PRODUCTS

### Description

FLEX-AUTO	ZooMe Auto-Flex Folding Travel Scooter, EA
-----------	--

### SPECIFICATIONS

### Description

Anti-Tip Wheels	2" x .5"
Base Weight	55 lbs.
Batteries (pair)	25.9V x 10AH Lithium
Battery Charger	2A Offboard
Battery Weight (pair)	5 lbs.
Brakes	Electromagnetic
Climbing Angle	6°
Controller	PG S Drive 45A
Estimated Range	13 miles
Freewheel Mode	Yes
Front Wheels	6" x 1.5"
Ground Clearance	1.6"
Heaviest Piece Weight	55 lbs.
Motor	24V x 120W
Overall Dimensions	37" (L) x 17" (W) x 32.5" (H)
Overall Dimensions Folded	17" (L) x 17" (W) x 23.5" (H)
Rear Wheels	7" x 2.35"
Seat Dimensions	16" (W) x 14" (D)

**SPECIFICATIONS****Description**

Top Speed	4 mph
Total Weight (excl. hand control)	60 lbs.
Turning Radius	47.2"
Warranty Motor   Electronics   Battery	24   24   12 Months
Warranty on Frame	Limited Lifetime
Weight Capacity	300 lbs
Weight Capacity (13046   13046-HD)	300 lbs.