



EW-M49 Folding Smart Wheelchair

The EW-M49 Smart Folding Power Wheelchair is one of the most advanced fully adjustable folding wheelchairs on the market today! Equipped with both air and spring suspension this wheelchair is the ultimate in comfort. Additionally the fully adjustable backrest and neck rest provide a custom fitted experience. The digital Smart control panel provides the rider with all needed information during operation. This top of the line EWheels Medical wheelchair is an absolute market leader.

Key Specs:

- Digital Smart Control Panel
- Air Spring Backrest Suspension
- Adjustable Backrest Maximum 135 degree
- Ergonomic Armrest Adjustable Height, Forward/Backward movement
- Front Fork suspension
- Adjustable Neck Rest Height and Forward/Back Rotation

General Specifications:

- Swing-Away joystick: Yes
- Folded Size: 18.9"L x 24.4"W x 34.6"H
- Unfolded Size: 40"L x 24.4"W x 40"H
- Lap Belt: Yes
- Storage Area: Yes, Conveniently Located
- Motor: Two 250W Brushless Hub Motors
- Brake: Electromagnetic
- Back Seat Pocket: Yes

Front Suspension: Yes
Optional Left or Right Hand Joystick: Yes
Frame Material: Aluminum
Top Speed: 3.7 mph
Weight Capacity: 330 lbs.
Heaviest Piece: 65 lbs.
Battery Charge Distance: 9.7 mph
Climbing Slope: 12 Degree
Rear Tire Splash Guard: Yes
Drive Type: Rear Wheel Drive
Anti-Tipping Rear Wheels: Yes
Foldable: Yes
Footplate: Large Aluminum
Seat Width: 18"
Seat Depth: 17"
Seat to Floor Height: 20.5"
Backrest Height: 21.6"
Overall Width: 24.4"
Overall Length: 40"
Joystick Options: Right or Left
Frame Type: Aluminum Transportable
Battery: 24V 12AH Lithium Battery
Battery Included: Yes
Battery Charger Type: Off Board
Turning Radius: 32"
Rear Wheel Size: 12"
Front Wheel Size: 8"
Flat Free Tires: Yes
Battery Pack Weight: 5 lbs.
Weight Without Batteries: 65 lbs.