

	<i>EW-10</i>
<i>Warranty</i>	3 Year Limited Warranty
<i>Colors</i>	Red, Blue
<i>Weight Capacity</i>	400 lbs
<i>Fully Assembled</i>	Yes
<i>3 or 4 Wheel</i>	3 Wheels
<i>FWD or RWD</i>	Rear Wheel Drive
<i>Max Passengers</i>	2
<i>Airline Approved</i>	No
<i>Seat Dimensions (LxWxH)</i>	14" x 26"
<i>Type of Seat</i>	Extended Saddle Seat
<i>Type of Battery</i>	Deep Cycle Sealed Lead Acid
<i>Overall length</i>	70"
<i>Overall width</i>	28"
<i>Overall height</i>	43"
<i>Armrests</i>	No
<i>Scooter Weight</i>	209 lbs
<i>Shipping Method</i>	LTL Freight
<i>Top speed (Up To)</i>	15 mph
<i>Distance (Up To)</i>	40 miles
<i>Light(s)</i>	Turn signals, brake lights, high/low beam, headlight rear light
<i>Turning Radius</i>	66"
<i>Braking Type</i>	Hand Brakes
<i>Left Hand Control Option</i>	Yes
<i>Horn</i>	Yes
<i>Incline capability</i>	Up to 12 Degrees
<i>Ground Clearance</i>	6.5"
<i>Suspension</i>	Front and Rear
<i>Type of Tires</i>	Tubeless
<i>Battery Charger</i>	48 volt offboard
<i>Speed Control</i>	Low-Med-High
<i>Basket</i>	No
<i>Forward/Reverse</i>	Yes
<i>Throttle Type</i>	Twist
<i>Alarm</i>	Yes
<i>Front Tire Size</i>	15" x 3"
<i>Rear Tire Size</i>	15" x 3"
<i>Cell Phone Holder</i>	No
<i>USB Port</i>	No
<i>Voltage</i>	48 volts
<i>Motor Wattage</i>	500 watt
<i>Motor Type</i>	Brushless Motor with Differential
<i>Storage</i>	Lockable Rear and Under Seat
<i>Rear View Mirror</i>	Yes
<i>Extras</i>	Bag hook, rear reflectors,