



EW-M39

Features include a maximum Speed of 5mph, a 300 lb. capacity, a range of 14.5 miles (18AH), wireless design for trouble free dismantling, an on board battery charger, and delta handlebar.

EW-M39 Specifications

Weight Capacity:
300 lbs

Range:
15 miles

Front LED Light

Overall Length
40 "

Overall Width
20"

Overall Height
34"

Wheels:

Front 8"

Rear 8"

Easy Remove, Quick Connect Battery Pack

Weight w/ Batteries

93 lbs

Disassembles Into 5 Pieces:

Heaviest Piece 30 lbs

Weight Capacity

300 lbs

Ground Clearance

1.6"

Grade Climbable

8 degree

Turning Radius

44 "

Brake

Electromagnetic

Seat Type

Swivel Padded Foldable

Seat Width

17"

Motor Size

200 W

Battery Size

Powerful 18AH

Battery Weight
20lbs

Electronics
On/Off Key Switch, Battery Indicator, Speed Control Knob