

WEIGHS ONLY 19 lbs.!

**NEW
Lighter
Weight!**



EXP19LTRD



**HANDY LEVER ON CHAIR
BACK MAKES IT EASY
AND CONVENIENT TO
FOLD BACK DOWN!**



**THE EXPEDITION FOLDS
FLAT FOR EASY TRANSPORT!**

Lightweight Expedition Aluminum Transport Chair with Loop Locks and 12" Rear "Flat Free" Wheels


EXP19LTBL Blue, 1/cs

EXP19LTRD Red, 1/cs

FEATURES:

- Weighs only 19 lbs.
- Available in blue or red frames.
- Back folds down for storage and transport.
- Comes with seat belt.
- Companion activated wheel locks.
- 12" rear "flat free" wheels.
- Limited Lifetime Warranty.

Available
1 Pack



Blue

Red

SPECIFICATIONS

Item #	Overall Width (open)	Overall Width (closed)	Overall Height	Seat To-Floor Height	Back Height	Width Seat	Depth Seat	Overall Length w/Riggings	Product Weight	Carton Shipping Weight	Weight Limit
EXP19LTBL, RD	22"	10.75"	38"	19"	18"	19"	16"	38"	19 lbs.	24 lbs.	300 lbs.

transport chairs

To order call toll free: 877.224.0946

drive