

"I sell several [scooter lines], but the scooter that I have used for the past twelve years is a PaceSaver".
Bernie Webb, Sec./Tres., Mobility Plus, Inc. (dealer)

Specifications	Espree Premier 9"/10"	Espree Titan	Espree Atlas
Overall Length	45.5"	45.5"	45.5"
Overall Width - 9"	24"	N/A	N/A
Overall Width - 10"	24" or 24.75"	24.75"	24.75"
Wheel Base	33.5"	33.5"	33.5"
Weight			
Rear Drive Unit	40/41lbs.	41 lbs.	41 lbs.
Base & Tiller	36.5"/37 lbs.	37 lbs.	45 lbs.
Seat Section	24.5 lbs.	24.5 lbs.	29.75 lbs.
Shroud	3.50 lbs.	3.50 lbs.	3.50 lbs.
Maximum Speed	4.30/5.25 MPH	5.40 MPH	5.40 MPH
Range (approximate)	20 - 30 miles	20 - 30 miles	20 - 30 miles
Ground Clearance	4.25"/4.50"	4.50"	4.50"
Turning Radius	34.5"/35"	35"	35"
Pneumatic Tire Size - 9" (front/back)	2.5" x 9"/3" x 9"	N/A	N/A
Pneumatic Tire Size - 10" (front/back)	3" x 9"/3" x 10"	3" x 9"/3" x 10"	3" x 9"/3" x 10"
Seat Dimensions			
Width	18"	18"	21"
Depth	18"	18"	21"
Back Rest Height	20"	20"	22.5"
Seat Height From Floorpan	16.75" - 18.75"	16.75" - 18.75"	16.75" - 18.75"
Seat Height With Power Seat	18" - 23"	18" - 23"	18" - 23"
Electronic Brakes			
Regenerative, Dynamic	Yes	Yes	Yes
Disc Braking	"Posi-Lock" Electric	"Posi-Lock" Electric	"Posi-Lock" Electric
Rear Wheel Drive	Transaxle	HD Transaxle	HD Transaxle
Power Source	1 HP	1 1/4 HP	1 1/4 HP
Structural Capacity	300 lbs	350 lbs	450 lbs

PaceSaver Extra Degree of Protection™ -			
Maximum Incline Rating	12° (22% grade)	12° (22% grade)	12° (22% grade)

Specifications subject to change. Incline stability is determined per US Ansi/Resna Dynamic Stability test 8.2. Stability increases by leaning forward or seat not in highest position.

Five-Year Limited Warranty
Five-year prorated limited transaxle warranty. Three-year prorated limited warranty on Main Frame components. One-year prorated limited warranty covering virtually all internal components including the electronics.

Lifetime Extended Warranty
Lifetime warranty on frame components. Two-year prorated limited parts warranty. Please consult the Owners Manual for eligibility requirements and details of the warranties. Batteries are not part of either warranty.

Battery Charger - Battery chargers are standard equipment with PaceSavers and are a heavy duty, Micro-Light 24 volt constant current fully automatic charger specially manufactured for Leisure-Lift.

Batteries - PaceSavers are designed to utilize size U-1 sealed deep cycle type batteries. These batteries, supplied by your dealer, may be purchased from Leisure-Lift or supplied by a local source. Leisure-Lift only provides warranty coverage on its own batteries. Because they are 24 volts, PaceSavers require two (2) batteries. Please note: all batteries require some form of maintenance. Deep cycle batteries must be kept charged! Batteries left for some time without using and/or charging will go bad.

Your Professional Local PaceSaver Dealer

Color Choices

- Ocean Blue
- Midnight Black
- Candy Apple Red

Velvet Seats

- Black
- Crimson
- Gray

Plush Vinyl Seats

- Charcoal
- Tan
- Gray

Accessories & Options

- U1 Sealed Batteries
- Walker Holder
- Sliding Seat
- Wide Seat
- Admiral Seat
- Flat-Free Tires
- Pelvic Strap
- Stretch Frame



Please see the Accessories & Options brochure for a complete listing.

PaceSaver® Scooters



- Top Reliability*
- Rock-Solid Stability*
- Beautiful Styling*
- Real Transportability*
- 5 - Year Warranty*
- Maximum Comfort*
- Long-Term Value*

Our most advanced scooters

PaceSaver® Scooters! Celebrating 21 years of legendary performance.

PaceSaver Scooters!
Celebrating 21 years of legendary performance

Set Your Own Pace!

PaceSaver®

OUR HISTORY

In 1967, Bill and Zella Burke risked all they had to invent and manufacture healthcare products. Mr. Burke said “[I will] dedicate myself to inventing home care products for the infirm”. The resulting family of products allows people to maintain their independence so they can stay in their home close to friends and relatives.

The company that Mr. and Mrs. Burke started employed thirty people; twenty were seniors or disabled and they received national awards for hiring veterans. In our 3rd generation, we continue Mr. Burke’s dedication to building reliable long lasting products and supporting America.

WHAT AMERICAN MADE MEANS TO YOU

Today while many scooters and power chairs are made in low production cost countries like China, PaceSaver continues the 30 year tradition of building to a higher standard of durability and reliability. We are the **only** manufacturer of scooters and power chairs to build all our models in the U.S.A. Because our products are built in the U.S.A., the FDA inspects our plant for design controls, good manufacturing practices and quality systems on a regular basis while foreign plants are not. Industry tests of foreign products reveal poor quality components and designs that do not work as claimed or worse are dangerous if used per the instructions in the owners manual. American made gives you a higher level of Quality Assurance.


AMERICAN MADE MEANS U.S. STOCKED PARTS

U. S. stocked parts keeps your scooter running for years. We stock millions of dollars worth of parts in Kansas City, whether for new product or decades old. Our parts are “over here” not “over seas”! So over 90% of our replacement parts are shipped within one working day.

AMERICAN MADE KEEPS AMERICANS WORKING

Leisure-Lift keeps all of America working! We buy over 5,000 parts manufactured by 580 vendors located in 44 different states. The efforts of people from all over the United States go into our mobility products.

The New Com-For-Back™ Bucket Seat for all day support & comfort.



All PaceSaver products are built in America.

Espree™

The spirit of freedom planted in tradition.



Bill Burke, PaceSaver



Extra Degree of Protection™

- Incline stability is determined per US Ansi/Resna Dynamic Stability Test 8.2.
- A confident ride - indoors & out.

Mirror Glaze™ Lasting Beauty

- Rich “depth of color” reflective finish.
- Color “through and through” ABS substrate.
- Won’t chip, peel or fade.
- “Scratch guard” easy maintenance.
- Deep clear coat hardened acrylic layered.
- Simply wipe with damp rag or wax.
- No paint, No ugly touch ups.

Choices

Our most luxurious and advanced flagship scooter, all Espree’s feature the Soft-touch tiller for **shoulder comfort** and **easy access**, luxurious “Com-For-Back” high backed seat for all day comfort and support, stylish mag wheels, lights, basket, “Mirror-Glaze” color for **long lasting beauty**. Espree performance features the **industry’s best incline stability** rating, the **most powerful, easiest to drive, safer state-of-the-art standard controller**, quiet powerful 3 brake transaxle drive system, **shock absorbing resilient standard suspension**, **easy no tool separation** for transport, **on board worldwide charger** and the industry’s longest warranty. Enjoy a comfortable, confident ride on America’s most luxurious scooter -- *Espree*.

Espree Premier

Drive with elegant luxurious full size comfort and confidence. *Espree* features excellent indoor maneuverability and abundant leg room, 300 lb. capacity, brisk 5.25 mph. speed, 20 - 30 mile range and the industry’s best stability and warranty.

Espree Titan

The brawny *Espree Titan*, our maximum power Espree design, for long hill climbs and demanding terrain. Titan adds a rugged 4 Pole 4” motor for heavy duty efficiency and 25% more power with swift 5.4 mph. speed and robust 350 lb. capacity. Espree Titan is the American made luxurious performance scooters.

Espree Atlas

Our largest capacity Espree at 450 lb., *Atlas* features a heavy duty reinforced frame with side anti-tips for stability, plush 21” wide seat, height, width and length adjustable arm-rests and a 450 lb. power seat option. Espree Atlas’ distinctive style, maneuverability and transportability with luxurious comfort.

“I just love my new PaceSaver [Espree] Atlas. It drives like a dream and is so good looking. Oregon weather being what it is, it will get much use”. “Thanks”.

Wendy Turner, Portland, OR.



Digital Electronics

- Powerful 110 amp. controller for more power.
- *Push-Too-Fast™* eliminates disengage runaways.
- *Power Boost™* technology for inclines.
- *Battery Saver™* for extended range.
- Audio-visual indicators provide constant status.



American made controller is enclosed & isolated. VA & CE moisture tested.

- *Look under the hood!*
- Encased digital controller.
- Isolated from road moisture.
- Double layered moisture protected throttle.

Smooth, Reassuring Control

- Advance Controller Programs for smooth transition through all speed ranges.
- Meets FDA required radio wave runaway protection of 20 volts per meter.

Real Transportability

- Quick & easy single pull-pin separation.
- Separates into 4 easy pieces.
- Fits in most cars, vans and SUV’s.
- Or use our Magnum Power Lift.



Easily Transportable