



ORDER FORM

Everyday Activity Seat

Activity chair for children or young adults with moderate needs who still require pelvic stability during the day.

Customer			
Account			
Account #		Date	
<input type="checkbox"/> Quote	<input type="checkbox"/> Order	PO #	
Buyer			
Mark For			

Shipping Address			
Name			
Address			
City		State	
ZIP		Phone	
Email / Fax			

The HCPCS CODES herein are based on PDAC verification or interpretation of Medicare definitions and guidelines. Non-Medicare payers may accept alternative HCPCS CODES, including misc. codes to ensure access for their enrollees. The use of HCPCS CODES does not ensure coverage or payment.

Everyday Activity Seat

Size 1 Weight Capacity: 66 lbs | Size 2 Weight Capacity: 110 lbs | Size 3 Weight Capacity: 220 lbs

✓ OPTION	DESCRIPTION	HCPCS	RETAIL
Base Frame			\$
Selection Required.			
<input type="checkbox"/> 151-6102	Size 1 Seat Base, Backrest and Hi-Low Chassis with Gas Spring	E1399	\$3,118
<input type="checkbox"/> 151-6202	Size 2 Seat Base, Backrest and Hi-Low Chassis with Gas Spring	E1399	\$3,236
<input type="checkbox"/> 151-6302	Size 3 Seat Base, Backrest and Hi-Low Hydraulic Assist Chassis	E1399	\$5,998
Chair Cover			\$
Selection Required. Must match base frame size.			
<input type="checkbox"/> 151-6104-03VIN	Size 1 Vinyl Seat and Back Covers - Blue		\$374
<input type="checkbox"/> 151-6206-03VIN	Size 2 Vinyl Seat and Back Covers - Blue		\$423
<input type="checkbox"/> 151-6304-03VIN	Size 3 Vinyl Seat and Back Covers - Blue		\$702
Lateral Support			\$
Must match base frame size.			
<input type="checkbox"/> 151-6106-07VIN	Size 1 Rigid Chest Laterals - Black Vinyl (Pair)		\$392
<input type="checkbox"/> 151-6210-07VIN	Size 2 Rigid Chest Laterals - Black Vinyl (Pair)		\$392
<input type="checkbox"/> 151-6306-07VIN	Size 3 Rigid Chest Laterals - Black Vinyl (Pair)		\$535
<input type="checkbox"/> 151-6108-07VIN	Size 1 Swing-Away Chest Laterals - Black Vinyl (Pair)		\$488
<input type="checkbox"/> 151-6212-07VIN	Size 2 Swing-Away Chest Laterals - Black Vinyl (Pair)		\$493
<input type="checkbox"/> 151-6308-07VIN	Size 3 Swing-Away Chest Laterals - Black Vinyl (Pair)		\$647
<input type="checkbox"/> 151-6218-07VIN	Size 1 and 2 Seat Protraction Pads - Black Vinyl (Pair)		\$475
<input type="checkbox"/> 151-6314-07VIN	Size 3 Seat Protraction Pads - Black Vinyl (Pair)		\$512
<input type="checkbox"/> 137-808-05	Chest Support - Size 1 & 2 (Requires rigid or swing-away laterals)		\$187
<input type="checkbox"/> 133-625	Chest Support - Size 3 (Requires rigid or swing-away laterals)		\$187
Harness			\$
<input type="checkbox"/> 155-N001-L600	Butterfly Harness (Small)		\$406
<input type="checkbox"/> 155-N001-L601	Butterfly Harness (Medium)		\$406
<input type="checkbox"/> 155-N001-L602	Butterfly Harness (Large)		\$406
Upper Extremity Positioning			\$
Must match base frame size. Tray requires armrests.			
<input type="checkbox"/> 151-6117	Size 1 Armrests (Pair)		\$407
<input type="checkbox"/> 151-6226	Size 2 Armrests (Pair)		\$407
<input type="checkbox"/> 151-3747	Size 3 Armrests (Pair)		\$573
<input type="checkbox"/> 151-1644	Size 1 Tray - Plastic Gray		\$407
<input type="checkbox"/> 151-2644	Size 2 Tray - Plastic Gray		\$407
<input type="checkbox"/> 151-3767	Size 3 Tray - Plastic Black		\$510

Everyday Activity Seat

✓ OPTION		DESCRIPTION	RETAIL
Head Supports			\$
<input type="checkbox"/>	151-6208-07VIN	Size 1 and 2 Flat Headrest Support and Cushion - Black Vinyl Must match base frame size	\$381
<input type="checkbox"/>	151-6209-07VIN	Size 1 and 2 Laterals for Flat Headrest - Black Vinyl (Pair) Must match base frame size. Includes width adjustable lateral pads & covers. Requires Flat Headrest Support and Cushion.	\$321
<input type="checkbox"/>	151-3751-05	Size 3 Contoured Headrest - Black Neoprene Must match base frame size	\$495
<input type="checkbox"/>	151-1756	Size 1 and 2 Contoured Headrest Adapter Bracket	\$85
<input type="checkbox"/>	120-859	Size 1 and 2 Contoured Headrest Support Hardware	\$385
<input type="checkbox"/>	137-800-05	Size 1 and 2 Contoured Headrest Cushion - Black	\$176
Foot Positioning			\$
<input type="checkbox"/>	151-1745	Size 1 One Piece Footplate - Metal	\$407
<input type="checkbox"/>	151-2745	Size 2 One Piece Footplate - Metal	\$407
<input type="checkbox"/>	151-3745	Size 3 One Piece Footplate - Metal	\$692
<input type="checkbox"/>	152-1600	Small Sandals w/ Straps (Pair)	\$273
<input type="checkbox"/>	152-2600	Medium Sandals w/ Straps (Pair)	\$298
<input type="checkbox"/>	152-3600	Large Sandals w/ Straps (Pair)	\$321
<input type="checkbox"/>	151-1760	Size 1 and 2 Sandal Raiser 1" (Each) Requires Sandal Selection, must match Sandal size.	\$60
<input type="checkbox"/>	151-3760	Size 3 Sandal Raiser 1" (Each) Requires Sandal Selection, must match Sandal size.	\$69
	QTY	<input type="text"/>	
	QTY	<input type="text"/>	
Pelvic Fixation			\$
<input type="checkbox"/>	151-1660	Size 1 Pelvic Belt 4-Point	\$149
<input type="checkbox"/>	151-2660	Size 2 Pelvic Belt 4-Point	\$160
<input type="checkbox"/>	151-3660	Size 3 Pelvic Belt 4-Point	\$172
<input type="checkbox"/>	151-6115-07VIN	Size 1 Side Pads - Black Vinyl (Pair)	\$529
<input type="checkbox"/>	151-6222-07VIN	Size 2 Side Pads - Black Vinyl (Pair)	\$629
<input type="checkbox"/>	151-6317-07VIN	Size 3 Side Pads - Black Vinyl (Pair)	\$808
<input type="checkbox"/>	151-6224-07VIN	Size 2 Adduction Side Pad - Black Vinyl	\$303
<input type="checkbox"/>	151-6320-07VIN	Size 3 Adduction Side Pad - Black Vinyl	\$303
<input type="checkbox"/>	151-6120-07VIN	Size 1 Wedge Shape Pommel - Black Vinyl	\$214
<input type="checkbox"/>	151-6229-07VIN	Size 2 Wedge Shape Pommel - Black Vinyl	\$256
<input type="checkbox"/>	151-6334-07VIN	Size 3 Wedge Shape Pommel - Black Vinyl	\$266

Technical Specifications - Imperial

Frame Sizing Specification	Size 2	Size 3
Max Load	110 lbs	220 lbs
Weight [Seat and Hi-Low Chassis]	75 lbs	108 lbs
Seat Height [Hi-Low Chassis]	12.6" - 23.6"	17.7" - 23.6"
Seat Width	10.2" - 15"	13.4" - 19.7"
Seat Depth	9" - 15"	13" - 18.9"
Lateral Distance [Chest]	6.7" - 9.4"	9.4" - 14.6"
Back Height	15.7" - 20.5"	17.7" - 24.8"
Back Angle	30° Supine / 10° Prone	30° Supine / 10° Prone
Armrest Height	7.9" - 11.8"	7" - 12.7"
Lower Leg Length	9.8" - 17.7"	9.8" - 17.7"
Lower Leg Angle	-10° to +50°	-10° to +50°
Footplate Angle	-15° to +10°	-15° to +10°
Tilt-In-Space	30° Supine / 10° Prone	30° Supine / 10° Prone

Technical Specifications - Metric

Frame Sizing Specification	Size 2	Size 3
Max Load	50 kgs	100 kgs
Weight [Seat and Hi-Low Chassis]	34 kgs	49 kgs
Seat Height [Hi-Low Chassis]	320 mm - 600 mm	450 mm - 600 mm
Seat Width	260 mm - 380 mm	340 mm - 500 mm
Seat Depth	230 mm - 380 mm	330 mm - 480 mm
Lateral Distance [Chest]	170 mm - 240 mm	240 mm - 370 mm
Back Height	400 mm - 520 mm	450 mm - 630 mm
Back Angle	30° Supine / 10° Prone	30° Supine / 10° Prone
Armrest Height	200 mm - 300 mm	180 mm - 325 mm
Lower Leg Length	250 mm - 450 mm	250 mm - 450 mm
Lower Leg Angle	-10° to +50°	-10° to +50°
Footplate Angle	-15° to +10°	-15° to +10°
Tilt-In-Space	30° Supine / 10° Prone	30° Supine / 10° Prone

Sunrise Medical (US) LLC

2842 Business Park Avenue
 Fresno, CA 93727 USA
 Tel.: 1-800-333-4000
 Fax: 1-800-300-7502
 E-Mail: webmaster@sunmed.com



If you are visually impaired, this document can be viewed in PDF format at www.SunriseMedical.com



MK-129846 Rev. F ©01/2024



800.333.4000

www.SunriseMedical.com