

Bill To (Ultimate Medical Supplier Only)

Purchase Order #: _____

Quote #: _____

Account #: _____

Ordered By: _____

Supplier Company: _____

Address: _____

City, State, Zip: _____

Phone: _____ Email: _____

Payer Source: _____

Name of ATP (Required): _____

Ship To

Name: _____

Address: _____

City, State, Zip: _____

Phone: _____

Tag For: _____

Evolv XT User Specifications

Approximate Height Range 6'0"-6'10"

Weight Limit 350 lbs.

Seat Depth Range (seat pivot to flat back) 20"-25"

Seat to Foot Plate Range (from seat pivot) 17"-23"

Base *Note: The DEFAULT items will be added if no selection is made.*

■ **EasyStand[®] Evolv[®] XT** \$3,631

Price Includes:

PNG50209-1 EasyStand Evolv XT Base \$2,259

PNG50377 Extra Tall Configuration \$1,372

Lift Mechanisms-must select one

☐ PNG50320 Manual Hydraulic Actuator **DEFAULT** N/C

Manual Hydraulic Actuator Handle Location

☐ PNG50059 Right Side **DEFAULT** N/C

☐ PNG50126 Left Side N/C

Manual Hydraulic Actuator Handle Accessories

☐ PNG50328 Adjustable/Removable Actuator Handle \$236

☐ PNG50044 Quad Grip Handle Extension \$159

☐ PNG50070 T-Style Handle Extension \$159

☐ PNG50321 Pow'r Up Lift (battery powered) \$1,265

Front Frames and Trays**MUST SELECT ONE FRAME AND ONE TRAY IN THE SAME FRAME SECTION****Shadow Tray Frame-back required***

☐ PNG51054 Column Front for Shadow Tray with Wheels \$382

☐ PNG51056 Column Front for Shadow Tray with Swivel Casters \$568

Shadow Tray*

☐ PNG50515 Black Molded Adjustable Shadow Tray-Right \$1,465
(7"-12"H from seat, 4"-20"D from back)

☐ PNG50537 Black Molded Adjustable Shadow Tray-Left \$1,465
(7"-12"H from seat, 4"-20"D from back)

☐ PNG50517 Clear Adjustable Shadow Tray-Right \$1,635
(6"-11"H from seat, 4"-20"D from back)

☐ PNG50538 Clear Adjustable Shadow Tray-Left \$1,635
(6"-11"H from seat, 4"-20"D from back)

Column Frame*

☐ PNG51055 Column with Wheels **DEFAULT** \$382

☐ PNG51057 Column with Swivel Casters \$568

Column Tray*

☐ PNG50247-1 Black Molded Tray **DEFAULT** \$606

☐ PNG50248-1 Clear Tray \$780

☐ PNG50194-1 Oversized Clear Tray \$924

☐ PNG50368-1 Black Molded Angle Adjustable Tray \$967

☐ PNG50366-1 Clear Angle Adjustable Tray \$1,134

Mobile Frame

☐ PNG51058 Column with Mobile Option \$1,903

Mobile Tray

☐ PNG50372-1 Narrow Clear Tray \$780

☐ PNG50327-1 Narrow Clear Angle Adjustable Tray \$1,134

Tray Accessories

☐ PNG50042 Hand Grips (for clear tray only, pair)* \$164

☐ PNG30295 Padded Tray Cover (not available with oversized clear trays or mobile trays)* \$234

Chest Pad-must select one

☐ PNG30001 Medium Chest Pad **DEFAULT** N/C

☐ PNG30000 Large Contoured Chest Pad \$147

☐ PNG50331 Large Contoured Adjustable Chest Pad \$212

Knee Pads-must select one

☐ PNG50460-1 1-Piece Knee Pad (pair) **DEFAULT** \$340

☐ PNG50384-1 Independent Knee Pads (pair) \$843

☐ PNG50385-1 Independent Roho[®] Knee Pads (pair) \$969

* not available with mobile frame

Foot Plates-must select one

- ☐ PNG50278-1 Multi-Adjustable Foot Plates (pair) **DEFAULT** \$342

Foot Plate Accessories

- ☐ PNG30030 Foot Straps-15"L (length over top of foot, pair) \$50
- ☐ PNG50505 Foot Straps-21"L (length over top of foot, pair) \$50
- ☐ PNG30031 Secure Foot Straps-15"L (length over top of foot, two pair) \$126
- ☐ PNG50506 Secure Foot Straps-21"L (length over top of foot, two pair) \$126

Seats-must select one

- ☐ PNG30006-1 Planar Seat **DEFAULT** \$258
- ☐ PNG30270 Transfer Seat (replaces standard planar seat) \$402
- ☐ PNG50318 Rotating Seat (replaces standard planar seat) \$508
- ☐ PNG30023 Seat with Roho[®] Insert (replaces standard planar seat) \$631

Seat Accessories

- ☐ PNG50500 Easy-Adjust Seat Depth-CANNOT be ordered off-unit \$112
(roller bearing allows for smoother seat depth adjustment)

Pelvic Support

- ☐ PNG50161 Hip Supports (10"-15" range, pair) \$376
- ☐ PNG50198 Hip Supports (12"-17" range, pair, limited adjustment with shadow tray) \$376
- ☐ PNG50518 Hip Supports (17"-22" range, pair, limited adjustment with shadow tray) \$376
- ☐ PNG30029 Velcro[®] Positioning Belt-26"-55" (back required) \$65
- ☐ PNG30306 Velcro[®] Positioning Belt-38"-67" (back required) \$65
- ☐ PNG30028 Positioning Belt with Airline Style Buckle-26"-55" (back required) \$159

Backs

- ☐ PNG50066 Flat Back-15"H \$689
- ☐ PNG50067 Flat Back-19"H \$727
- ☐ PNG50062 Contoured Back-15"H \$758
- ☐ PNG50068 Contoured Back-19"H \$796
- ☐ PNG50171 Removable Flat Back-15"H \$790
- ☐ PNG50173 Removable Flat Back-19"H \$828
- ☐ PNG50172 Removable Contoured Back-15"H \$860
- ☐ PNG50174 Removable Contoured Back-19"H \$899

Back Accessories

- ☐ P600002 Back Angle Adjustment Wheel (back required) \$108
- ☐ P80246 Angle Locator \$40

Notes:**Chest Straps-back required**

- ☐ PNG50192 Chest Strap-30"-46" \$159
- ☐ PNG50193 Chest Strap-34"-54" \$159
- ☐ PNG50504 Chest Strap-46"-66" \$159

Upper Body Support-back required**ATTENTION - ALL OPTIONS BELOW REQUIRE THE ACCESSORIES MOUNTING BRACKET**

- ☐ PNG50034 Accessories Mounting Bracket \$71

Chest Vests-back and accessories mounting bracket required

- ☐ PNG11212 High Mount Chest Vest Bracket \$81
- ☐ PNG30170 X-Style Chest Vest-11"Lx9.5"W \$297
- ☐ PNG30056 X-Style Chest Vest-14.5"Lx10.5"W \$297
- ☐ PNG30171 Y-Style Chest Vest-17.5"Lx11.5"W \$297
- ☐ PNG30021 Y-Style Chest Vest-20.5"Lx12.5"W \$297
- ☐ PNG30318 Non-Stretch X-Style Chest Vest-11"Lx9.5"W \$345
- ☐ PNG30322 Non-Stretch X-Style Chest Vest-14.5"Lx10.5"W \$345
- ☐ PNG30320 Non-Stretch Y-Style Chest Vest-17.5"Lx11.5"W \$345
- ☐ PNG30324 Non-Stretch Y-Style Chest Vest-20.5"Lx12.5"W \$345

Lateral Supports-back and accessories mounting bracket required

- ☐ PNG50535 Lateral Supports with Flat Pads-8"-16.5"W (pair) \$376
- ☐ PNG50536 Lateral Supports with Flat Pads-11"-19.5"W (pair) \$376
- ☐ PNG50224 Lateral Supports with Curved Pads-8"-16.5"W (pair) \$376
- ☐ PNG50176 Lateral Supports with Curved Pads-11"-19.5"W (pair) \$376

Head Supports-back and accessories mounting bracket required

- ☐ PNG50240 Head Support-Short \$388
- ☐ PNG50037 Head Support-Tall \$388

Other Accessories

- ☐ PNG50040 Push Handles (back required, pair) \$159
- ☐ P80809 Tool Pouch \$27

* not available with mobile

(A) Seat Depth Range: 20"-25"
(seat pivot to flat back)

Measure the distance between the back of the buttocks and the back of the knee.

(B) Seat to Foot Plate Range: 17"-23"
(from seat pivot)

Measure the distance between the foot/shoe and the knee bend.

**AltimateMedical**

Ph: 800.342.8968 or 507.697.6393
 Fax: 877.342.8968 or 507.697.6900
 info@easystand.com easystand.com

SubmitEmail form to webquotes@easystand.com