

Bill To (Ultimate Medical Supplier Only)

Purchase Order #: _____

Quote #: _____

Account #: _____

Ordered By: _____

Supplier Company: _____

Address: _____

City, State, Zip: _____

Phone: _____ Email: _____

Payer Source: _____

Name of ATP (Required): _____

Ship To

Name: _____

Address: _____

City, State, Zip: _____

Phone: _____

Tag For: _____

Evolv Large User Specifications

Approximate Height Range 5'0"-6'2"

Weight Limit 280 lbs.

Seat Depth Range (seat pivot to flat back) 18"-23"

Seat to Foot Plate Range (from seat pivot) 14"-21"

Base *Note: The DEFAULT items will be added if no selection is made.*

■ **PNG50084-1 EasyStand[®] Evolv[®] Large** \$2,182

Lift Mechanisms-must select one

☐ PNG50009 Manual Hydraulic Actuator **DEFAULT** N/C

Manual Hydraulic Actuator Handle Location

☐ PNG50059 Right Side **DEFAULT** N/C

☐ PNG50126 Left Side N/C

Manual Hydraulic Actuator Handle Accessories

☐ PNG50328 Adjustable/Removable Actuator Handle \$228

☐ PNG50044 Quad Grip Handle Extension \$153

☐ PNG50070 T-Style Handle Extension \$153

☐ PNG50010 Pow'r Up Lift (battery powered) \$1,222

Front Frames and Trays**MUST SELECT ONE FRAME AND ONE TRAY IN THE SAME FRAME SECTION****Shadow Tray Frame-back required***

☐ PNG51042R Swing-Away Front-Right for Shadow Tray with Wheels \$580

☐ PNG51042L Swing-Away Front-Left for Shadow Tray with Wheels \$580

☐ PNG51044R Swing-Away Front-Right for Shadow Tray with Swivel Casters \$759

☐ PNG51044L Swing-Away Front-Left for Shadow Tray with Swivel Casters \$759

☐ PNG51050 Column Front for Shadow Tray with Wheels \$369

☐ PNG51052 Column Front for Shadow Tray with Swivel Casters \$548

Shadow Tray*

☐ PNG50515 Black Molded Adjustable Shadow Tray-Right \$1,415
(7"-12"H from seat, 4"-20"D from back)

☐ PNG50537 Black Molded Adjustable Shadow Tray-Left \$1,415
(7"-12"H from seat, 4"-20"D from back)

☐ PNG50517 Clear Adjustable Shadow Tray-Right \$1,579
(6"-11"H from seat, 4"-20"D from back)

☐ PNG50538 Clear Adjustable Shadow Tray-Left \$1,579
(6"-11"H from seat, 4"-20"D from back)

Swing-Away Frame*

☐ PNG51046R Swing-Away Front-Right Mounted Trays with Wheels \$639

☐ PNG51046L Swing-Away Front-Left Mounted Trays with Wheels \$639

☐ PNG51048R Swing-Away Front-Right Mounted Trays with Swivel Casters \$818

☐ PNG51048L Swing-Away Front-Left Mounted Trays with Swivel Casters \$818

Swing-Away Tray*

☐ PNG50418-1 Black Molded Tray \$585

☐ PNG50456-1 Clear Tray \$753

☐ PNG50458-1 Oversized Clear Tray \$892

☐ PNG50490-1 Black Molded Angle Adjustable Tray \$934

☐ PNG50489-1 Clear Angle Adjustable Tray \$1,095

☐ PNG50491-1 Oversized Angle Adjustable Tray \$1,232

Column Frame*

☐ PNG51051 Column with Wheels **DEFAULT** \$369

☐ PNG51053 Column with Swivel Casters \$548

Column Tray*

☐ PNG50247-1 Black Molded Tray **DEFAULT** \$585

☐ PNG50248-1 Clear Tray \$753

☐ PNG50194-1 Oversized Clear Tray \$892

☐ PNG50368-1 Black Molded Angle Adjustable Tray \$934

☐ PNG50366-1 Clear Angle Adjustable Tray \$1,095

Mobile Frame

☐ PNG51059 Column with Mobile Option \$1,838

Mobile Tray

☐ PNG50372-1 Narrow Clear Tray \$753

☐ PNG50327-1 Narrow Clear Angle Adjustable Tray \$1,095

Tray Accessories

☐ PNG50042 Hand Grips (for clear tray only, pair)* \$158

☐ PNG30295 Padded Tray Cover (not available with oversized clear trays or mobile trays)* \$226

Chest Pad-must select one

☐ PNG30001	Medium Chest Pad DEFAULT	N/C
☐ PNG30000	Large Contoured Chest Pad	\$142
☐ PNG50331	Large Contoured Adjustable Chest Pad	\$204

Knee Pads-must select one

☐ PNG50460-1	1-Piece Knee Pad (pair) DEFAULT	\$328
☐ PNG50384-1	Independent Knee Pads (pair)	\$814
☐ PNG50385-1	Independent Roho® Knee Pads (pair)	\$936
☐ PNG50424-1	1-Piece Knee Pad for Swing-Away Front (pair) *	\$328
☐ PNG50425-1	Independent Knee Pads for Swing-Away Front (pair) *	\$814
☐ PNG50426-1	Independent Roho® Knee Pads for Swing-Away Front (pair) *	\$936

Foot Plates-must select one

☐ PNG50391-1	Height-Adjustable Foot Plates (pair)	\$297
☐ PNG50015-1	Multi-Adjustable Foot Plates (pair) DEFAULT	\$330

Foot Plate Accessories

☐ PNG30030	Foot Straps-15"L (length over top of foot, pair)	\$48
☐ PNG50505	Foot Straps-21"L (length over top of foot, pair)	\$48
☐ PNG30031	Secure Foot Straps-15"L (length over top of foot, two pair)	\$121
☐ PNG50506	Secure Foot Straps-21"L (length over top of foot, two pair)	\$121

Seats-must select one

☐ PNG30006-1	Planar Seat DEFAULT	\$249
☐ PNG30270	Transfer Seat (replaces standard planar seat)	\$388
☐ PNG50318	Rotating Seat (replaces standard planar seat)	\$490
☐ PNG30023	Seat with Roho® Insert (replaces standard planar seat)	\$609

Seat Accessories

☐ PNG50499	Easy-Adjust Seat Depth-CANNOT be ordered off-unit (roller bearing allows for smoother seat depth adjustment)	\$108
-----------------------------------	--	-------

Pelvic Support

☐ PNG50161	Hip Supports (10"-15" range, pair)	\$363
☐ PNG50198	Hip Supports (12"-17" range, pair, limited adjustment with shadow tray)	\$363
☐ PNG50518	Hip Supports (17"-22" range, pair, limited adjustment with shadow tray)	\$363
☐ PNG30029	Velcro® Positioning Belt-26"-55" (back required)	\$62
☐ PNG30306	Velcro® Positioning Belt-38"-67" (back required)	\$62
☐ PNG30028	Positioning Belt with Airline Style Buckle-26"-55" (back required)	\$153

Backs

☐ PNG50066	Flat Back-15"H	\$665
☐ PNG50067	Flat Back-19"H	\$702
☐ PNG50062	Contoured Back-15"H	\$732
☐ PNG50068	Contoured Back-19"H	\$769
☐ PNG50171	Removable Flat Back-15"H	\$763
☐ PNG50173	Removable Flat Back-19"H	\$800
☐ PNG50172	Removable Contoured Back-15"H	\$830
☐ PNG50174	Removable Contoured Back-19"H	\$868

Back Accessories

☐ P600002	Back Angle Adjustment Wheel (back required)	\$104
☐ P80246	Angle Locator	\$38

Chest Straps-back required

☐ PNG50192	Chest Strap-30"-46"	\$153
☐ PNG50193	Chest Strap-34"-54"	\$153
☐ PNG50504	Chest Strap-46"-66"	\$153

Upper Body Support-back required**ATTENTION - ALL OPTIONS BELOW REQUIRE THE ACCESSORIES MOUNTING BRACKET**

☐ PNG50034	Accessories Mounting Bracket	\$68
-----------------------------------	------------------------------	------

Chest Vests-back and accessories mounting bracket required

☐ PNG11212	High Mount Chest Vest Bracket	\$78
☐ PNG30170	X-Style Chest Vest-11"Lx9.5"W	\$286
☐ PNG30056	X-Style Chest Vest-14.5"Lx10.5"W	\$286
☐ PNG30171	Y-Style Chest Vest-17.5"Lx11.5"W	\$286
☐ PNG30021	Y-Style Chest Vest-20.5"Lx12.5"W	\$286
☐ PNG30318	Non-Stretch X-Style Chest Vest-11"Lx9.5"W	\$333
☐ PNG30322	Non-Stretch X-Style Chest Vest-14.5"Lx10.5"W	\$333
☐ PNG30320	Non-Stretch Y-Style Chest Vest-17.5"Lx11.5"W	\$333
☐ PNG30324	Non-Stretch Y-Style Chest Vest-20.5"Lx12.5"W	\$333

Lateral Supports-back and accessories mounting bracket required

☐ PNG50535	Lateral Supports with Flat Pads-8"-16.5"W (pair)	\$363
☐ PNG50536	Lateral Supports with Flat Pads-11"-19.5"W (pair)	\$363
☐ PNG50224	Lateral Supports with Curved Pads-8"-16.5"W (pair)	\$363
☐ PNG50176	Lateral Supports with Curved Pads-11"-19.5"W (pair)	\$363

Head Supports-back and accessories mounting bracket required

☐ PNG50240	Head Support-Short	\$374
☐ PNG50037	Head Support-Tall	\$374

Other Accessories

☐ PNG50040	Push Handles (back required, pair)	\$153
☐ P80809	Tool Pouch	\$26

* not available with mobile

(A) Seat Depth Range: 18"-23"
(seat pivot to flat back)

Measure the distance between the back of the buttocks and the back of the knee.

(B) Seat to Foot Plate Range: 14"-21"
(from seat pivot)

Measure the distance between the foot/shoe and the knee bend.

**AltimateMedical**

Ph: 800.342.8968 or 507.697.6393

Fax: 877.342.8968 or 507.697.6900

info@easystand.com easystand.com

Email form to webquotes@easystand.com