



HEAVY DUTY PVC POOL ACCESS CHAIR



SPECIFICATIONS

THE HEAVY DUTY PVC POOL ACCESS CHAIR SHALL BE CONSTRUCTED OF $\phi 1 \frac{3}{8}$ " [$\phi 4.3$ cm O.D.] SCHEDULE 40 FURNITURE GRADE PVC PIPE AND FITTINGS. THE CHAIR SHALL HAVE A WEIGHT CAPACITY OF 400lbs [181 kg]. THE CHAIR SHALL COME EQUIPPED WITH AN EXTENDABLE FOOTREST. THE SEAT SHALL BE 24" [61 cm] WIDE AND CONSTRUCTED OF MILDEW RESISTANT PVC MESH. THE CHAIR SHALL COME EQUIPPED WITH DUAL LOCKING BRAKES AND A SWING-OUT ARM REST TO ALLOW FOR EASY TRANSFERS. THE REAR WHEELS SHALL BE $\phi 24$ " [61 cm] URETHANE NON-SKID WHEELS WITH USER GRIPS. THE FRONT WHEELS SHALL BE 3" [7.6 cm] RUBBER WHEELS WITH STAINLESS STEEL SWIVELS. THE CHAIR SHALL HAVE AN ANTI-TIP DEVICE.

FRAME: $\phi 1 \frac{3}{8}$ " [$\phi 4.3$ cm O.D.] FURNITURE GRADE PVC PIPE & FITTINGS
 UPHOLSTERY: PVC MESH - MEETS CALIFORNIA 117 BULLETIN FOR FIRE RETARDENCY.
 REAR WHEEL: 24" X 1" [61 X 2.5 cm] URETHANE TIRE ON PLASTIC 'MAG' STYLE WHEEL.
 FRONT CASTER: $\phi 3$ " [$\phi 7.6$ cm] RUBBER WITH STAINLESS STEEL SWIVELS.
 ARM RESTS: SWING-OUT ARMREST FOR TRANSFERS
 FOOT RESTS: SLIDING PVC FOOTREST SLIDES IN FOR CONVIENCE
 HARDWARE: STAINLESS STEEL




Aqua Creek Products
 The Total Solution for Pool and Spa Access

9889 GARRY MORE LANE
 MISSOULA, MT 59808
 (406) 549-0769
 FAX (406) 549-2602

NOTE:
 *SPECIFICATIONS ARE SUBJECT TO CHANGE WITHOUT NOTICE
 *ALL DIMENSIONS ARE IN INCHES

PRODUCT:	F-400PWC		
DWG #:	F-400PWC-1		
SCALE:	1:16	DATE:	12/3/15
DRAWN BY:	EF	REVISION:	E