

FST-300 Free Standing Track

Take Your Track Anywhere

Are you heading on vacation and need to bring along a lift? Or just want the flexibility to switch rooms at home? Handicare's FST-300 Freestanding track system is both lightweight and durable. For use with our reliable P-440 portable lift, the FST-300 is designed with both caregiver and user in mind. Ready out of the box for lifts up to 300 lb, it can be moved and assembled by one person in minutes for easy set-up in any space.

FST-300 Free Standing Track Specifications

Part Numbers

341400: FST-300

341401: Roller bag

Material

Aluminum

Safe Working Load (SWL)

300 lb (136 kg)

Dimensions (L x W x H)

Length: 96 in (2440 mm) Width, base: 36 in (910 mm)

Height, minimum: 60 in (1520 mm) Height, maximum: 90 in (2290 mm) Height, adjustable: 2 in (51 mm) increments

Weight

Unit: 45 lb (20.4 kg), including trolley

Track: 15 lb (6.8 kg)

Post and leg assembly: 12 lb (5.4 kg)

Safety


Designed and tested to be compliant to ISO 10535:2006

Accessories

Roller bag: components fold easily to fit

For use with a Handicare brand portable lift that has a safe working load (SWL) of up to 300 lb, or with a lift that has a higher safe working load, but the user must obey the SWL of the FST-300.

 ϵ



Are you ready to take mobility to the next level?

Call 1-888-637-8155 or visit www.handicareusa.com

