

AMP HANDRIM



Mark For: _____	QUOTE <input type="checkbox"/>	ORDER <input type="checkbox"/>	# _____
Acct #: _____	Date: _____		
Bill To Name: _____	PO#: _____		
Billing Address: _____	Ship to: _____		
City, State: _____	Address: _____		
Confirmation Email: _____	City, State: _____		
Contact : _____	Zip: _____		

END USER WHEELCHAIR SPECIFICATIONS

Wheelchair Type: <input type="checkbox"/> Folding <input type="checkbox"/> Rigid	Brand: _____	Model: _____
1. Front Seat Height in inches <input style="width: 80px;" type="text"/>	2. Rear Seat Height in inches <input style="width: 80px;" type="text"/>	3. Inside Tire Distance in inches <input style="width: 80px;" type="text"/>

1 Wheelchair front seat height must be a minimum 1/2" greater than rear; to maintain same incline when mounted to AMP add an additional 1/2" in front height of wheelchair as docking increases rear seat height 1/2". 2 To avoid wheelchair interference with the AMP Module the rear lower frame must have 8.5" ground clearance. 3 At bottom of chair inside to inside distance between tires.

AMP HANDRIM & OPTIONS

- | | | | | |
|--|------------|-------------|--|---------------|
| <input type="checkbox"/> | SKU | HCPS | Description | Retail |
| <input type="checkbox"/> | AMP | E0986 | AMP Handrim Power Module w/ Battery & Charger | \$7,995 |
| (Must choose one from the size matrix to the right:) | | | | |
| <input type="checkbox"/> | 380AHRS | | AMP Small Base 14" - 16" Wide | |
| <input type="checkbox"/> | 380AHRM | | AMP Medium Base 16" - 18" Wide | |
| <input type="checkbox"/> | 380AHLR | | AMP Large Base 18" - 20" Wide | |
| <input type="checkbox"/> | 180WEXT | | Width Extension Kit (Adds 2" to Wheel Width Fit) | \$275 |
| Joystick Input | | | | |
| <input type="checkbox"/> | 380JIO | | Joystick Input | \$595 |
| Joystick Mount Side (Must choose one if ordering Joystick option) | | | | |
| <input type="checkbox"/> | 380JIR | | Joystick Input Right | |
| <input type="checkbox"/> | 380JIL | | Joystick Input Left | |
| Swing Away/Attendant Control (Must choose at least one if ordering Joystick option) | | | | |
| <input type="checkbox"/> | 380SAD | | Swing Away Mount Dual Post | \$350 |
| <input type="checkbox"/> | 380SAT | | Swing Away Mount Tubular | \$350 |
| Joystick Handles | | | | |
| <input checked="" type="checkbox"/> | 380SH | | Standard Handle (When ordering Joystick input) | Included |
| <input type="checkbox"/> | 380BPU | | Bodypoint® 3" U-Shaped Handle | \$135 |
| <input type="checkbox"/> | 380BPGM | | Bodypoint® Grooved Mushroom Handle | \$135 |
| Battery/ Charger | | | | |
| <input checked="" type="checkbox"/> | 380LPSTC | | LiFePO4 Battery/ Charger (6.5 Hours) | Included |
| <input type="checkbox"/> | 380USC | | Upgrade Super Charger (2.5 hours) | \$125 |
| Battery Options | | | | |
| <input type="checkbox"/> | 380BULP | | Back Up LiFePO4 Battery | \$995 |
| <input type="checkbox"/> | 380BUSC | | Spare Battery Charger (6.5 Hours) | \$225 |
| <input type="checkbox"/> | 380SSC | | Spare Super Charger (2.5 Hours) | \$335 |
| Military Branch Logo (VA ONLY) | | | | |
| <input type="checkbox"/> | 380ABL | | Army | NC |
| <input type="checkbox"/> | 380NBL | | Navy | NC |
| <input type="checkbox"/> | 380AFBL | | Air Force | NC |
| <input type="checkbox"/> | 380MBL | | Marines | NC |
| <input type="checkbox"/> | 380CGBL | | Coast Guard | NC |

Sub Totals

AMP SIZE MATRIX

NOTE: Choose AMP size based on chair type, along with camber, which determines the AMP size recommended. Each size is adjustable in width, but limited based on camber or overall wheel base width.

CHECK ONE BELOW

FOLDING/RIGID 0/3 CAMBER		WHEEL WIDTH RANGE
CHAIR SEAT WIDTH	AMP SIZE (growth)	INSIDE TO INSIDE OF TIRES
14"	<input type="checkbox"/> SMALL (14"-16")	17 1/4"-19 1/2"
15"	<input type="checkbox"/> SMALL (14"-16")	
16"	<input type="checkbox"/> SMALL (14"-16")	
16"	<input type="checkbox"/> MEDIUM (16"-18")	19" - 21 1/4"
17"	<input type="checkbox"/> MEDIUM (16"-18")	
18"	<input type="checkbox"/> MEDIUM (16"-18")	
18"	<input type="checkbox"/> LARGE (18"-20")	20 1/2" - 23"
19"	<input type="checkbox"/> LARGE (18"-20")	
20"	<input type="checkbox"/> LARGE (No Growth)	

FOLDING/RIGID UP TO 5 CAMBER		WHEEL WIDTH RANGE
CHAIR SEAT WIDTH	AMP SIZE (growth)	INSIDE TO INSIDE OF TIRES
14"	<input type="checkbox"/> SMALL (14"-15")	17 1/4"-19 1/2"
15"	<input type="checkbox"/> MEDIUM (15"-17")	
16"	<input type="checkbox"/> MEDIUM (15"-17")	
17"	<input type="checkbox"/> LARGE (17"-19")	19" - 21 1/4"
18"	<input type="checkbox"/> LARGE (17"-19")	
19"	<input type="checkbox"/> LARGE (No Growth)	
20"	N/A	

For quotes and orders EMAIL to:
CustomerCare@methodmobility.com
 or FAX to: 877-398-4545

Shipping	
Total w/Assess.,Delivery & Shipping	
Credit Card Fee	
Total with Fees (Pay this amount)	



Method Mobility VA Registration: UVRSAU9M5AE7/CAGE #9JCF9

www.methodmobility.com
 888-88AMP88 (882-6788)

MM-AMPOF 11.5.1