



**CORPORATE HEADQUARTERS**

99 Seaview Blvd., Port Washington, NY 11050

Toll Free: 877.224.0946 • 516.998.4600

Fax: 516.708.1104 • [www.drivemedical.com](http://www.drivemedical.com)

Medical Depot Inc. dba Drive DeVilbiss Healthcare  
Inovo, Inc.  
DeVilbiss Healthcare, LLC  
Medical Depot, Inc. Distribution centers NJ CA IN GA

June 5, 2018

### **Battery Regulation for Flight**

To whom it may concern:

Medical Depot, Inc. confirms that the product - FLEX-AUTO - is a medical device.

The Travel Scooter should be stored in the cabin for flight. If the Scooter needs to be stowed with the luggage, the battery needs to be removed and placed with carry-on luggage in the cabin.

The battery type and Safety Data Sheet information for the product is listed below.

Item:	FLEX-AUTO
Battery:	Cylindrical Lithium Ion Battery, Model 24/10
Voltage:	25.9V
SDS Information:	Dong Guan Li Long Battery Technology Co. Ltd
SDS #:	STRD1701027S
Labeling:	UN3480/1
E-Mail:	lxq@chinallb.com
24 HR Emergency Phone:	18929251199

Any questions may be directed to the number above if damage occurs.

Regards,

[Redacted Signature]

[Redacted Contact Information]