

# FLEXION FOOTBOX™ - SINGLE



**HCPCS CODING:** K0108  
WHEELCHAIR ACCESSORY

## PRODUCT FEATURES

- + Available in four sizes: 5"W, 6"W, 7"W, and 8"W.
- + Easy angle adjustment from 90° to 160° with a single screw.
- + Lightweight, durable aluminum shell.
- + After drilling mounting holes in footbox, bolt on to chair for permanent placement.
- + All surfaces and edges are foam padded for added comfort.
- + Comes standard with Comfort-Tek cover.



## STANDARD COVER



**COMFORT-TEK FABRIC:** the fluid proof exterior has supple leathery properties that provide an easily cleaned surface and protection for the inner cushion. Comfort-Tek has multi-dimensional stretch and helps alleviate pressure on bony prominences.

## ANGLE ADJUSTMENT

Easy angle adjustment from 90° to 160° with a single screw.



## FLEXION FOOTBOX™ SINGLE

| PART #       | WIDTH X HEIGHT X DEPTH |
|--------------|------------------------|
| <b>FBS-5</b> | 5.25"W x 11"H x 11"D   |
| <b>FBS-6</b> | 6.25"W x 11"H x 11"D   |
| <b>FBS-7</b> | 7.25"W x 11"H x 11"D   |
| <b>FBS-8</b> | 8.25"W x 11"H x 11"D   |

NOTES

# FLEXION FOOTBOX™ - DOUBLE

**HCPCS CODING:** K0108  
WHEELCHAIR ACCESSORY



## PRODUCT FEATURES

- + Available in four sizes: 14"W, 16"W, 18"W, and 20" chairs.
- + Easy angle adjustment from 90° to 160° with a single screw.
- + Lightweight, durable aluminum shell.
- + After drilling mounting holes in footbox, bolt on to chair for permanent placement.
- + All surfaces and edges are foam padded for added comfort.
- + Comes standard with Comfort-Tek cover.



### FLEXION FOOTBOX™ DOUBLE

| PART #        | CHAIR SIZE | WIDTH X HEIGHT X DEPTH |
|---------------|------------|------------------------|
| <b>FBD-14</b> | 14"        | 12"W x 11"H x 11"D     |
| <b>FBD-16</b> | 16"        | 14"W x 11"H x 11"D     |
| <b>FBD-18</b> | 18"        | 16"W x 11"H x 11"D     |
| <b>FBD-20</b> | 20"        | 18"W x 11"H x 11"D     |

### STANDARD COVER



**COMFORT-TEK FABRIC:** the fluid proof exterior has supple leathery properties that provide an easily cleaned surface and protection for the inner cushion. Comfort-Tek has multi-dimensional stretch and helps alleviate pressure on bony prominences.

### ANGLE ADJUSTMENT

Easy angle adjustment from 90° to 160° with a single screw.



NOTES