

THE FIRST EVER FOLDABLE WHEELCHAIR WHEELS



Photo: Colin Ross



MORPHTM
WHEELS
OPEN POSSIBILITIESTM

OPEN POSSIBILITIES™



- Fail-safe design prevents wheels from folding when in use
- Travel in smaller cars with less storage space – even in taxis
- Reduce the painful strain on your shoulders when transferring into cars
- Fit into overhead bins, closets and under your seat during flights
- Limited spoke design enables you to easily reach under your chair
- Compatible with all manual wheelchairs that use quick release axles



Wheel size: 24"

Dimensions when folded: 32" x 12½"

Wheel set weight capacity: 265 lbs.

Weight per wheel: 7.5 lbs.

Material: glass-filled nylon

Hand rim material: polypropylene

Tires: solid

Passed all RESNA/ANSI impact, curb drop and double drum testing standards


MORPH
WHEELS
morphwheels.com