





Invacare® Top End® Force™ K Handcycle

The **Invacare® Top End® Force™ K Handcycle** is a kneeling position handcycle available in two custom, made-to-measure aerodynamic designs. The Invacare Top End Force K Kneeler is used by the most able athletes, while the Top End Force K Double Amputee Kneeler is custom designed for the double above-the-knee amputee. What makes the Force K Handcycle different from the rest of the Invacare Top End Force Series handcycles? The Force K Handcycle is a kneeling-only, non-adjustable position handcycle that requires dialed-in measurements and engineering drawings—all of which are easy to find on the website: www.topendwheelchair.com. Can you force the pace? Get the Force K Handcycle.

Features

- 7005 aluminum custom welded frame
- 7005 aluminum aerodynamic fork
- Breathable mesh upholstery
- Seat and back cushions
- Tapered crank arms
- Frictionless pedal bearings
- Aluminum ergonomic hand pedals
- SRAM® drive train
- Chainrings: 53-39 or 48-39-24
- Cassettes: SRAM® 11-36
- Tektro® Caliper Brake
- Fork mounted locking second brake
- Steering dampener
- 26" high performance wheels
- 3°, 6°, 9° or 11° of camber
- Threaded axles
- Drafting bumper
- Parking brake
- Safety flag
- Chain guard
- Many paint choices
- Reflective stickers
- 250 lb. weight capacity

Model #: FRCK



1. Top End Carbon Wheel Upgrade is stiff, light and aerodynamic
2. SRAM components allow shifting and braking at your fingertips
3. Adjustable depth seat upholstery in short position
4. Adjustable depth seat upholstery in long position

Options

- Carbon fiber wheels
- Component, brake and cassette upgrades
- SRM™ Power Meter measuring tool integrated directly into Invacare Top End cranks

Frame Styles:

- The Force K Kneeler is designed for the most able athletes with a custom knee/seat bucket
- The Force K Double Amputee Kneeler is designed for the double above-the-knee amputee utilizing a custom cage

See the Force K Handcycle order form for all of the great options available!

Invacare Corporation
www.invacare.com

USA

One Invacare Way
Elyria, Ohio
44035-4190
(800) 333-6900

Canada

570 Matheson Blvd. E.,
Unit 8
Mississauga, Ontario
L4Z 4G4 Canada
(800) 668-5324

Top End

www.topendwheelchair.com
Sports and Recreation
Products
4501 63rd Circle North
Pinellas Park, Florida
33781-5914
(800) 532-8677
Fax: (727) 522-1007

©2018

Invacare Corporation.
All rights reserved.
Trademarks are identified
by the symbols ™ and ®.
All trademarks are owned
by or licensed to Invacare
Corporation unless
otherwise noted. SRM is a
trademark of SRM Service
Center Inc. USA. SRAM
is a registered trademark
of SRAM Corp. Tektro is
a registered trademark of
Tektro USA. Specifications
are subject to change
without notification.
Form No. 18-320
181122

 **TOP END®**