



ORDER FORM

Freestander

Upright accessible standing frame
from Early Intervention to adult.

Customer			
Account			
Account #		Date	
☐ Quote	☐ Order	PO #	
Buyer			
Mark For			

Shipping Address			
Name			
Address			
City		State	
ZIP		Phone	
Email / Fax			

The HCPCS CODES herein are based on PDAC verification or interpretation of Medicare definitions and guidelines. Non-Medicare payers may accept alternative HCPCS CODES, including misc. codes to ensure access for their enrollees. The use of HCPCS CODES does not ensure coverage or payment.

Freestander

Size 1 Weight Capacity: 132 lbs | Size 2 Weight Capacity: 176 lbs | Size 3 Weight Capacity: 198 lbs

✓ OPTION	DESCRIPTION	HCPCS	RETAIL
Freestander			\$
Selection Required. Includes Chest and Hip Pad, Flexible Lateral, PU Knee Supports, Footplate, Tray Cover and Activity Bowl.			
☐ LFS/1-MW	Size 1 Freestander - Meadow	E0638	\$2,818
☐ LFS/2-MB	Size 2 Freestander - Mondeo Blue	E0638	\$3,053
☐ LFS/3-MB	Size 3 Freestander - Mondeo Blue	E0638	\$3,453
Additional Positioning			\$

If Solid Pelvic Plate selected, pelvic band selection required.

☐ LFS/1/B-MW	Size 1 Hip Belt - Meadow		\$189
☐ LFS/2/B-MB	Size 2 Hip Belt - Mondeo Blue		\$213
☐ LFS/3/B-MB	Size 3 Hip Belt - Mondeo Blue		\$231
☐ LFS/1/L-MW	Size 1 Solid Pelvic Plate - Meadow		\$331
☐ LFS/2/L-MB	Size 2 Solid Pelvic Plate - Mondeo Blue		\$374
☐ LFS/3/L-MB	Size 3 Solid Pelvic Plate - Mondeo Blue		\$418
☐ LPS/1/B-MW	Size 1 Pelvic Band - Meadow		\$194
☐ LPS/2/B-MB	Size 2 Pelvic Band - Mondeo Blue		\$213
☐ LPS/3/B-MB	Size 3 Pelvic Band - Mondeo Blue		\$231

Additional Knee Supports			\$
☐ LPS/1/K-MW	Size 1 Upholstered Knee Cups - Meadow (Pair)		\$132
☐ LPS/2/K-MB	Size 2 Upholstered Knee Cups - Mondeo Blue (Pair)		\$132
☐ LPS/3/K-MB	Size 3 Upholstered Knee Cups - Mondeo Blue (Pair)		\$132
☐ LPS/4/K-MB	Size 4 Upholstered Knee Cups - Mondeo Blue (Pair)		\$132

Additional Foot Positioning			\$
-----------------------------	--	--	----

Sandal selection required with Sandal Raiser.

☐ 152-1600	Small Sandals with Straps (Pair)		\$273
☐ 152-2600	Medium Sandals with Straps (Pair)		\$298
☐ 152-3600	Large Sandals with Straps (Pair)		\$321
☐ 117-787	QTY Sandal Raiser 1" (Each) Requires Sandal Selection.		\$95

Technical Specifications - Imperial

Frame Sizing Specification	Size 1	Size 2	Size 3
Age (approximate)	1 1/2 - 5 years	4 - 10 years	9 - 18 years
User Body Height (Min / Max)	31.5" / 45.3"	43.3" / 57"	55.1" / 72.8"
Maximum User Weight	132 lbs	176 lbs	198 lbs
Chest Pad Height (Min / Max)	20.7" / 33.5"	26.8" / 43.3"	38.4" / 56.1"
Width Between Chest Supports (Min / Max)	6" / 10"	8" / 12"	10" / 14"
Activity Tray	16" x 21"	18" x 23"	20" x 25"

Technical Specifications - Metric

Frame Sizing Specification	Size 1	Size 2	Size 3
Age (approximate)	1 1/2 - 5 years	4 - 10 years	9 - 18 years
User Body Height (Min / Max)	800 mm / 1151 mm	1100 mm / 1448 mm	1400 mm / 1849 mm
Maximum User Weight	60 kgs	80 kgs	90 kgs
Chest Pad Height (Min / Max)	470 mm / 680 mm	680 mm / 975 mm	975 mm / 1350 mm
Width Between Chest Supports (Min / Max)	152 mm / 254 mm	203 mm / 305 mm	254 mm / 356 mm
Activity Tray	400 mm x 505 mm	450 mm x 575 mm	500 mm x 625 mm

Sunrise Medical (US) LLC

2842 Business Park Avenue
Fresno, CA 93727 USA
Tel.: 1-800-333-4000
Fax: 1-800-300-7502
E-Mail: webmaster@sunmed.com



If you are visually impaired, this document can be viewed in PDF format at www.SunriseMedical.com



MK-129847 Rev. E ©01/2024



800.333.4000

www.SunriseMedical.com

