


450 lb Capacity Fusion™

The Ultimate Big & Tall Scooter



Extraordinary Indoor Maneuverability
33 3/4" Turning Radius.

Exceptional Outdoor Capability
10° Stability With 450 lb Rider.

Awesome All Day Power.
Climb Tortuous 140 yd Long,
7-9° Hills 14 Times With a
450 lb Rider.

Luxurious Comfortable Seating
21" Wide Seating Plus 3 "
Extra Leg Room.

Separates For Car Transport.

Top Rated Reliability and
American Made With Kansas City
Stocked Parts and Longest
Warranties.



PaceSaver®
Scooters!

Set Your Own Pace!

PaceSaver® Fusion™ 450

**ALWAYS
AMERICAN MADE**

Specifications	Fusion™ 450
Overall Length	48"
Overall Width (with 18" seat)	24"
(with 21" seat)	26"
Weight	
Rear Drive Unit	56 lbs.
Base & Tiller	39 lbs.
Seat Section	27 lbs. (30 lbs. with 21" seat)
Batteries	27 lbs. (each)
Shroud	3.5 lbs.
Maximum Speed	4.5 mph
Range	450 lb. rider-up to 16 miles
Ground Clearance	2"
Turning Radius	33 3/4"
Tire Size	
Front	3" x 9"
Rear	3" x 9"
Seat Dimensions	
Width	18" or 21"
Depth	18" or 21"
Back Rest Height	16" std.
Seat Height (From Floorpan)	16.75" to 17.75"
Electronic Brakes (Regenerative & Dynamic with "Posi-Lock" Disk Brake)	Yes
Drive Type	Transaxle
Motor Size	1 1/4 HP Motor
Incline Stability	10° (17.6% grade)
Rated Structural Capacity	450 lbs.

Specifications tested to US Ansi/Resna procedures. Incline stability rating by US Ansi/Resna 8.4. Information current at time of publication. We reserve the right to change specifications without public notice. Speed and range may vary with rider weight, battery condition and terrain.

Lifetime Limited Warranty on main frame.

Five Year Limited Warranty on transaxle.

Two Year Limited Warranty on components.

Battery Warranty by battery supplier.

See Owners Manual for complete warranty details.

PaceSaver®
American Made

Color Choices

Ocean Blue

Midnight Black

Candy Apple Red

Velvet Seats

Black

Gray

Plush Vinyl Seats

Charcoal

Gray

Accessories & Options



Flag



Oxygen
Holder



Cup
Holder



Saddle
Bag



Front
Bag

- On-Board Charger
- Cane Holder
- Left Hand Controls
- Lap Belt
- Rearview Mirror
- Cane Seating
- Flat-Free Tires

**Reclining Admiral Seat
with Adjustable Headrest**

Ocean Blue



Please see the Accessories & Options brochure for a complete listing.

Celebrating 23 years of legendary performance. PaceSaver® Scooters!

PaceSaver® is a registered trademark and
Plus III™ is a trademark of Leisure-Lift, Inc.
© 2005 Leisure-Lift, Inc. All rights reserved.

www.pacesaver.com

800-255-0285

leisure-lift@kc.rr.com

Leisure-Lift, Inc.
1800 Merriam Lane
Kansas City, KS 66106
Fusion450