

G5 Pro



39lb Remote Controlled Electric Wheelchair with Suitcase Handle




Battery	Range	Total Weight	Weight w/o Batteries	Transport Weight	Capacity	Motor	Seat Width	Chair Width	Airline Cruise	Max Speed	Max Climb	Turn Radius
10ah Lithium Ion	Up to 11 mi*	46 lbs.	42 lbs.	38 lbs.	265 lbs	2 x 250W Brushless	18.5"	23"	Approved with 10ah	4 mph	9°	47"

- ✓ Lightweight aluminum frame + 265lb capacity + 18.5in seat space
- ✓ Dual 250W brushless whisper quiet motors + airline approved
- ✓ A sleek LCD display joystick
- ✓ Handheld remote control for caregivers



FRAME



MOTOR & ELECTRONICS

LCD Display



Remote Control

