

[Safety Data Sheet Documents](#)

No SDS content is available for this product.

Tub Grab Bars
BAR, GRAB, BATHTUB, BI-LEVEL, GUARDIAN

- Tub grab bar has a step-through clamp design that provides unobstructed entry/exit path for better hand positioning
- Tool-free installation with rubber-lined pads that won't mar tub surface
- Adjusts to fit tubs 2-3/4- 6-1/4" (7-16 cm)

Packaging	1 Each / Each
Manuf / Supplier	Medline
Manuf / Supplier #	G98006

SKU **G98006H****Specifications**

Specifications	Values
BRAND	Guardian Signature
Finish	White
HPIS CLASSIFICATION	630_20_30_0
Latex Free	Yes
Length	11
Mount Type	Bath Tub
PRODUCT TYPE	Bath Tub Grab Bar
Texture	Smooth
UNSPSC	42211608
Weight Capacity	250 lb

Unit of Measure	Conversion	Net / Gross Weight (lbs)	Volume (cubic ft)	Dimensions (inch) L x W x H	GTIN
Each (EA)	1.0 Each	0.0 / 4.4	0.791	9.6 x 7.3 x 19.5	10080196176007 / 80196176000
Case (CS)	3.0 Each	0.0 / 27.0	2.707	22.0 x 10.5 x 20.25	40080196176008 / 80196176000

HCPCS Information**HCPCS Code: E0241**

HCPCS Disclaimer: HCPCS codes and Home Health Consolidated Billing codes provided by Medline are intended as general guidelines only. Medline does not guarantee coverage or reimbursement of any products. You must address all coverage and reimbursement issues (including the correctness and accuracy of codes) with your individual payers. It is your responsibility to ensure the accuracy and appropriateness of each claim you submit, in accordance with all applicable payer requirements.

For more information on HCPCS codes, [Click here](#)