



Standard Colors:



Velvet Red
(Installed)



Coastal Blue
(Included in Box)

Color Through Shrouds!

ALL NEW Sporty Blacked-Out Styling!

MicroTurn™ Technology!

Buzzaround 4-Wheel GB150



The next generation of the popular Buzzaround XL 4-Wheel scooter is now available! The new GB150 offers improved standard features to make your riding experience better than ever. It features a large wraparound Delta tiller, higher ground clearance of 3.5". The XL also offers a new stylish, contoured stadium style mesh seat for more comfort and increased airflow. The LED dash is full color with speed and battery indicators, and the basket includes a built-in drink holder. The GB150 offers our new MicroTurn™ technology with an impressive 39" turning radius to navigate even the tightest spaces with ease!

Standard Features:

- Higher 325 lb. Weight Capacity
- Improved Drivability
- Turning Radius of a 3-W with Stability of 4-W
- Color LED Console with Battery Meter & Speedometer
- Higher Ground Clearance
- Front & Rear LED Lighting
- Front & Rear Suspension
- Width & Angle Adjustable Flip-up Armrests
- USB A&C Ports on the Tiller
- Phone Holder
- Front Basket with Carry Handle & Drink Holder
- Quick & Easy Disassembly



FDA Class II Medical Devices designed by Golden Technologies are intended to assist persons with chronic physical disabilities.





LIGHTWEIGHT PORTABLE SCOOTER

Item	GB150
Description	XL
Medicare Code	K0801
Weight Capacity	325 lbs.
Drive Wheels	Rear
Maximum Speed	4.3 mph
Operating Range up to ¹	8 miles
Under Deck Clearance	3.5"
Turning Radius	39"
Type Batteries	2-12AH
Colors	Red, Blue
Freewheel Mode	Yes
Electronic Speed Control	Yes
Electro-Mechanical Brakes	Yes
Charger	Off-board
Controller	50 Amp

Dimensions

Length	42"
Width	22"
Height (ground to top of back of seat)	35.5"-37.5"
Ground to Top of Seat	20.5"-22.5"
Ground to Top of Deck	6"
Deck to Top of Seat	14.5"-16.5"
Number of Seat Height Positions	3
Size of Increments	1"
Front Axle to Rear Axle	30.6"
Weight of Unit: (assembled)³	110 lbs.
Front Body	42.5 lbs.
Rear Body	28 lbs.
Seat w/Arms	16 lbs.
Batteries	18 lbs.

Tires:

PU Foam	
Front	8"
Rear	8"
Rear Anti-Tip Wheels	2"

Standard Seat

Stadium Mesh Back

Seat Back Height - no headrest	17"
Width x Depth	16"W x 15.5"D
Color	Black

Warranty:

Frame	Lifetime Warranty
Drive Train	2 Year Warranty
Electronics ²	2 Year Warranty
Batteries	1 Year Warranty

¹ Battery range at 200 lbs. and will vary due to rider weight, drive surface, terrain, and battery type.

² Electronics warranty excludes batteries.



microTURN

NEW MicroTurn
feature on GB150
offers an impressive
39" turning radius!



Disassembles into 5 Pieces!



Full-color LED Console with
Battery & Speed Indicators



Phone: 800-624-6374

Fax: 800-628-5165

www.goldentech.com



www.facebook.com/goldentechnologies