



**GO LIGHT.  
FOLD FAST.  
SPEND LESS.**

ONLY  
**25.8 lbs.\***

**300 lb.**  
CAPACITY

Carbon fiber construction.  
FREE Memory foam comfort.  
Side-fold design that stays on its wheels for effortless loading, travel, and everyday freedom.



FDA Class II Medical Device

*\*Without cushion, battery or drink holder.*



**LIGHTWEIGHT  
CARBON FIBER**



**300 LB.  
WEIGHT CAPACITY**



**FREE MEMORY  
FOAM CUSHION**



**AIRLINE  
FRIENDLY**





# Designed for Everyday Freedom

**Travel lighter. Fold easier. Move with confidence.**

Lightweight carbon fiber construction, innovative side-fold convenience, memory foam comfort, and thoughtful everyday features come together in a power wheelchair designed for effortless mobility.



Folded  
Front View  
with Base Cushion

### TRAVEL READY

Compact 12" side-fold design stays on its wheels for easier loading, transport and storage.

**Airline-approved side-mounted lithium-ion battery.**



**PORTABLE POWER**  
FREE battery carry bag for convenient storage.



**EVERYDAY COMFORT**  
FREE memory foam cushion and integrated drink holder.



**EASY TRANSFERS**  
Flip-up armrests simplify transfers.

### SPECIFICATIONS

Model Number	GP309
Medicare Code	K0899
Weight Capacity	300 lbs.
Drive Wheels	Rear
Maximum Speed	3.7 mph
Operating Range <sup>1</sup>	9.5 miles
Ground Clearance	1.25"
Turning Radius	29.5"
Type Batteries	10AH 24V Lithium-Ion
Color	Black Carbon Fiber
Freewheel Mode	Yes
Electronic Speed Control	Yes
Electro-Mechanical Brakes	Yes
Charger	Off Board
Controller	35 Amp

### DIMENSIONS

Length	36"
Width	24"
Height (ground to top of back of seat)	33"
Ground to Top of Seat	19"

### WEIGHT

Weight of Unit: (assembled w/ memory foam cushion)	31.3 lbs.
Weight of Unit: (assembled w/ base cushion)	29.5 lbs.
Weight of Unit: (minus parts)*	25.8 lbs.
Battery	3 lbs.

### TIRES

Tires	Flat Free
Drive Wheels (Rear)	11"
Casters (Front)	6.5"
Rear Anti-Tip Wheels	2"

### STANDARD SEAT

Back Height - no headrest	16.4"
Width x Depth	18" x 18"
Color	Black

### WARRANTY

Frame	5 Years
Drive Train	13 Months
Electronics <sup>2</sup>	13 Months
Battery	1 Year

<sup>1</sup> Battery range tested at 200 lbs. and will vary due to rider weight, drive surface, terrain and battery type.

<sup>2</sup> Electronics warranty excludes batteries.

\* Excludes battery, cushion, and drink holder.



www.goldentech.com

Literature is current at time of printing.  
Golden Technologies reserves the right to make changes to the product or literature at any time.