

# R82 Gazelle - E0641

Quote

Order

## USD Retail Price List/ Order Form


Account #: \_\_\_\_\_ Purchase Order No.: \_\_\_\_\_ Date: \_\_\_\_\_  
 Contact: \_\_\_\_\_ Mark For: \_\_\_\_\_

Bill To:	Ship To:
Address:	Address:
City/State/Zip:	
Phone:	City/State/Zip:
Email	Phone:

### US Configuration

Effective: 01.01.22 - Rev 7

- ☐ Gazelle
- ☐ Head support
- ☐ Trunk support
- ☐ Hip/side support
- ☐ Tray



A (Total Height) \_\_\_\_\_  
 B (Shoulder Height) \_\_\_\_\_  
 C (Chest Height) \_\_\_\_\_  
 D (Buttock Height) \_\_\_\_\_  
 E (Chest Width) \_\_\_\_\_  
 F (Shoulder Width) \_\_\_\_\_  
 Age \_\_\_\_\_  
 Weight \_\_\_\_\_

### Measurements

	Size 1 Inches	Size 2 Inches	Size 3 Inches
Overall width	19"	23"	26"
Length, total	30"	35"	39"
Length, foot-axil	18-37"	24-41"	28-49"
Overall height, horizontal	30"	31½"	35"
Folded height	11"	12"	12"
Folded length	37"	44½"	51"
Height of buttock	14½"	18"	22"
Angle, board	90° to -15°	90° to -15°	90° to -15°
Total mass	44 lbs.	48 lbs.	52 lbs.
Max load	88 lbs.	132 lbs.	176 lbs.

### NOTES:

## Gazelle Standard Features



- Prone, upright and supine in one unit
- Multi adjustable foot and knee (reversible) support
- Foldable bracket for back support
- Removable back support
- Trunk safety belt
- Locking castors
- Manual leg abduction
- Gas spring assisted tilt from vertical to horizontal
- Adjustable trunk or chest support
- Up to 30° of Hip Abduction

Description	Size	Part No.	USD Price	HPCPS
<input type="checkbox"/> Gazelle PS, Yellow	1	88010-1	\$3,561	E0641
<input type="checkbox"/> Gazelle PS, Green	2	88010-2	\$3,619	E0641
<input type="checkbox"/> Gazelle PS, Blue	3	88010-3	\$3,723	E0641

Accessories	Description	Measurements	Size	Part No.	USD Price	HPCPS
-------------	-------------	--------------	------	----------	-----------	-------



<input type="checkbox"/> Trunk Support Extension	4 x 9¾"	1	88106-1	\$219
<input type="checkbox"/> Trunk Support Extension	4 x 11"	2	88106-2	\$219
<input type="checkbox"/> Trunk Support Extension	4 x 12½"	3	88106-3	\$219

Description	Measurements	Size	Part No.	USD Price	HPCPS
-------------	--------------	------	----------	-----------	-------



<input type="checkbox"/> Tray, Plastic, w/ Bowl & Top Plate 24 x 22½"	1, 2	88108-0	\$527
<input type="checkbox"/> Tray, Plastic, w/ Bowl & Top Plate 24 x 22½"	3	88108-1	\$527



<input type="checkbox"/> Armrest for Plastic Tray	One Size	861131	\$167
---	----------	--------	-------



<input type="checkbox"/> Toybar for Plastic Tray	One Size	861135	\$250
--	----------	--------	-------

<input type="checkbox"/> Grip Handle	One Size	8910750	\$422
--------------------------------------	----------	---------	-------

Description	Measurements	Size	Part No.	USD Price	HPCPS
-------------	--------------	------	----------	-----------	-------



<input type="checkbox"/> Lower Leg Angle Adjustment Bracket	1	88165-1	\$506
<input type="checkbox"/> Lower Leg Angle Adjustment Bracket	2	88165-2	\$506
<input type="checkbox"/> Lower Leg Angle Adjustment Bracket	3	88165-3	\$506

<input type="checkbox"/> Strap for Foot Support	One Size	85096	\$38
---	----------	-------	------

Description	Measurements	Size	Part No.	USD Price	HPCPS
-------------	--------------	------	----------	-----------	-------

<input type="checkbox"/> Lateral Supports	5¾ x 5¾"	Small	88055-1	\$230
<input type="checkbox"/> Lateral Supports	7 x 5"	Large	88055-2	\$230






<input type="checkbox"/> Lateral Supports (Exterior) ♦	5¾ x 5¾"	Small	88065-1	\$230
<input type="checkbox"/> Lateral Supports (Exterior) ♦♦	7 x 5"	Large	88065-2	\$230

♦ Extends width by 4"




♦♦ Custom order; Approx. lead time: 6-8 weeks.


<input type="checkbox"/> Lateral Supports (Interior) ♦	5¾ x 5¾"	Small	88075-1	\$230
<input type="checkbox"/> Lateral Supports (Interior) ♦	7 x 5"	Large	88075-2	\$230

♦ Reduces width by 4"

Accessories	Description	Measurements	Size	Part No.	USD Price	HCPCS
	<input type="checkbox"/> Straps for Side Supports		1	88250-1	\$69	
	<input type="checkbox"/> Straps for Side Supports		2	88250-2	\$69	
	<input type="checkbox"/> Straps for Side Supports		3	88250-3	\$69	
	<input type="checkbox"/> Mounting for Head Support ♦			88130	\$214	
	♦ <b>Required when ordering any head support</b>					
	<input type="checkbox"/> Head Supp. w/ Facial Aperture •• 8¼ x 6¼"		1	88117-1	\$245	
	<input type="checkbox"/> Head Supp. w/ Facial Aperture •• 10 x 8¼"		2	88117-2	\$245	
	<input type="checkbox"/> Head Supp. w/ Facial Aperture •• 11½ x 10¼"		3	88117-3	\$245	
	•• <b>Custom order; Approx. lead time: 6-8 weeks.</b>					
	<input type="checkbox"/> Head Support No. 1, Black PU		One Size	6171550	\$214	
	<input type="checkbox"/> Head Support No. 2, Black PU		One Size	6172550	\$214	
	<input type="checkbox"/> Head Support No. 3, Black PU		One Size	6173550	\$214	
	<input type="checkbox"/> Anatomic Head Support ♦	A:5½ x B:5½ x C:9½"	Small	5099001	\$266	E0955
	<input type="checkbox"/> Head Support Cover (Required)		Small	9923-83	\$0	
	<input type="checkbox"/> Anatomic Head Support ♦	A:6 x B:5½ x C:10"	Medium	5099002	\$266	E0955
	<input type="checkbox"/> Head Support Cover (Required)		Medium	9945-83	\$0	
	<input type="checkbox"/> Anatomic Head Support ♦	A:6½ x B:5½ x C:12½"	Large	5099003	\$266	E0955
	<input type="checkbox"/> Head Support Cover (Required)		Large	9967-83	\$0	
	♦ <b>Turtle Bar System Short or Long Required; 98619-25 or 98619-26</b>					
	<b>Swing-away Mounting Brackets - Height, Depth, Angle and Lateral Adjustable</b>					
	<input type="checkbox"/> Turtle Bar System ♦	A: 5", B: 4-8½"	Short	98619-25	\$240	E1028
	<input type="checkbox"/> Turtle Bar System ♦	A: 10", B: 4-8½"	Long	98619-26	\$240	E1028
	♦ <b>Required with Anatomic Head Support; 5099001, 5099002 or 5099003</b>					
	<input type="checkbox"/> Universal Headrest Mount		One Size	87590	\$64	E1028
	♦ <b>Required when using other manufactures Head Support</b>					
	<b>Allow most Whitmyer, Stealth, Otto Bock headrest, I2I to mount to our products (x:panda, Serval, etc.). Please contact Customer Service for confirmation</b>					

**Other types of headrest available, please see Head Supports Order Form**

Spare Parts	Description	Measurements	Size	Part No.	USD Price	HCPCS
	<input type="checkbox"/> Back Support		1	88151-1	\$432	
	<input type="checkbox"/> Back Support		2, 3	88151-2	\$432	
	<input type="checkbox"/> Knee Support (Pr.)		1	88161-1	\$584	
	<input type="checkbox"/> Knee Support (Pr.)		2	88161-2	\$584	
	<input type="checkbox"/> Knee Support (Pr.)		3	88161-3	\$584	
	<input type="checkbox"/> Foot Support (Pr.)		1	88145-1	\$427	
	<input type="checkbox"/> Foot Support (Pr.)		2	88145-2	\$584	
	<input type="checkbox"/> Foot Support (Pr.)		3	88145-3	\$584	

Spare Parts	Description	Measurements	Size	Part No.	USD Price
	<input type="checkbox"/> Trunk Strap (Ea.)		1	88200-1	\$140
	<input type="checkbox"/> Trunk Strap (Ea.)		2	88200-2	\$140
	<input type="checkbox"/> Trunk Strap (Ea.)		3	88200-3	\$140

*HCPCS codes are provided for your convenience and are NOT to be considered billing or legal advice. Providers are responsible for insuring that all appropriate billing codes are used when submitting claims. Proper use of HCPCS codes does not insure coverage or payment. Contact the appropriate payer for coverage information.*

**Prices in USD & is subject to change without notice**

**For Quotes and Orders please send completed Order Form to: [sales.us@etac.com](mailto:sales.us@etac.com) or fax: 1-800-760-0537**

Etac  
13137 Bleinheim Lane  
Matthews, North Carolina 28105  
+1 844-US-MOBILITY (844-876-6245)  
[www.etac.us.com](http://www.etac.us.com)

OF5400-ENG-REV.07  
  
Creating Possibilities