

GO GO[®]

ELITE TRAVELLER GEN 3



FDA Class II
Medical Device¹



NEW deck design that offers more leg room and increased seating comfort.



Intuitive control panel framed by a soft-touch delta tiller handle.



Tiller angle adjustment knob.

AN ELITE TRAVEL MOBILITY EXPERIENCE

Thanks to Pride Mobility's patented, intelligent-turning EZ Turn technology, the new **Go Go Elite Traveller[®] Gen 3** is a 4-wheel scooter that boasts a turning radius of a typical 3-wheel. You'll experience tight turning without compromising 4-wheel stability.

Perfect for all-day travel or activities around the house, the Go Go Elite Traveller[®] Gen 3 features a modern deck design that offers more leg room and the most comfortable seating position available!



SHOWN IN **CARDINAL RED**

GO GO ELITE TRAVELLER® GEN 3

FEATURES

- **New sloped deck** design for improved comfort and **leg room**
- **Front and rear pathway lights**
- **EZ Turn** technology provides unrivaled tight turning without sacrificing the stability of a 4-wheel scooter
- Feather-touch **disassembly** allows for easy separation into 5 lightweight pieces for transport or storage
- Angle-adjustable tiller mechanism
- Interchangeable shroud color panels (comes standard with Cardinal Red)
- Large front basket (standard)



FDA Class II Medical Device¹

ACCESSORIES

- 2.1 amp XLR USB-A Charger
- 2.1 amp XLR USB-C Charger
- Cell phone holder
- Cup holder
- Lap belt
- Military patches
- Oxygen holder
- Rear basket
- Rear view mirror
- Saddle bag
- Safety flag
- Single/double rear mount 3" cane/crutch holder
- Under seat storage
- Walker holder
- Weather cover

SPECIFICATIONS

WEIGHT CAPACITY	300 lbs.	
MAXIMUM SPEED²	Up to 4.23 mph	
DIMENSIONS⁴,⁵	Width: 19.45"	Length: 41.5"
TURNING RADIUS⁴	39.2"	
GROUND CLEARANCE	1.29" at motor, 3.5" under deck	
WEIGHT	83.6 lbs. without batteries 103 lbs. with 12 Ah batteries 111 lbs. with 18 Ah batteries	
RANGE ²,³	18 Ah: Up to 14.1 miles at 220 lbs. 18 Ah: Up to 11.3 miles at 300 lbs. 12 Ah: Up to 8.9 miles at 220 lbs. 12 Ah: Up to 7.1 miles at 300 lbs.	

WARRANTY

FRAME: Lifetime limited
ELECTRONICS: 2-year limited

DRIVETRAIN: 2-year limited
BATTERY: 6 months



FDA Class II Medical Device¹

INCLUDED COLOR PANELS (MATTE)



Cardinal Red

ACCESSORY COLOR PANELS (MATTE)



Black Iceberg Blue Robin's Egg Blue Sugar Plum Chell-o

HCPCS Code: K0800

- 1 Pride® FDA Class II Medical Devices are designed to aid individuals with mobility impairments.
- 2 Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition, and tire condition. This specification can be subject to a variance of (+ or -) 10%.
- 3 Tested in accordance with ANSI/RESNA, MC Vol2, Section 4 & ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance. Test conducted at maximum weight capacity.
- 4 Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or -) 3%.
- 5 Dimensions listed are for scooter only and do not include seats or arms. Width may vary depending on seat options and accessories.

HCPCS codes provided should not be considered as legal advice and do not guarantee reimbursement. DME providers are responsible for determining appropriate billing codes when submitting for insurance reimbursement.

NOTE: All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3 %. Special order seats or non-standard dimensions will not be fold-flat and depth-adjustable.

AUTHORIZED DEALER