





TURNS LIKE A 3-WHEEL STABLE LIKE A 4-WHEEL

No four-wheel scooter can turn like the Go Go Elite Traveller® 2 4-Wheel. Thanks to Pride's patented, intelligent-turning, iTurn Technology, it boasts the tightest turning radius of any available 4-wheel scooter at just 37", which means you have 3-wheel scooter agility without compromising 4-wheel scooter stability.

Perfect for all-day travel or activities around the house, the Go Go Elite Traveller 2 4-Wheel has a 300-lb. weight capacity, a top speed of up to 4.23 mph, and a battery range of up to 13.8 miles.⁵





Scan to pridemobility.co





SC442E



Cardinal Red (matte)

Black (matte)

HCPCS Code: K0800

- 1 Pride® FDA Class II Medical Devices are designed to aid individuals with mobility impairments.
- 2 Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition, and tire condition. This specification can be subject to a variance of (+ or -) 10%.
- 3 Tested in accordance with ANSI/RESNA, MC Vol2, Section 4 & ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance. Test conducted at maximum weight capacity.
- 4 Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or -) 3%.
- 5 Range of 13.8 miles tested with 18 Ah battery and 220-lb. max user weight.

HCPCS codes provided should not be considered as legal advice and do not guarantee reimbursement. DME providers are responsible for determining appropriate billing codes when submitting for insurance reimbursement.

NOTE: All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of \pm or - 3 %. Special order seats or non-standard dimensions will not be fold-flat and depth-adjustable.

FEATURES

- 300-lb. weight capacity
- Maximum speed of up to 4.23 mph
- Improved turning radius of 37"
- Delta tiller with adjustment knob and charging port
- Feather-touch disassembly allows the unit to be easily separated into 5 lightweight pieces
- Front frame-mounted seat post for optimal stability
- Dual voltage charger can charge battery pack on-board or off-board
- Removable footboard basket for additional storage
- Front basket (standard)
- 2 sets of removable color-through shrouds included (Cardinal Red and Black)

ACCESSORIES

- 2.1A XLR USB charger
- Cup holder
- Forearm crutch holder
- Military patches
- Oxygen holder
- Rear basket
- Rear mount 3" double cane/crutch holder
- · Rear view mirror

- Saddle bag
- Safety flag
- Single rear mount 3" cane/crutch holder
- Under seat storage
- Walker holder
- Weather cover
- Wishbone crutch holder
- Additional shroud panels

ACCESSORY SHROUD PANEL COLORS (additional charge applies)



SPECIFICATIONS

WEIGHT CAPACITY	300 lbs.
MAXIMUM SPEED ²	Up to 4.23 mph
DIMENSIONS ⁴	Width: 19.45" Length: 39.37"
TURNING RADIUS4	37"
GROUND CLEARANCE	1.34" at motor 3.15" at mid-deck
PRODUCT WEIGHT	83 lbs. without batteries 102.5 lbs. with 12 Ah batteries 112.6 lbs. with 18 Ah batteries
RANGE PER CHARGE ^{2,3}	18 Ah: Up to 13.8 miles at 220 lbs. Up to 10.4 miles at 300 lbs. 12 Ah: Up to 8.8 miles at 220 lbs. Up to 6.6 miles at 300 lbs.

WARRANTY

FRAME: Lifetime limited D ELECTRONICS: 2-year limited B.

DRIVETRAIN: 2-year limited **BATTERY:** 6-month limited



