


S I E S T A S E R I E S



**ALL NEW!** Power Headrest & Power Lumbar Standard!  
USB Port on Every Hand Control!



Ash

PR445  
shown in  
Hickory

**PR445** Introducing the new Siesta Series power lift recliner! Power headrest and power lumbar come standard, along with a USB port on the soft-feel hand control. Luxurious fabrics feature designer contrast stitching to stand out in any décor. It features an independent back and foot rest and a motion seat for ease of movement. Recline and relax in style with the new PR445!



# PR445

Model Number	PR445
Size	Medium
Weight Capacity	375 LBS
Recline Positions	Infinite
Standard Fabrics	3
Fabric Options	None
Chaise Seat	YES
Coil Springs	YES
Back Type	Seam
Head Rest Cover	No
Arm Covers	No
Battery Back Up	YES
Overall Width	35"
Overall Height	45.5"
Width Between Arms	20.5"
Floor to Top of Seat	20"
Seat Depth	20"
Seat to Top of Back	31.5"
Seat to Top of Arm	6"
Seat Angle Raised	15°
Floor to Top of Foot Rest Reclined	19"
Floor to Top of Back Reclined	30"
Back Angle Closed	105°
Back Angle Reclined	152°
Foot Rest Angle Closed	90°
Foot Rest Angle Opened	173°
Distance Required from Wall Reclined	21"
Extension reclined from Back of Seat to end of Foot Rest	42"
Life Time Warranty, Chair Frame, Lift Frame, Recline Mechanism	YES
Warranty on Electrical Parts	3 Yrs.
Pro Rated Years 4 Thru 7	YES
Warranty on Mechanical Labor	3 Yrs.
Warranty on Electrical Labor	1 Yr.
Removable Back Option Available	YES
Smart Tek System Standard	YES
Customizable	NO
Left Hand Control Option Available	NO
Extra Pocket Option Available	Standard
Heat & Massage Option Available	NO
Gel Memory Foam Option Available	NO
Power Pillow Option Available	Standard
Power Adjustable Lumbar Option Available	Standard
Foot Rest Extension Option Available	NO
Straight Lift Option Available	NO
USB Port	Standard



Shown in Timber

## Features & Benefits

- Independent Back & Footrest
- Motion Seat Reduces Shear
- Formed Head Pillow on back
- Power Headrest & Power Lumbar Standard
- USB port on every hand control
- Soft-feel hand control
- Luxurious fabric
- Contemporary Contrast Stitching



Contrast Stitching



Power Lumbar



Power Headrest with Formed Head Pillow

Sizes may vary depending on fabric, filling material, upholstery and carpet thickness. All measurements taken on concrete floor using levels and metal rulers.



### Smart Tek

In the event of a malfunction, our patented SmartTek™ diagnostic system enables you to get your lift chair up and running quickly.



### Battery Backup System

Batteries not included. Battery Backup System is designed to bring the chair from the reclined to the seated position in the event of a power outage.

Literature is current at time of printing. Golden Technologies reserves the right to make changes to the product or literature at any time.



Complies with  
UL 73 and  
CSA C22.2 No. 68

