



## Features:

- First travel scooter that disassembles with no wires or clips.
- Disassembles into 4 pieces with the heaviest piece being only 22 lbs!



- Includes both a basket and a padded vinyl seat with flip-up arms.
- Maximum speed of 4.0 mph
- Drive range of up to 10 miles per charge
- Height adjustable seat

**Warranty:** Frame - 5 years  
Controller – 2 years  
Motor and Gearboxes – 1 year

Specifications	TREK 3	TREK 4
Overall Length:	38"	40.5
Overall Width:	20"	20"
Per Charge Range* (Max):	Up to 10 miles	Up to 10 miles
Maximum Speed:	Up to 4	Up to 4
Turning Radius:	32"	44.5"
Total Weight Incl. Batteries:	81 lbs	90 lbs
Wheels:	(3) 2" x 7" (front & rear)	(4) 2" x 7" (front & rear)
Ground Clearance:	2 .5"	2 .5"
Battery:	12V 12Ah x 2	12V 12Ah x 2
Drive Train:	24V DC Motor sealed mini-transaxle	24V DC Motor sealed mini-transaxle
Battery Weight:	22 lbs	22 lbs
Battery Charger:	2 Amp (Off-Board)	2 Amp (Off-Board)
Weight Capacity:	250 lbs	250 lbs
Seat Type:	Foldable Molded Foam	Foldable Molded Foam
Standard Seat Weight:	14.5 lbs	14.5 lbs
Body Colors:	Maroon or Blue	Maroon or Blue
Heaviest Piece Weight:	22 lbs	22 lbs
Disassembly:	4 pieces	4 pieces