

SL350AC

Comfort and Convenience



OPTIONAL ADD-ONS



Optional: Grocery Basket



Flip-up arms

The SL350AC Stair Lift is engineered with a carriage that is pre-installed on a short section of the durable extruded aluminum track. The backrest has three different settings, and the fold-up footrest adjusts to two different heights for the utmost in custom comfort. The spacious, wide seat features thick padding and attractive upholstery.

FEATURES & BENEFITS

- Designed for long-term indoor use
- Constant-pressure rocker controls
- Operates on standard 115V AC household current
- Seat and footrest fold to save space
- Durable off-white powder coat finish

COMFORT AND CONVENIENCE

SL350AC



This lift is there when I need it. We have an old two story home and surgery required us to install a lift so that I can use the upstairs. This lift has been dependable and enabling.

20

Feet per Minute

350

*Pound Capacity

LOW

Maintenance

14"

Folding Width from Wall

SL350AC PRODUCT SPECIFICATIONS

Weight Capacity	350 lb* (159 kg)
Speed	20 fpm
Folded Width	14" (35.6 cm)
Drive System	Cable Drive
Power	24V AC
Rail Length	16' (4.9 m) standard
Safety Features	Limit switches, Slack cable brake with switch, Footrest obstruction switch, Swivel seat, Cut-off switch, Constant pressure controls, Seat belt
Warranty	3-Year Parts
Option(s)	Retractable seat belt, Key lock, Extra rail

* 250 lb (113.4) capacity for rail length 16' - 25' (4.9 m - 7.6 m)



Harmar.com | 800.833.0478



22NOV19 | MKT-000066 REV A