



Helio A6

ORDER FORM - USA

P.O. number:* _____ Date : _____
Company Name:* _____ Phone:* _____
Ordered by:* _____
Clinician/ATP/RTS: _____
Email: _____
Marked for: _____

Ship to:* _____
Name: _____
Address: _____
City: _____ State: _____
Zip Code: _____
Phone: _____

* Required

A6MA082 ☐

English Manual USA (STD)

A6MA089 ☐

Spanish Manual USA

Taxes not included

MODEL									USD			
1	A6WC19	Ultralight folding Helio A6 Wheelchair - 6000 series aluminum (265lb weight capacity)								K0005	\$2 095	☐
2	A6WC09	HD Heavy Duty package (350 lb weight capacity)									\$500	☐
SEAT WIDTH ³												
3	A6SW	Seat width:	14" ☐	15" ☐	16" ☐	17" ☐	18" ¹ ☐	19" ¹ ☐			NC	☐
4	A6SW	Seat width:		20" ¹ ☐	21" ² ☐	22" ² ☐				E2201	\$495	☐
SEAT DEPTH												
5	A6SD	Seat depth:	14" ☐	15" ☐	16" ☐	17" ☐	18" ☐	19" ☐			NC	☐
6	A6SD	Seat depth:		20" ☐						E2203	\$475	☐

1 - Available with HD Kit (350 lb)

2 - Only with HD Kit (350 lb)

[3. Width calculator](#)

A6 FRONT Seat-to-floor Height - Please select your choice

Fork/Position		3" caster			4" caster			5" caster			6" caster			7" caster			8" caster		
		Stem bolt Length			Stem bolt length			Stem bolt length			Stem bolt length			Stem bolt length			Stem bolt length		
		STD	+1"	+2"	STD	+1"	+2"	STD	+1"	+2"	STD	+1"	+2"	STD	+1"	+2"	STD	+1"	+2"
4"	P1	12	13	14	12 1/2	13 1/2	14 1/2												
	P2	12 1/2	13 1/2	14 1/2	13	14	15	13 1/2	14 1/2	15 1/2									
5"	P1																		
	P2				13	14	15	13 1/2	14 1/2	15 1/2									
	P3	13	14	15	13 1/2	14 1/2	15 1/2	14	15	16	14 1/2	15 1/2	16 1/2						
	P4	13 1/2	14 1/2	15 1/2	14	15	16	14 1/2	15 1/2	16 1/2	15	16	17	15 1/2	16 1/2	17 1/2			
7"	P1										15 1/2	16 1/2	17 1/2	16	17	18			
	P2							15 1/2	16 1/2	17 1/2	16	17	18	16 1/2	17 1/2	18 1/2			
	P3							16	17	18	16 1/2	17 1/2	18 1/2	17	18	19	17 1/2	18 1/2	19 1/2
	P4	15 1/2	16 1/2	17 1/2	16	17	18	16 1/2	17 1/2	18 1/2	17	18	19	17 1/2	18 1/2	19 1/2	18	19	20
	P5	16	17	18	16 1/2	17 1/2	18 1/2	17	18	19	17 1/2	18 1/2	19 1/2	18	19	20	18 1/2	19 1/2	20 1/2

Height available +/- 1/4" (see stem bolt and fork combination at the end of the order form)

CASTERS												
7	A6CA01	Composite wheel w/ PU tire (1" width)	3"□	4"□	5"□	6"□	7"□	8"□			NC	□
8	A6CA10	Pneumatic casters 6" x 1 ¼"								E2214/E2215	\$160	□
9	A6CA11	Pneumatic casters 8" x 1 ¼"								E2214/E2215	\$160	□
10	A6CA14	Newton UltraCaster <i>Soft-Roll</i> Composite	4" x 1 ½ "								\$80	□
11	A6CA15		5" x 1 ½ "								\$80	□
12	A6CA16		6" x 1 ½ "								\$80	□
13	A6CA17	Newton UltraCaster <i>Soft-Roll</i> Aluminum	4" x 1 ½ "								\$175	□
14	A6CA18		5" x 1 ½ "								\$185	□
15	A6CA19		6" x 1 ½ "								\$195	□
CASTER OPTIONS												
16	A6CA12	Caster Pin Locks (<i>Available with +1" or +2" Stem bolt ONLY</i>)								K0073	\$295	□
17	A6CA13	Frog Legs suspension forks phase-one (<i>Not available with caster pin locks</i>)	User weight :							E1015	\$475	□

Helio A6 - REAR Seat-to-floor Height - Please select your choice

20"	22"	24"	25"	26"
12	13	13 3/4	14 1/4	14 3/4
12 1/2	13 1/2	14 1/4	14 3/4	15 1/4
13	14	14 3/4	15 1/4	15 3/4
13 1/2	14 1/2	15 1/4	15 3/4	16 1/4
14	15	15 3/4	16 1/4	16 3/4
14 1/2	15 1/2	16 1/4	16 3/4	17 1/4
15	16	16 3/4	17 1/4	17 3/4
15 1/2	16 1/2	17 1/4	17 3/4	18 1/4
16	17	17 3/4	18 1/4	18 3/4
16 1/2	17 1/2	18 1/4	18 3/4	19 1/4
17	18	18 3/4	19 1/4	19 3/4
17 1/2	18 1/2	19 1/4	19 3/4	20 1/4

Configuration used: 1" seat slope with a center of gravity at 2"

Note: Tire size, center of gravity setting and seat slope may result in a variation of +/- ½" of selected seat to floor height"

SWING-IN/SWING-OUT FOOTREST											
18	A6FR03	60° Swing-in/swing-out footrest	See line 21 (5-13.5") ☐	13.5 - 17" ☐	14.5 - 19.5" ☐	15.5 - 20.5" ☐	Set at: *			NC	☐
19	A6FR04	70° Swing-in/swing-out footrest	See line 21 (5-13.5") ☐	13.5 - 17" ☐	14.5 - 19.5" ☐	15.5 - 20.5" ☐	Set at: *			NC	☐
20	A6FR15	80° Swing-in/swing-out footrest	See line 21 (5-13.5") ☐	13.5 - 16" ☐	14.5 - 19.5" ☐	15.5 - 20.5" ☐	Set at: *			NC	☐
		FOOTREST OPTIONS									
21	A6FR09		High mount (footplate attached to hanger with clamps)			* Range guide: 5 - 13.5"	Set at: *			\$75	☐

* +/- 1/2"

90° SWING-IN/SWING-OUT FOOTREST WITH FRONT MOUNT (Footplate mounted to hanger tube)									
22	A6FR06	90° swing-in/swing-out with front mount - Standard (no extension) ¹	7.5 - 11"	Set at: *			NC	☐	
23	A6FR07	90° swing-in/swing-out with front mount - Short extension 2" ^{1, 2}	9.5 - 13"	Set at: *			NC	☐	
24	A6FR08	90° swing-in/swing-out with front mount - Medium extension 4" ^{1, 2}	9.5 - 14.5"	Set at: *			NC	☐	
25	A6FR99	Omit Footrest					NC	☐	

1 - Available only with adjustable footplates

2 - May interfere with front casters

* +/- 1/2"

ELEVATING LEGREST (<i>unit</i>)										The elevating legrest replace the footrest side selected			
26	A6FR10	Elevating legrest with calf pad, 14" - 20" Length adjustable	Left ☐	Right ☐	Set at:		E0990	\$240	☐				

FOOTPLATES						
27	A6FP01	Standard (Composite)			STD	☐
28	A6FP02	Newton adjustable angle (Composite)		K0040	\$210	☐
29	A6FP03	Newton adjustable angle (Aluminum)		K0040	\$230	☐
30	A6FP06	One piece adjustable angle flip-up (Aluminum)		K0040	\$190	☐

FOOTPLATE ACCESSORIES						
31	A6FP08	Length adjustable heel loop (<i>not available with one piece footplate</i>)		E0951	\$40	☐
32	A6FP09	Calf strap		K0038	\$35	☐
33	A6FP30	Bodypoint Aeromesh Padded calf strap		K0038	\$115	☐

RESIDUAL LIMB SUPPORT <i>(height and angle adjustable) (unit)</i>									
34	A6FR12	Swing away residual limb support - 8" x 11"	Left ☐ Right ☐			E1020	\$375	☐	
35	A6FR23	Swing away residual limb support - 7" x 9"	Left ☐ Right ☐			E1020	\$375	☐	
Back canes with integrated push handle									
FIXED ANGLE BACK									
36	A6BC58	Straight aluminum back cane	14" ☐	16" ☐	18" ☐ 20" ☐			NC	☐
37	A6BC01	8° bend aluminum back cane	16" ☐	18" ☐	20" ☐			NC	☐
Back canes, height adjustable, with or without push handle									
38	A6BC11	Straight aluminum back cane	12 - 15" ☐	15 - 18" ☐	18 - 21" ☐ STD			NC	☐
39	A6BC12	8° bend aluminum back cane	13.5 - 15" ☐	15 - 18" ☐	18 - 21" ☐ STD			\$50	☐
40	A6BC13	Adjustable angle, straight aluminum back cane	15 - 18" ☐	18 - 21" ☐ STD				\$200	☐
41	A6BC14	Adjustable angle, 8° bend aluminum back cane	15 - 18" ☐	18 - 21" ☐ STD				\$250	☐
42		Select back cane angle:	85° ☐ ¹	90° STD ☐	95° ☐ 100° ☐	105° ☐	110° ☐		
Set height to: ²		(measured from seat sling to top of back upholstery)							
43	A6PH01	Push handles						STD	☐
44	A6PH05	Folding push handles ^{3, 5}						\$120	☐
45	A6PH06	Clamp on height adjustable push handle ^{3, 4, 5}	9 - 20" STD ☐	9-23" ☐				\$300	☐
46	A6PH02	Without push handles						NC	☐
<i>1 - Not compatible with 8 degree bent back 2 - If not indicated, default height setting is at the lowest position 3 - 265 pounds weight limit</i>									
<i>4 - Not compatible with swing away armrests , 8° bend back canes and ampute axle plates.</i>									
<i>5 - Not compatible with stroller handle, stabilizer bar, headrest and attendant lock</i>									
STROLLER HANDLE AND STABILIZER BAR ¹ <i>(only available with push handle)</i>									
47	A6PH04	Stroller handle (not available with the HD option)				K0108	\$160	☐	
48	A6PH03	Newton folding stabilizer bar					\$250	☐	
<i>1 - Stroller handle and stabilizer bar can't be combined</i>									
BACK UPHOLSTERY									
49	A6BH05	Slip-on Nylon back upholstery					STD	☐	
50	A6BH06	Adjustable tension Nylon back upholstery				E2611	\$300	☐	
51	A6BH99	Omit Back upholstery					NC	☐	
SEAT CUSHION									
52	A6SC01	2" Cushion					\$165	☐	
53	A6SC02	3" Cushion					\$185	☐	
54	A6SC99	Omit seat cushion					NC	☐	
REAR WHEELS									
55	A6RW01	Mag ¹	20" ☐	22" ☐	24" ☐		\$50	☐	
56	A6RW02	Newton One - spoke wheel ²	20" ☐	22" ☐	24" ☐		NC	☐	
57	A6RW03	Newton Gravity - ultralight wheel	20" ☐	22" ☐	24" ☐ 25" ☐ 26" ☐		\$300	☐	
58	A6RW04	Spinergy Spox	22" ☐	24" ☐	25" ☐ 26" ☐		\$795	☐	
59		Spinergy Spox spokes colors ³	Blk ☐ STD	Wht ☐ Blue ☐ Red ☐ Yel ☐			NC	☐	
60	A6RW05	Spinergy LX	22" ☐	24" ☐	25" ☐ 26" ☐		\$995	☐	
61		Spinergy LX spokes colors ³	Blk ☐ STD	Pur ☐ Wht ☐ Blue ☐ Red ☐ Yel ☐ Pink ☐ Grn ☐ Org ☐			NC	☐	
62	A6RW07	Spinergy CLX Carbon Blade	24" ☐	25" ☐			\$1 480	☐	
63		Spinergy CLX spokes colors ³	Blk ☐ STD	Pur ☐ Wht ☐ Blue ☐ Red ☐ Yel ☐ Pink ☐ Grn ☐ Org ☐			NC	☐	

TIRES			Available				
64	A6TI01	Soft Urethane 1 1/4"	20", 22" 24"				STD ☐
65	A6TI02	Soft Urethane knobby 1 3/8 "	22", 24", 25", 26"				\$95 ☐
66	A6TI04	Soft Urethane - Shox 1"	22", 24", 25", 26"				\$115 ☐
67	A6TI05	Pneumatic 1 3/8 "	20", 22", 24", 25", 26"		E2211/E2212		NC ☐
68	A6TI08	Pneumatic with Airless insert 1 3/8"	22", 24"		E2211/E2213		\$80 ☐
69	A6TI06	Pneumatic HP - 1"	20", 22", 24", 25", 26"		E2211/E2212		\$95 ☐
70	A6TI07	Pneumatic HP - Schwalbe Marathon Plus 1"	24", 25", 26"		E2211/E2212		\$125 ☐

HANDRIM								
71	A6HR01	Aluminum anodized					STD	☐
72	A6HR07	Newton Air Grip ¹ Black ☐					\$95	☐
73	A6HR0	Newton Air Grip ¹ (08) Sunkissed Orange ☐ (09) Burgundy ☐ (10) Sapphire Blue ☐					\$240	☐
		(11) Monster Green ☐ (12) Charcoal ☐ (13) Fuchsia ☐ (14) Steel Blue ☐						
		(15) Ferrari Red ☐ (16) Acid Green ☐ (17) White ☐ (18) Sky Blue ☐						
		Assembly position for aluminum and Newton Air Grip handrim						
74	A6HR29	Standard Handrim position Full tab					STD	☐
75	A6HR30	Narrow Handrim position <i>(not available in 22")</i> Short tab 2 positions					\$25	☐
76	A6HR31	Super Narrow Handrim <i>(not available in 22")</i> Short tab 2 positions STD ☐ Short tab 1 position ☐					\$25	☐

HANDRIM						
77	A6HR02	Plastic coated ³			\$95	☐
78	A6HR03	Q-Grip High friction Neoprene coated ^{1, 2, 3}			\$350	☐
79	A6HR28	Surge ^{1, 2, 3} (1-1/2" Oval handrim w/ Gripton high friction strip)			\$375	☐
80	A6HR05	Surge LT ^{1, 2, 3} (1-1/8" Oval handrim w/ Gripton high friction strip)			\$405	☐
81	A6HR04	Natural Fit Standard Grip ^{1,2,3} (1-1/2" Oval handrim w/standard thumb piece)			\$395	☐
82	A6HR25	Natural Fit Super Grip ^{1,2,3} (1-1/2" Oval handrim w/super grip thumb piece)			\$450	☐
83	A6HR99	Omit handrim			NC	☐

1- Not available with Mag wheels 2 - Not available in 20" 3- Only standard handrim position available


ONE ARM DRIVE									
84	A6HR22	One arm drive ⁴	24"	left ☐	right ☐	New		\$1 395	☐

4 - Only available with Newton One wheel, aluminum or plastic handrims

AXLE						
85	A6AX01	Quick release axle			STD	☐
86	A6AX02	Permanent axle			NC	☐
87	A6AX03	Quad release axle		K0108	\$125	☐
AXLE PLATE						
88	A6AX04	Amputee adjustable axle plate <i>(may interfere with swing away armrest)</i>		E0959	\$195	☐

REAR WHEEL CAMBER							
89	A6CM00	0° camber				STD	☐
90	A6CM03	3° camber				NC	☐

CENTER OF GRAVITY



Choose position of rear wheel axle plate:

Measure distance from front of back post to center of axle receiver.

Default setting is at 1" 1/4"

Most stable						Most responsive
AMPUTÉ	3/4"	1" 1/4	1" 3/4	2" 1/4	2" 3/4	
☐	☐	☐	☐	☐	☐	

If a configuration issue arises due to the choice of the Center Of Gravity (COG) Motion Composites reserves the right to reduce or increase the COG to the closest position.

WHEEL LOCK						
91	A6WL08	Newton Wheel lock - push to lock			STD	☐
92	A6WL02	Pull to lock			NC	☐
93	A6WL04	Newton grade aid push to lock <i>(suggested with pneumatic or high pressure tire)</i>		E0974	\$150	☐
94	A6WL13 A6WL14	6" removable extension handle		E0961	\$75	☐
95	A6WL05	Scissor Lock			\$240	☐
96	A6WL06	Attendant Lock (additional set of wheel lock mounted to the back)			\$200	☐
97	A6WL99	Omit wheel locks			NC	☐

FLIP BACK ARMREST <i>(Convertible to "T" armrest)</i> ¹		Height Adjustment				
98	A6AR01	10" Desk length arm pad	8 - 12" STD ☐	10 - 14" ☐	Locked ☐ Unlocked ☐ <i>STD</i>	E0973 \$365 ☐
99	A6AR02	14" Full length arm pad	8 - 12" STD ☐	10 - 14" ☐	Locked ☐ Unlocked ☐ <i>STD</i>	E0973 \$365 ☐
"T" ARMREST ¹		Height Adjustment				
100	A6AR13	10" Desk length arm pad	8 - 12" STD ☐	10 - 14" ☐		E0973 \$325 ☐
101	A6AR14	14" Full length arm pad	8 - 12" STD ☐	10 - 14" ☐		E0973 \$325 ☐

1 - May have a 1/2 " adjustment restriction on both limits of the range

SWING AWAY ARMREST ²			Height Adjustment				
102	A6AR05	Tubular arm pad	8 - 11"		E0973	\$225	☐
103	A6AR06	10" Desk length arm pad	8 - 11"		E0973	\$225	☐
		GEL ARM PAD					
104	A6AR12	14" Full Length Arm Pad - Gel Ovation (pair)				\$150	☐

2 - Could interfere with some rigid back rests

105	A6AR99	Omit armrest			NC	☐
-----	--------	--------------	--	--	----	--------------------------

SIDE GUARDS <i>(not available with Flip Back and "T" Armrest)</i>						
106	A6SG01	Plastic side guards <i>(removable)</i>			\$195	☐
107	A6SG02	Carbon fiber side guards <i>(fixed)</i> ¹			\$250	☐
108	A6SG03	Carbon fiber side guards <i>(removable)</i>			\$290	☐
109	A6SG06	Plastic fender side guards <i>(removable)</i> ²			\$270	☐
110	A6SG04	Carbon fiber fender side guards <i>(removable)</i>			\$540	☐
111	A6SG99	Omit side guard			NC	☐

1- Not available with back upholstery

2- Not compatible with 20" and 22" rear wheels

POSITIONING BELT								
112	A6SB01	Velcro adjustable	Standard 35" ☐	Long 60" ☐		E0978	\$45	☐
113	A6SB02	Auto buckle	Standard 35" ☐	Long 60" ☐		E0978	\$85	☐
114	A6SB07	Hardware attachment for 2 point seatbelt (1 pair)					NC	☐
115	A6SB03	Additional hardware for 2 extra attachment points (needed for the 4 point seat belt) (1 pair)					\$15	☐
116	A6SB04	BodyPoint Hip Belt - 2 point not padded	Medium 56" ☐	Long 61" ☐		E0978	\$80	☐
117	A6SB05	BodyPoint Hip Belt - 2 point padded	Small 51" ☐	Medium 56" ☐	Long 61" ☐	E0978	\$170	☐
118	A6SB06	BodyPoint Hip Belt - 4 point padded	Small 51" ☐	Medium 56" ☐	Long 61" ☐	E0978	\$190	☐
119	A6SB08	BodyPoint EVOFLEX® Positioning Belt Rehab latch and Band Clamps	Medium 18" ☐	Long 24" ☐		E0978	\$285	☐
120	A6SB99	Omit positioning belt					STD	☐

HEADREST									
121	A6AC01	Profiled Headrest	Small (3.5"X6") ☐	Medium (4.5"X7") ☐	Large (5.5"X8") ☐		E0955	\$180	☐
122	A6AC02	Headrest mounting hardware with removable non folding stabilizer bar ¹					E1028	\$295	☐

1 - Only available with push handles - Not available with a stroller handle or a folding stabiliser bar

ACCESSORIES					
123	A6AC03	Anti tippers		E0971	\$140 ☐
124	A6AC27	Tip assist ¹	New		\$60 ☐
125	A6AC04	Transit Tie Down (<i>Unoccupied</i>)			NC ☐
126	A6AC06	Cane Carrier Left ☐ Right ☐		E2207	\$195 ☐
127	A6AC07	Oxygen tank holder - E Type Left		E2208	\$310 ☐
128	A6AC08	Spoke guards ² Available: 20", 22", 24", 26"		K0065	\$115 ☐
129	A6AC09	Footrest neoprene impact guard (pair) - Black		K0108	\$90 ☐
130	A6AC10	Seat pouch			\$75 ☐
131	A6AC26	BodyPoint Mobility Bag			\$125 ☐
132	A6AC23	Newton Backpack			\$120 ☐
133	A6AC12	Velcro Closing Strap			\$50 ☐
134	A6AC13	Tool kit			\$45 ☐
135	A6AC16	Motion Composites Polo Shirt W/S ☐ W/M ☐ W/L ☐ M/S ☐ M/M ☐ M/L ☐ M/XL ☐ M/XXL ☐ M/XXXL ☐ M/XXXXL ☐			\$32 ☐
136	A6AC15	Military Patches Army ☐ Navy ☐ Air Force ☐ Marines ☐ Coast Guard ☐			\$90 ☐
137	A6AC15	Combat Wounded Purple Heart Patch			\$90 ☐

1 - Cannot be combined with Anti Tippers

2 - Not available with Mag wheels

FRAME COLOR (with black hydroformed aluminum crossbrace)					
138	A6FC01	Black ☐			NC ☐
139	A6FC03	Sunkissed Orange ☐ Burgundy ☐ Sapphire Blue ☐ Monster Green ☐ Charcoal ☐ Fuchsia ☐ Steel Blue ☐ Ferrari Red ☐ Acid Green ☐			NC ☐
140	A6FC99	Custom Color Specify PPG Paintit color:			\$360 ☐

TOTAL	
-------	--

Fork and Stem bolt combinations

4"

STD

+1"

+2"

5"

STD

+1"

+2"

7"

STD

+1"

+2"

Notes :