



Helio A6

ORDER FORM - USA

P.O. number:* _____ Date: _____

Company Name:* _____ Phone:* _____

Ordered by:* _____

Clinician/ATP/RTS: _____

Email: _____

Marked for: _____

Ship to:* _____

Name: _____

Address: _____

City: _____ State: _____

Zip Code: _____

Phone: _____

* Required

A6MA082
English Manual USA

A6MA089
Spanish Manual USA

Taxes not included

MODEL									USD			
1	A6WC19	Ultralight folding Helio A6 Wheelchair - 6000 series aluminum (265lb weight capacity)								K0005	\$2 095	☐
2	A6WC09	HD Heavy Duty package (350 lb weight capacity)									\$500	☐
SEAT WIDTH ³												
3	A6SW	Seat width:	14" ☐	15" ☐	16" ☐	17" ☐	18" ¹ ☐	19" ¹ ☐			NC	☐
4	A6SW	Seat width:		20" ¹ ☐	21" ² ☐	22" ² ☐				E2201	\$495	☐
SEAT DEPTH												
5	A6SD	Seat depth:	14" ☐	15" ☐	16" ☐	17" ☐	18" ☐	19" ☐			NC	☐
6	A6SD	Seat depth:		20" ☐						E2203	\$475	☐

1 - Available with HD Kit (350 lb)

2 - Only with HD Kit (350 lb)

[3. Width calculator](#)

A6 FRONT Seat-to-floor Height - Please select your choice

Fork/Position		3" caster			4" caster			5" caster			6" caster			7" caster			8" caster		
		Stem bolt Length			Stem bolt length			Stem bolt length			Stem bolt length			Stem bolt length			Stem bolt length		
		STD	+1"	+2"	STD	+1"	+2"	STD	+1"	+2"	STD	+1"	+2"	STD	+1"	+2"	STD	+1"	+2"
4"	P1	12	13	14	12 1/2	13 1/2	14 1/2												
	P2	12 1/2	13 1/2	14 1/2	13	14	15	13 1/2	14 1/2	15 1/2									
5"	P1																		
	P2				13	14	15	13 1/2	14 1/2	15 1/2									
	P3	13	14	15	13 1/2	14 1/2	15 1/2	14	15	16	14 1/2	15 1/2	16 1/2						
	P4	13 1/2	14 1/2	15 1/2	14	15	16	14 1/2	15 1/2	16 1/2	15	16	17	15 1/2	16 1/2	17 1/2			
7"	P1										15 1/2	16 1/2	17 1/2	16	17	18			
	P2							15 1/2	16 1/2	17 1/2	16	17	18	16 1/2	17 1/2	18 1/2			
	P3							16	17	18	16 1/2	17 1/2	18 1/2	17	18	19	17 1/2	18 1/2	19 1/2
	P4	15 1/2	16 1/2	17 1/2	16	17	18	16 1/2	17 1/2	18 1/2	17	18	19	17 1/2	18 1/2	19 1/2	18	19	20
	P5	16	17	18	16 1/2	17 1/2	18 1/2	17	18	19	17 1/2	18 1/2	19 1/2	18	19	20	18 1/2	19 1/2	20 1/2

Height available +/- 1/4" (see stem bolt and fork combination at the end of the order form)

CASTERS												
7	A6CA01	Composite wheel w/ PU tire (1" width)	3" ☐	4" ☐	5" ☐	6" ☐	7" ☐	8" ☐			NC	☐
8	A6CA10	Pneumatic casters 6" x 1 ¼"								E2214/E2215	\$160	☐
9	A6CA11	Pneumatic casters 8" x 1 ¼"								E2214/E2215	\$160	☐
10	A6CA14	Newton UltraCaster <i>Soft-Roll</i> Composite	4" x 1 ½ "								\$80	☐
11	A6CA15		5" x 1 ½ "								\$80	☐
12	A6CA16		6" x 1 ½ "								\$80	☐
13	A6CA17	Newton UltraCaster <i>Soft-Roll</i> Aluminum	4" x 1 ½ "								\$175	☐
14	A6CA18		5" x 1 ½ "								\$185	☐
15	A6CA19		6" x 1 ½ "								\$195	☐
CASTER OPTIONS												
16	A6CA12	Caster Pin Locks (<i>Available with +1" or +2" Stem bolt ONLY</i>)								K0073	\$295	☐
17	A6CA13	Carbon fiber Frog Legs suspension forks (<i>Not available with caster pin locks</i>)								E1015	\$475	☐

Helio A6 - REAR Seat-to-floor Height - Please select your choice

20"	22"	24"	25"	26"
12	13	13 3/4	14 1/4	14 3/4
12 1/2	13 1/2	14 1/4	14 3/4	15 1/4
13	14	14 3/4	15 1/4	15 3/4
13 1/2	14 1/2	15 1/4	15 3/4	16 1/4
14	15	15 3/4	16 1/4	16 3/4
14 1/2	15 1/2	16 1/4	16 3/4	17 1/4
15	16	16 3/4	17 1/4	17 3/4
15 1/2	16 1/2	17 1/4	17 3/4	18 1/4
16	17	17 3/4	18 1/4	18 3/4
16 1/2	17 1/2	18 1/4	18 3/4	19 1/4
17	18	18 3/4	19 1/4	19 3/4
17 1/2	18 1/2	19 1/4	19 3/4	20 1/4

Configuration used: 1" seat slope with a center of gravity at 2"

Note: Tire size, center of gravity setting and seat slope may result in a variation of +/- ½" of selected seat to floor height"

SWING-IN/SWING-OUT FOOTREST									
18	A6FR03	60° Swing-in/swing-out footrest ¹	13.5 - 17" □	14.5 - 19.5" □	15.5 - 20.5" □	Set at: *		NC	□
19	A6FR04	70° Swing-in/swing-out footrest ¹	13.5 - 17" □	14.5 - 19.5" □	15.5 - 20.5" □	Set at: *		NC	□
20	A6FR15	80° Swing-in/swing-out footrest ¹	13.5 - 16" □	14.5 - 19.5" □	15.5 - 20.5" □	Set at: *		NC	□
		FOOTREST OPTIONS							
21	A6FR09		High mount <i>(footplate attached to the hanger with clamps)</i>		5 - 13.5" □	Set at: *		\$75	□

1 - For a 5.13" range, see line 21

* +/- 1/2"

90° SWING-IN/SWING-OUT FOOTREST WITH FRONT MOUNT <i>(Footplate mounted to hanger tube)</i>								
22	A6FR06	90° swing-in/swing-out with front mount - Standard (no extension) ²	7.5 - 11"	Set at: *			NC	☐
23	A6FR07	90° swing-in/swing-out with front mount - Short extension 2" ^{2, 3}	9.5 - 13"	Set at: *			NC	☐
24	A6FR08	90° swing-in/swing-out with front mount - Medium extension 4" ^{2, 3}	9.5 - 14.5"	Set at: *			NC	☐
25	A6FR99	Omit Footrest					NC	☐

2 - Available only with adjustable footplates

3 - May interfere with front casters

* +/- 1/2"

ELEVATING LEGREST <i>(unit)</i>			The elevating legrest replace the footrest side selected					
26	A6FR10	Elevating legrest with calf pad, 14" - 19" Length adjustable	Left ☐	Right ☐	Set at:	E0990	\$240	☐

FOOTPLATES						
27	A6FP01	Standard (Composite)			STD	☐
28	A6FP02	Newton adjustable angle (Composite)		K0040	\$210	☐
29	A6FP03	Newton adjustable angle (Aluminum)		K0040	\$230	☐
30	A6FP06	One piece adjustable angle flip-up (Aluminum)		K0040	\$190	☐

FOOTPLATE ACCESSORIES						
31	A6FP08	Length adjustable heel loop (<i>not available with one piece footplate</i>)		E0951	\$40	☐
32	A6FP09	Calf strap		K0038	\$35	☐
33	A6FP30	Bodypoint Aeromesh Padded calf strap		K0038	\$115	☐

RESIDUAL LIMB SUPPORT <i>(height and angle adjustable) (unit)</i>								
34	A6FR12	Swing away residual limb support - 8" x 11"	Left ☐	Right ☐		E1020	\$375	☐
35	A6FR23	Swing away residual limb support - 7" x 9"	Left ☐	Right ☐		E1020	\$375	☐
Back canes with integrated push handle								
FIXED ANGLE BACK								
36	A6BC58	Straight aluminum back cane	14" ☐	16" ☐	18" ☐	20" ☐	NC	☐
37	A6BC01	8° bend aluminum back cane	16" ☐	18" ☐	20" ☐		NC	☐
Back canes, height adjustable, with or without push handle								
FIXED ANGLE BACK								
38	A6BC11	Straight aluminum back cane	12 - 15" ☐	15 - 18" ☐	18 - 21" ☐	STD	NC	☐
39	A6BC12	8° bend aluminum back cane	13.5 - 15" ☐	15 - 18" ☐	18 - 21" ☐	STD	\$50	☐
ADJUSTABLE ANGLE BACK								
40	A6BC13	Adjustable angle, straight aluminum back cane	15 - 18" ☐	18 - 21" ☐	STD		\$200	☐
41	A6BC14	Adjustable angle, 8° bend aluminum back cane	15 - 18" ☐	18 - 21" ☐	STD		\$250	☐
42		Select back cane angle:	85° ¹ ☐	90° STD ☐	95° ☐	100° ☐	105° ☐	110° ☐
Set height to: ²						(measured from seat sling to top of back upholstery)		
PUSH HANDLE OPTION								
43	A6PH01	Push handles					STD	☐
44	A6PH05	Folding push handles ^{3,5}					\$120	☐
45	A6PH06	Clamp on height adjustable push handle ^{3,4,5}	9 - 20" STD ☐	9-23" ☐			\$300	☐
46	A6PH02	Without push handles					NC	☐
1 - Not compatible with 8 degree bent back 2 - If not indicated, default height setting is at the lowest position 3 - 265 pounds weight limit 4 - Not compatible with swing away armrests , 8° bend back canes and amputee axle plates. 5 - Not compatible with stroller handle, stabilizer bar, headrest and attendant lock								
STROLLER HANDLE AND STABILIZER BAR ¹ <i>(only available with push handle)</i>								
47	A6PH04	Stroller handle (not available with the HD option)				K0108	\$160	☐
48	A6PH03	Newton folding stabilizer bar					\$250	☐
1 - Stroller handle and stabilizer bar can't be combined								
BACK UPHOLSTERY								
49	A6BH05	Slip-on Nylon back upholstery					STD	☐
50	A6BH06	Adjustable tension Nylon back upholstery				Product improved	E2611	\$300 ☐
51	A6BH99	Omit Back upholstery					NC	☐
SEAT CUSHION								
52	A6SC01	2" Cushion					\$165	☐
53	A6SC02	3" Cushion					\$185	☐
54	A6SC99	Omit seat cushion					NC	☐
ONE ARM DRIVE ¹ <i>Only available with Newton One rear wheels with Aluminum, Plastic coated Handrims</i>								
55	A6HR22	One arm drive ^{2,3} 24"	Left ☐	Right ☐			1 395 \$	☐
56	A6HR98	Omit handrim on opposite side of the One Arm Drive (OAD)					NC	☐
1- Offered with handrims on each side of the wheelchair 2- Only standard handrim position available 3 - Not available with HD option								
REAR WHEELS								
57	A6RW01	Mag ¹	20" ☐	22" ☐	24" ☐		\$50	☐
58	A6RW02	Newton One - spoke wheel ²	20" ☐	22" ☐	24" ☐		NC	☐
59	A6RW03	Newton Gravity - ultralight wheel	20" ☐	22" ☐	24" ☐	25" ☐	26" ☐	\$300 ☐
60	A6RW04	Spinergy Spox	22" ☐	24" ☐	25" ☐	26" ☐	\$795	☐
61		Spinergy Spox spokes colors ³	Blk STD ☐	Wht ☐	Blue ☐	Red ☐	Yel ☐	NC ☐
62	A6RW05	Spinergy LX	22" ☐	24" ☐	25" ☐	26" ☐	\$995	☐
63		Spinergy LX spokes colors ³	Blk STD ☐	Pur ☐	Wht ☐	Blue ☐	Red ☐	Yel ☐
64	A6RW07	Spinergy CLX Carbon Blade	24" ☐	25" ☐				\$1 480 ☐
65		Spinergy CLX spokes colors ³	Blk STD ☐	Pur ☐	Wht ☐	Blue ☐	Red ☐	Yel ☐
1 - Not compatible with HP tires ; Only available with aluminum or plastic coated handrim in standard position 2 - Not compatible with Shox tires 3 - Colour spokes require an additionnal 10 days lead-time								

TIRES			Available				
66	A6TI01	Soft Urethane 1 1/4"	20", 22" 24"			STD	☐
67	A6TI02	Soft Urethane knobby 1 3/8 "	22", 24", 25", 26"			\$95	☐
68	A6TI04	Soft Urethane - <i>Shox 1"</i>	22", 24", 25", 26"			\$115	☐
69	A6TI05	Pneumatic 1 3/8 "	20", 22", 24", 25", 26"		E2211/E2212	NC	☐
70	A6TI08	Pneumatic with Airless insert 1 3/8"	22", 24"		E2211/E2213	\$80	☐
71	A6TI06	Pneumatic HP - 1"	20", 22", 24", 25", 26"		E2211/E2212	\$95	☐
72	A6TI07	Pneumatic HP - <i>Schwalbe Marathon Plus 1"</i>	24", 25", 26"		E2211/E2212	\$125	☐


HANDRIM							
73	A6HR01	Aluminum anodized				STD	☐
74	A6HR07	Newton Air Grip ¹ Black ☐			WAMO	\$150	☐
75	A6HR08	Newton Air Grip ¹ (08) Sunkissed Orange ☐ (09) Burgundy ☐ (10) Sapphire Blue ☐			WAMO	\$330	☐
		(11) Monster Green ☐ (12) Charcoal ☐ (13) Fuchsia ☐ (14) Steel Blue ☐					
		(15) Ferrari Red ☐ (16) Acid Green ☐ (17) White ☐ (18) Sky Blue ☐					
		Assembly position for aluminum and Newton Air Grip handrim					
76	A6HR29	Standard Handrim position Full tab				STD	☐
77	A6HR30	Narrow Handrim position Short tab 2 positions				\$25	☐
78	A6HR31	Super Narrow Handrim Short tab 2 positions STD ☐ Short tab 1 position ☐				\$25	☐

HANDRIM						
79	A6HR02	Plastic coated ³		WAMO	\$160	☐
80	A6HR03	Q-Grip High friction Neoprene coated ^{1, 2, 3}		WAMO	\$430	☐
81	A6HR28	Surge ^{1, 2, 3} (1-1/2" Oval handrim w/ Gripton high friction strip)			\$470	☐
82	A6HR05	Surge LT ^{1, 2, 3} (1-1/8" Oval handrim w/ Gripton high friction strip)			\$500	☐
83	A6HR04	Natural Fit ^{1, 2, 3}			\$310	☐
84	A6HR99	Omit handrim			NC	☐
Thumb piece ⁴ selection for Natural Fit handrims						
85	A6HR23	Standard grip			\$85	☐
86	A6HR27	Super grip			\$140	☐
87	A6HR099	Omit thumb piece			NC	☐

1 - Not available with Mag wheels or OAD 2 - Not available in 20" 3- Only standard handrim position available 4- Compatible with Surge, Surge LT handrims

AXLE						
88	A6AX01	Quick release axle			STD	☐
89	A6AX02	Permanent axle			NC	☐
90	A6AX03	Quad release axle		K0108	\$125	☐
AXLE PLATE						
91	A6AX04	Amputee adjustable axle plate (may interfere with swing away armrest)		E0959	\$195	☐

REAR WHEEL CAMBER						
92	A6CM00	0° camber			STD	☐
93	A6CM03	3° camber			NC	☐

CENTER OF GRAVITY																							
																							
Choose position of rear wheel axle plate:																							
Measure distance from front of back post to center of axle receiver.																							
Default setting is at 1" 1/4"																							
<table border="1"> <thead> <tr> <th>Most stable</th> <th></th> <th></th> <th></th> <th></th> <th>Most responsive</th> </tr> </thead> <tbody> <tr> <td>AMPUTÉ</td> <td>3/4"</td> <td>1" 1/4</td> <td>1" 3/4</td> <td>2" 1/4</td> <td>2" 3/4</td> </tr> <tr> <td>☐</td> <td>☐</td> <td>☐</td> <td>☐</td> <td>☐</td> <td>☐</td> </tr> </tbody> </table>						Most stable					Most responsive	AMPUTÉ	3/4"	1" 1/4	1" 3/4	2" 1/4	2" 3/4	☐	☐	☐	☐	☐	☐
Most stable					Most responsive																		
AMPUTÉ	3/4"	1" 1/4	1" 3/4	2" 1/4	2" 3/4																		
☐	☐	☐	☐	☐	☐																		
If a configuration issue arises due to the choice of the Center Of Gravity (COG) Motion Composites reserves the right to reduce or increase the COG to the closest position.																							

WHEEL LOCK						
94	A6WL08	Newton Wheel lock - push to lock			STD	☐
95	A6WL02	Pull to lock			NC	☐
96	A6WL04	Newton grade aid push to lock <i>(suggested with pneumatic or high pressure tire)</i>		E0974	\$150	☐
97	A6WL13 A6WL14	6" removable extension handle		E0961	\$75	☐
98	A6WL05	Scissor Lock			\$140	☐
99	A6WL06	Attendant Lock (additional set of wheel lock mounted to the back)			\$200	☐
100	A6WL99	Omit wheel locks			NC	☐

FLIP BACK ARMREST <i>(Convertible to "T" armrest)</i> ¹			Height Adjustment						
101	A6AR01	10" Desk length arm pad	8 - 12" STD ☐	10 - 14" ☐	Locked ☐	Unlocked ☐ STD	E0973	\$365	☐
102	A6AR02	14" Full length arm pad	8 - 12" STD ☐	10 - 14" ☐	Locked ☐	Unlocked ☐ STD	E0973	\$365	☐
"T" ARMREST ¹			Height Adjustment						
103	A6AR13	10" Desk length arm pad	8 - 12" STD ☐	10 - 14" ☐			E0973	\$325	☐
104	A6AR14	14" Full length arm pad	8 - 12" STD ☐	10 - 14" ☐			E0973	\$325	☐

1 - May have a 1/2 " adjustment restriction on both limits of the range

SWING AWAY ARMREST ²			Height Adjustment						
105	A6AR05	Tubular arm pad	8 - 11"				E0973	\$225	☐
106	A6AR06	10" Desk length arm pad	8 - 11"				E0973	\$225	☐
		GEL ARM PAD							
107	A6AR12	14" Full Length Arm Pad - Gel Ovation (pair)						\$150	☐
108	A6AR99	Omit armrest						NC	☐

2 - Could interfere with some rigid back rests

SIDE GUARDS <i>(not available with Flip Back and "T" Armrest)</i>									
109	A6SG01	Plastic side guards <i>(removable)</i>						\$195	☐
110	A6SG02	Carbon fiber side guards <i>(fixed)</i> ¹						\$250	☐
111	A6SG03	Carbon fiber side guards <i>(removable)</i>						\$290	☐
112	A6SG06	Plastic fender side guards <i>(removable)</i> ²						\$270	☐
113	A6SG04	Carbon fiber fender side guards <i>(removable)</i>						\$540	☐
114	A6SG99	Omit side guard						NC	☐

1- Not available with back upholstery

2- Not compatible with 20" and 22" rear wheels

POSITIONING BELT									
115	A6SB01	Velcro adjustable	Standard 35" ☐	Long 60" ☐			E0978	\$45	☐
116	A6SB02	Auto buckle	Standard 35" ☐	Long 60" ☐			E0978	\$85	☐
117	A6SB07	Hardware attachment for 2 point seatbelt (1 pair)						NC	☐
118	A6SB03	Additional hardware for 2 extra attachment points (needed for the 4 point seat belt) (1 pair)						\$15	☐
119	A6SB04	BodyPoint Hip Belt - 2 point not padded		Medium 56" ☐	Long 61" ☐		E0978	\$80	☐
120	A6SB05	BodyPoint Hip Belt - 2 point padded	Small 51" ☐	Medium 56" ☐	Long 61" ☐		E0978	\$170	☐
121	A6SB06	BodyPoint Hip Belt - 4 point padded	Small 51" ☐	Medium 56" ☐	Long 61" ☐		E0978	\$190	☐
122	A6SB08	BodyPoint EVOFLEX® Positioning Belt Rehab latch and Band Clamps		Medium 18" ☐	Long 24" ☐		E0978	\$285	☐
123	A6SB99	Omit positioning belt						STD	☐

HEADREST									
124	A6AC01	Profiled Headrest	Small (3.5"x6") ☐	Medium (4.5"x7") ☐	Large (5.5"x8") ☐		E0955	\$180	☐
125	A6AC28	Headrest mounting hardware with removable non folding stabilizer bar ¹					E1028	\$295	☐

1 - Only available with push handles - Not available with a stroller handle or a folding stabiliser bar

ACCESSORIES					
126	A6AC03	Anti tippers			E0971 \$140 ☐
127	A6AC27	Tip assist ¹			\$60 ☐
128	A6AC04	Transit Tie Down (<i>Unoccupied</i>)			NC ☐
129	A6AC05	WC-19 Transit Tie Down w/Q'straint belt ² Short belt 20-41" ☐ Long belt 20-46" ☐			\$320 ☐
130	A6AC06	Cane Carrier Left ☐ Right ☐			E2207 \$195 ☐
131	A6AC07	Oxygen tank holder - E Type Left			E2208 \$310 ☐
132	A6AC08	Spoke guards ³ Available: 20", 22", 24", 26"			K0065 \$115 ☐
133	A6AC09	Footrest neoprene impact guard (pair) - Black			K0108 \$90 ☐
134	A6AC10	Seat pouch			\$75 ☐
135	A6AC26	BodyPoint Mobility Bag			\$125 ☐
136	A6AC23	Newton Backpack			\$120 ☐
137	A6AC12	Velcro Closing Strap			\$50 ☐
138	A6AC13	Tool kit			\$45 ☐
139	A6AC16	Motion Composites Polo Shirt W/S ☐ W/M ☐ W/L ☐ M/S ☐ M/M ☐ M/L ☐ M/XL ☐ M/XXL ☐ M/XXXL ☐ M/XXXXL ☐			\$32 ☐
140	A6AC15	Military Patches Army ☐ Navy ☐ Air Force ☐ Marines ☐ Coast Guard ☐			\$90 ☐
141	A6AC15	Combat Wounded Purple Heart Patch			\$90 ☐

1 - Cannot be combined with Anti Tippers 2 - May interfere with swing away armrest and Clamp on height adjustable push handle; Not available with HD option

3 - Not available with Mag wheels and Spinergy

FRAME COLOR (with black hydroformed aluminum crossbrace)					
142	A6FC01	Black ☐			NC ☐
143	A6FC03	Sunkissed Orange ☐ Burgundy ☐ Sapphire Blue ☐ Monster Green ☐ Charcoal ☐ Fuchsia ☐ Steel Blue ☐ Ferrari Red ☐ Acid Green ☐			NC ☐
144	A6FC99	Custom Color Specify PPG Paintit color:			additional 10 days lead-time \$360 ☐

TOTAL	
-------	--

Fork and Stem bolt combinations

4"

STD

+1"

+2"

5"

STD

+1"

+2"

7"

STD

+1"

+2"

Notes :