

**HELIO****ORDER FORM - USA**

Account: _____ PO: _____
 Vendor: _____
 Date: _____
 Ordered by: _____
 Email: _____
 Marked for: _____

Ship to: _____
 Name: _____
 Address: _____
 City: _____ State: _____
 Zip Code: _____
 Phone: _____

MODEL		USD	
1	Ultralight folding HELIO wheelchair (250 lbs. weight capacity) A7 package: sideframes, symmetrical hydroformed crossbrace and back canes made of 7000 series aluminum	\$2 295	<input type="checkbox"/>
2	C2 carbon fiber upgrade: sideframes, symmetrical molded crossbrace and back canes made of ultralight carbon fiber	\$1 495	<input type="checkbox"/>
3	HD kit (350 lbs. weight capacity)	\$500	<input type="checkbox"/>
4	E-Motion adaptator Kit (4 adaptors)	\$196	<input type="checkbox"/>
FRAME COLOR			
5	Basalt Black/Carbon Fiber* <input type="checkbox"/> Sunkissed Orange <input type="checkbox"/> Burgundy <input type="checkbox"/> Sapphir Blue <input type="checkbox"/> Monster Green <input type="checkbox"/> Charcoal <input type="checkbox"/> Fushia <input type="checkbox"/> Steel Blue <input type="checkbox"/> Ferrari Red <input type="checkbox"/> Acid Green <input type="checkbox"/>	NC	<input type="checkbox"/>
6	Custom Color	\$360	<input type="checkbox"/>

* Basalt Black color replaced for Carbon fiber finish with the C2 upgrade

SEAT WIDTH			
6	Seat width: 14" <input type="checkbox"/> 15" <input type="checkbox"/> 16" <input type="checkbox"/> 17" <input type="checkbox"/> 18" <input type="checkbox"/> ¹ 19" <input type="checkbox"/> ¹	NC	<input type="checkbox"/>
7	Seat width: 20" <input type="checkbox"/> ¹ 21" <input type="checkbox"/> ² 22" <input type="checkbox"/> ²	\$495	<input type="checkbox"/>
SEAT DEPTH			
8	Seat depth: 14" <input type="checkbox"/> 15" <input type="checkbox"/> 16" <input type="checkbox"/> 17" <input type="checkbox"/> 18" <input type="checkbox"/> 19" <input type="checkbox"/>	NC	<input type="checkbox"/>
9	Seat depth: 20" <input type="checkbox"/>	\$475	<input type="checkbox"/>

1. Available with HD Kit (350 lb)

2. Only with HD Kit (350 lb)

NEWTON ULTRALIGHT PACKAGE			
10	Includes: Carbon fiber angle adjustable footplates + Newton Gravity ultralight rear wheels	Save \$205	\$480 <input type="checkbox"/>

STRAIGHT BACK			
11	Height adjustable back cane 12 - 15" <input type="checkbox"/> 15 - 18" <input type="checkbox"/> 18 - 21" <input type="checkbox"/>	NC	<input type="checkbox"/>
Set height to: (measured from seat sling to top of back upholstery / see drawing at the end of the order form)			
8° BEND BACK			
12	Height adjustable back cane 13 ½ - 15" <input type="checkbox"/> 15 - 18" <input type="checkbox"/> 18 - 21" <input type="checkbox"/>	\$50	<input type="checkbox"/>
Set height to: (measured from seat sling to top of back upholstery / see drawing at the end of the order form)			
ADJUSTABLE ANGLE BACK			
13	Height adjustable back cane straight 15 - 18" <input type="checkbox"/> 18 - 21" <input type="checkbox"/>	\$200	<input type="checkbox"/>
14	Height adjustable back cane 8° bend 15 - 18" <input type="checkbox"/> 18 - 21" <input type="checkbox"/>	\$250	<input type="checkbox"/>
15	Select back cane angle: 85° <input type="checkbox"/> 90° <input type="checkbox"/> 95° <input type="checkbox"/> 100° <input type="checkbox"/> 105° <input type="checkbox"/> 110° <input type="checkbox"/>		
Set height to: (measured from seat sling to top of back upholstery / see drawing at the end of the order form)			
PUSH HANDLES			
16	Push handles, 0 - 3" height adjustable	NC	<input type="checkbox"/>
17	Straight extensions, 0 - 3" height adjustable	NC	<input type="checkbox"/>
18	Stroller handle (only available with push handle)	\$160	<input type="checkbox"/>
19	Folding Stabilizer bar (only available with push handle)	\$250	<input type="checkbox"/>

SEATING CUSHION					
20	2" Cushion			\$165	<input type="checkbox"/>
21	3" Cushion			\$185	<input type="checkbox"/>
22	nXt FIT cushion series (See nXt FIT pricelist)	See pricelist			<input type="checkbox"/>
BACK UPHOLSTERY					
23	Slip on back upholstery			NC	<input type="checkbox"/>
24	Adjustable tension back upholstery			\$300	<input type="checkbox"/>
25	Omit Back upholstery			NC	<input type="checkbox"/>
RIGID BACK					
26	nXt back series (See nXt back pricelist)	See pricelist			<input type="checkbox"/>

FLIP BACK ARMREST <i>(Convertible to "T" armrest)</i>		Height Adjustment				
27	Desk Length	8 - 12" <input type="checkbox"/>	10 - 14" <input type="checkbox"/>		\$350	<input type="checkbox"/>
28	Full Length	8 - 12" <input type="checkbox"/>	10 - 14" <input type="checkbox"/>		\$350	<input type="checkbox"/>
"T" ARMREST		Height Adjustment				
29	Desk Length	8 - 12" <input type="checkbox"/>	10 - 14" <input type="checkbox"/>		\$325	<input type="checkbox"/>
30	Full Length	8 - 12" <input type="checkbox"/>	10 - 14" <input type="checkbox"/>		\$325	<input type="checkbox"/>
SWING AWAY ARMREST		Height Adjustment				
31	Tubular arm pad	8 - 11"			\$225	<input type="checkbox"/>
32	Desk length arm pad	8 - 11"			\$225	<input type="checkbox"/>

SIDE GUARDS					
33	Plastic side guards <i>(removable)</i>			\$195	<input type="checkbox"/>
34	Carbon fiber side guards <i>(fixed)</i>			\$250	<input type="checkbox"/>
35	Carbon fiber mud guards <i>(removable)</i>			\$540	<input type="checkbox"/>

SWING-IN/SWING-OUT FOOTREST					
Set length to:					
36	60° Swing-in/swing-out footrest	13.5 - 21"		NC	<input type="checkbox"/>
37	70° Swing-in/swing-out footrest	13.5 - 21"		NC	<input type="checkbox"/>
38	90° Swing-in/swing-out footrest	13 - 18"		NC	<input type="checkbox"/>
39	90° Sing-in/swing-out footrest front mount <i>(Footplate mounted to hanger tube)</i>	8"-14"		NC	<input type="checkbox"/>
40	Elevating legrest with calf pad	14 - 20"	Left <input type="checkbox"/>	Right <input type="checkbox"/>	\$240 <input type="checkbox"/>
41	Length modification outside of specified range	Set length to:		\$75	<input type="checkbox"/>

STUMP SUPPORT <i>(Height and angle adjustable)</i>					
42	Swing away stump support:	8" x 11"	Left <input type="checkbox"/>	Right <input type="checkbox"/>	\$375 <input type="checkbox"/>

FOOTPLATES					
43	Standard <i>(Composites)</i>			STD	<input type="checkbox"/>
44	Adjustable angle <i>(Aluminum)</i>			\$195	<input type="checkbox"/>
45	Adjustable angle <i>(Carbon fiber with rubber inserts)</i>			\$385	<input type="checkbox"/>
46	One piece adjustable angle flip-up <i>(Aluminum)</i>			\$190	<input type="checkbox"/>
47	One piece adjustable angle flip-up <i>(Carbon fiber with rubber inserts)</i>			\$435	<input type="checkbox"/>
48	Length adjustable heel loop <i>(not available with one piece footplate)</i>			\$40	<input type="checkbox"/>
49	Calf strap			\$35	<input type="checkbox"/>
50	Padded calf strap			\$90	<input type="checkbox"/>

CASTERS										
51	Composite wheel w/ PU tire (1" width)	3" <input type="checkbox"/>	4" <input type="checkbox"/>	5" <input type="checkbox"/>	6" <input type="checkbox"/>	7" <input type="checkbox"/>	8" <input type="checkbox"/>			NC <input type="checkbox"/>
52	Pneumatic	6" x 1 1/4" <input type="checkbox"/>			8" x 1 1/4" <input type="checkbox"/>					\$160 <input type="checkbox"/>
FROG LEGS CASTERS										
53		4" x 1 1/2" "								\$80 <input type="checkbox"/>
54	Composite wheel Soft-Roll® by Frog Legs ¹	5" x 1 1/2" "								\$90 <input type="checkbox"/>
55		6" x 1 1/2" "								\$100 <input type="checkbox"/>
56		4" x 1 1/2" "								\$175 <input type="checkbox"/>
57	Aluminum wheel Soft-Roll® by Frog Legs ¹	5" x 1 1/2" "								\$185 <input type="checkbox"/>
58		6" x 1 1/2" "								\$195 <input type="checkbox"/>
59	Ultra light aluminum court wheel by Frog Legs	4" x 1" <input type="checkbox"/>			5" x 1" <input type="checkbox"/>					\$185 <input type="checkbox"/>
CASTER OPTIONS										
60	Caster Pin Locks (Available with +1" or +2" Stem bolt ONLY)									\$295 <input type="checkbox"/>
61	Frog Legs suspension forks (Not available with caster pin locks)				Specify weight of user ² :					\$475 <input type="checkbox"/>

1: Soft-Roll® are casters designed to absorb vibrations

2: Max weight: 220 lbs

FRONT Seat-to-floor Height - Please circle your choice

Forks		3" caster			4" caster			5" caster			6" caster			7" caster			8" caster		
		Stem bolt length			Stem bolt length			Stem bolt length			Stem bolt length			Stem bolt length			Stem bolt length		
		STD	+1"	+2"	STD	+1"	+2"	STD	+1"	+2"	STD	+1"	+2"	STD	+1"	+2"	STD	+1"	+2"
4"	P1	13	14	15	13 1/2"	14 1/2"	15 1/2"												
	P2	13 1/2"	14 1/2"	15 1/2"	14	15	16	14 1/2"	15 1/2"	16 1/2"									
5"	P1																		
	P2	13 3/4"	14 3/4"	15 3/4"	14 1/4"	15 1/4"	16 1/4"	14 3/4"	15 3/4"	16 3/4"									
	P3	14 1/4"	15 1/4"	16 1/4"	14 3/4"	15 3/4"	16 3/4"	15 1/4"	16 1/4"	17 1/4"	15 3/4"	16 3/4"	17 3/4"						
	P4	14 3/4"	15 3/4"	16 3/4"	15 1/4"	16 1/4"	17 1/4"	15 3/4"	16 3/4"	17 3/4"	16 1/4"	17 1/4"	18 1/4"	16 3/4"	17 3/4"	18 3/4"			
7"	P1										16 1/4"	17 1/4"	18 1/4"	16 3/4"	17 3/4"	18 3/4"			
	P2										16 3/4"	17 3/4"	18 3/4"	17 1/4"	18 1/4"	19 1/4"			
	P3										16 3/4"	17 3/4"	18 3/4"	17 1/4"	18 1/4"	19 1/4"	17 3/4"	18 3/4"	19 3/4"
	P4				16 3/4"	17 3/4"	18 3/4"	17 1/4"	18 1/4"	19 1/4"	17 3/4"	18 3/4"	19 3/4"	18 1/4"	19 1/4"	20 1/4"	18 3/4"	19 3/4"	20 3/4"
	P5	16 3/4"	17 3/4"	18 3/4"	17 1/4"	18 1/4"	19 1/4"	17 3/4"	18 3/4"	19 3/4"	18 1/4"	19 1/4"	20 1/4"	18 3/4"	19 3/4"	20 3/4"	19 1/4"	20 1/4"	21 1/4"

REAR WHEELS															
62	Mag			20" <input type="checkbox"/>		22" <input type="checkbox"/>		24" <input type="checkbox"/>		\$50 <input type="checkbox"/>					
63	Newton One - spoke wheel					22" <input type="checkbox"/>		24" <input type="checkbox"/>		NC <input type="checkbox"/>					
64	Newton Gravity - ultralight wheel			20" <input type="checkbox"/>		22" <input type="checkbox"/>		24" <input type="checkbox"/>	25" <input type="checkbox"/>	26" <input type="checkbox"/>	\$300 <input type="checkbox"/>				
65	Spinergy Spox					22" <input type="checkbox"/>		24" <input type="checkbox"/>	25" <input type="checkbox"/>	26" <input type="checkbox"/>	\$795 <input type="checkbox"/>				
66	Spinergy Spox spokes colors*				Blk <input type="checkbox"/>			Wht <input type="checkbox"/>	Blue <input type="checkbox"/>	Red <input type="checkbox"/>	Yel <input type="checkbox"/>	NC <input type="checkbox"/>			
67	Spinergy LX					22" <input type="checkbox"/>		24" <input type="checkbox"/>	25" <input type="checkbox"/>	26" <input type="checkbox"/>	\$995 <input type="checkbox"/>				
68	Spinergy LX spokes colors *				Blk <input type="checkbox"/>			Wht <input type="checkbox"/>	Blue <input type="checkbox"/>	Red <input type="checkbox"/>	Yel <input type="checkbox"/>	Pink <input type="checkbox"/>	Grn <input type="checkbox"/>	Org <input type="checkbox"/>	NC <input type="checkbox"/>

* Black spokes standard leadtime / Colour spokes up to 2 weeks extra leadtime

TIRES											
69	Soft Urethane 1 3/8 "			20" <input type="checkbox"/>		22" <input type="checkbox"/>		24" <input type="checkbox"/>		STD <input type="checkbox"/>	
70	Soft Urethane knoby 1 3/8 "							24" <input type="checkbox"/>		\$95 <input type="checkbox"/>	
71	Soft Urethane - Shox 1"							24" <input type="checkbox"/>	25" <input type="checkbox"/>	\$115 <input type="checkbox"/>	
72	Pneumatic 1 3/8 "			20" <input type="checkbox"/>		22" <input type="checkbox"/>		24" <input type="checkbox"/>	25" <input type="checkbox"/>	26" <input type="checkbox"/>	NC <input type="checkbox"/>
73	Pneumatic HP - Schwalbe Right Run 1"			20" <input type="checkbox"/>		22" <input type="checkbox"/>		24" <input type="checkbox"/>	25" <input type="checkbox"/>	26" <input type="checkbox"/>	\$95 <input type="checkbox"/>
74	Pneumatic HP - Schwalbe Marathon Plus 1"							24" <input type="checkbox"/>	25" <input type="checkbox"/>	26" <input type="checkbox"/>	\$125 <input type="checkbox"/>

REAR Seat-to-floor Height

WHEEL SIZES	Height available +/- 1/4"	Specify desired height: _____
20" wheel	12 1/2" to 17"	
22" wheel	13 1/2" to 18"	
24" wheel	14 1/2" to 19"	
25" wheel	15" to 19 1/2"	
26" wheel	15 1/2" to 20"	

HANDRIM		Assembly	Available				
75	Aluminum anodized	regular <input type="checkbox"/> narrow* <input type="checkbox"/>	20", 22", 24", 25", 26"			STD	<input type="checkbox"/>
76	Plastic coated		20", 22", 24", 25", 26"			\$95	<input type="checkbox"/>
77	Q-Grip High friction plastic coated*		22", 24", 26"			\$240	<input type="checkbox"/>
78	Natural Fit- Standard grip* (<i>oval aluminum handrim w/thumb piece</i>)		22", 24", 26"			\$375	<input type="checkbox"/>
79	Surge LT handrim* (<i>oval aluminum handrim w/ a rubber high friction strip</i>)	regular <input type="checkbox"/> narrow* <input type="checkbox"/>	22", 24"			\$375	<input type="checkbox"/>
80	Omit handrim					NC	<input type="checkbox"/>

* Not available with Mag wheels

AXLE							
81	Quick release axle					STD	<input type="checkbox"/>
82	Permanent axle					NC	<input type="checkbox"/>
83	Quad release axle					\$125	<input type="checkbox"/>
84	Amputee adjustable axle plate (<i>may interfere with swing away armrest</i>)					\$195	<input type="checkbox"/>

REAR WHEEL CAMBER ANGLE							
85	0° camber					STD	<input type="checkbox"/>
86	3° camber					NC	<input type="checkbox"/>
87	6° camber (<i>not available with amputee axle plate</i>)					\$150	<input type="checkbox"/>

WHEEL LOCK							
88	Push to lock					STD	<input type="checkbox"/>
89	Pull to lock					NC	<input type="checkbox"/>
90	6" removable extension handle					\$70	<input type="checkbox"/>
91	Newton Grade aid push to lock (<i>suggested with pneumatic or high pressure tire</i>)					\$150	<input type="checkbox"/>
92	Scissor Lock					\$240	<input type="checkbox"/>
93	Attendant Lock (<i>additional set of wheel lock mounted to the back</i>)					\$200	<input type="checkbox"/>


POSITIONING BELT (<i>2 points seatbelt attachment hardware is included with the helio C2</i>)							
94	Velcro adjustable					\$45	<input type="checkbox"/>
95	Auto buckle					\$85	<input type="checkbox"/>
96	Additional attachment hardware for 4 point seatbelt					\$15	<input type="checkbox"/>
97	BodyPoint Hip Belt - 2 points non padded	M <input type="checkbox"/>	L <input type="checkbox"/>			\$80	<input type="checkbox"/>
98	BodyPoint Hip Belt - 2 points padded	S <input type="checkbox"/>	M <input type="checkbox"/>	L <input type="checkbox"/>		\$150	<input type="checkbox"/>
99	BodyPoint Hip Belt - 4 points padded	M <input type="checkbox"/>	L <input type="checkbox"/>			\$190	<input type="checkbox"/>

ACCESSORIES							
100	Anti tippers					\$140	<input type="checkbox"/>
101	Unoccupied transit tie down					NC	<input type="checkbox"/>
102	Cane Carrier	Left <input type="checkbox"/>	Right <input type="checkbox"/>			\$195	<input type="checkbox"/>
103	Spoke guards*	Available: 20", 22", 24", 26"				\$115	<input type="checkbox"/>
104	Footrest neoprene impact guard					\$90	<input type="checkbox"/>
105	Seat pouch					\$75	<input type="checkbox"/>
106	Backpack					\$95	<input type="checkbox"/>
107	Velcro Closing Strap					\$50	<input type="checkbox"/>
108	Helio T-Shirt	S <input type="checkbox"/>	M <input type="checkbox"/>	L <input type="checkbox"/>	XL <input type="checkbox"/>	XXL <input type="checkbox"/>	\$25 <input type="checkbox"/>
109	Army Patches	Army <input type="checkbox"/>	Navy <input type="checkbox"/>	Air Force <input type="checkbox"/>	Marines <input type="checkbox"/>	\$90	<input type="checkbox"/>
110	Combat Wounded Purple Heart Patch					\$90	<input type="checkbox"/>

* Not available with Mag wheels

TOTAL							
--------------	--	--	--	--	--	--	--

CENTER OF GRAVITY



Choose position of rear wheel axle plate:


Measure distance X from front of back post to center of axle receiver.

X =	Most stable									Most responsive				
	3/4"	1"	1/4"	1/2"	3/4"	2"	1/4"	1/2"	3/4"	3"	1/4"	1/2"	3/4"	4"
	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

Possible conflict between a seat depth under 16" and a center of gravity over 2.5".

Possible conflict between a 26" rear wheel and a center of gravity over 2.5".

BACK UPHOLSTERY HEIGHT



3" of adjustability

The slip-on or adjustable tension back upholstery is attached to the push handle offering 3" of height adjustment with increment of 1/2".

Back upholstery height

Rigid backs are attached to the back cane and the height of the push handle is adjusted independently.

Fork and Stem bolt combinations

