

AL435T

TAILGATER LIFT



PENDANT CONTROLLER

- Easy-to-grip base with soft, non-slip rubber coating
- Backlit buttons for night use
- Built-in LED flash light to assist with visibility while loading
- On-Off button with shut-off timer to prevent unintended activation
- Magnetic backing for convenient stowage — elevated from water in trucks
- Attractive, modern aesthetics

400 LB LIFTING WITH POWER EXTENSION

Harmar's Tailgater is engineered with versatile power extension to safely load scooters and wheelchairs from the curbside or over an open tailgate into the truck bed.

FEATURES:

- Folding arm
- Unique offset arm with power extension
- Reaches 43" to curb or over 24" tailgate
- Passenger-side install only
- 3-year transferrable warranty

AL435T PRODUCT SPECIFICATIONS

LIFTING CAPACITY	400 lb	INSTALLED WEIGHT	81 lb
ADJUSTABLE LENGTH	31" - 42"	ADJUSTABLE HEIGHT	32" - 43"
SECUREMENT	Docking device <i>(not included)</i>	WARRANTY	<ul style="list-style-type: none">- 3-year- Transferable
APPLICATION	<ul style="list-style-type: none">- Manual wheelchairs- Power wheelchairs- Scooters	POWERED	<ul style="list-style-type: none">- Power lifting- Power rotation- Power extension
INCLUDED ACCESSORIES	<ul style="list-style-type: none">- Folding arm- Universal base adapter	AVAILABLE OPTIONS	<ul style="list-style-type: none">- Battery pack- Custom adapters- Wireless remote