

IQ7000 AUTOFOLD

CHAIR LINK



VIDEO



Battery



12AH or 20AH
lithium-ion

Motor



2 * 250 watts

Airline



Approved *

Cruise



Approved *

Net Weight



67 lbs

Range



Up to 17 miles*

Max Speed



4 mph

Max Slope



6 degrees

Width



26"

Capacity



350 lbs

FITS ON THE LEFT OR RIGHT ARMREST




360° Joystick



Climbing Gradient 6°

Customize Your Wheelchair



Airline & Cruise Approved

COLOR OPTIONS





FEATURES	IQ - 7000	IQ - 8000	IQ - 9000
Auto Fold	Yes (optional)	No	No
Auto Recline	No	No	Yes (optional)
Remote Controlled	Yes	Yes	Yes
Battery	12AH / 20AH* Lithium-Ion	12AH / 20AH* Lithium-Ion	2 x 7.5AH* Lithium-Ion
Motor	2*250W	2*250W	2*250W Brushless
Range	Up to 17 miles*	Up to 19 miles*	Up to 19 miles*
Max Speed	4 mph	4 mph	4 mph
Max Slope	Up to 6 degrees	Up to 6 degrees	Up to 10 degrees
Width	26"	23.5" or 26"	23.5" or 26"
Capacity	350 lbs	330 lbs	310 lbs
Travel	Airline/Cruise Approved**	Airline/Cruise Approved**	Airline/Cruise Approved**

*Driving habits, temperature, weather, loading weight, road conditions, and other factors may affect the overall range.

**Only some variants are approved.

Contact Us For Further Information

Phone: +1 (844) 755-5055

Email: info@comfygomobility.com

Address: 8746 Kimball Avenue, Chino, CA 91708