




youtube video

**WITH DUAL BATTERY****Battery**2 * 7.5AH
lithium-ion**Motor**2 * 250 watts
brushless**Airline**


Approved *

Cruise

Approved *

Net Weight

58 lbs

Range

Up to 19 miles*

Max Speed

4 mph

Max Slope

10 degrees

Width


23.5" or 26"

Capacity

310 lbs

FITS ON THE LEFT OR RIGHT ARMREST**360° Joystick**intelligent
electromagnetic
brake

Climbing Gradient 10°

**Customize Your Wheelchair**Seating Width**
17.5" or 20"Max Loading Weight
310 lbs2 * 7.5AH lithium-ion battery
Fits in most trunks & through most doorways**Airline & Cruise Approved**

COLOR OPTIONS

