

Invacare® Perfecto₂™ V Oxygen Concentrator

The new Invacare® Perfecto₂™ V Oxygen Concentrator is an economical addition to the premier Invacare stationary concentrator line. This 5-liter concentrator is small, quiet and lightweight and comes standard with the SensO₂® oxygen sensor. The Perfecto₂ V unit incorporates many of the best features of the market-leading Invacare® Perfecto₂™ concentrator, including its contemporary, patient-preferred external design and color.

Features

Top handle	For ease of delivery
Servicing parts compatibility	Compatible with many Platinum™ and Perfecto ₂ ™ concentrator parts
Alarms and safety systems	Battery-free power-loss alarm; Sieve-GARD™ monitors sieve performance; SensO ₂ ® oxygen monitor; Compressor 35 psi pressure-relief valve, low flow alarm
Pressure-based system	Superior to timing-based systems
Self-diagnostic electronics	Minimizes troubleshooting time
Simplified P.M. schedule	Three year oxygen purity check, replace P.M. components as needed



Specifications

Dimensions		Flow range	1/2 - 5 L/min.
Width	13"	Electrical requirements	120 VAC +/- 10%, 60Hz
Height	23"	Rated current input	3 A
Depth	11.5"	Pressure relief mechanism operational	35 psi +/- 3.5 psi
Weight	39 lb.	Filters	Cabinet, outlet HEPA, compressor inlet
Shipping weight	45 lb.	Altitude	8,000 ft.
Average power consumption	325 watts	Altitude to fill HomeFill® oxygen system	6,000 ft.
Number per pallet	27 units	Limited warranty	3 years
Oxygen output concentration levels	95.6% to 87% at all flow rates		
Maximum outlet pressure	5 +/- 0.5 psi		

Model Selection Guide

Model No.	Units/Case	Limited Warranty	Description
IRC5P02V	1	3 Years	Perfecto ₂ V 5-liter concentrator with SensO ₂
IRC5P	1	3 Years	Perfecto ₂ 5-liter concentrator
IRC5P02	1	3 Years	Perfecto ₂ 5-liter concentrator with SensO ₂
IOH200	1	3 Years	HomeFill® Oxygen Compressor
IOH270	1	1 Year	HomeFill® Ready-Rack (for use with all Perfecto ₂ models)



HomeFill® oxygen system connection.



Invacare Corporation
www.invacare.com

USA

One Invacare Way
Elyria, Ohio
44035-4190
(800) 333-6900

Canada

570 Matheson Blvd. E.
Unit 8
Mississauga, Ontario
L4Z 4G4 Canada
(800) 668-5324

© 2009

Invacare Corporation.
All rights reserved. All trademarks are identified by the symbols ™ and ®. All trademarks are owned by or licensed to Invacare Corporation unless otherwise noted. Specifications are subject to change without notification.

Form No. 09-064



INVACARE®

Perfecto₂TM V

Oxygen Concentrator



Yes, you can.®