

JAY[®]

GO[™] BACK

≡ A comfortable, wipeable and adjustable general use backrest



JAY® GO™ BACK

≡ **The JAY® GO™ Back was specially designed for comfort and adjustability.** The JAY Go Back features wide Velcro® straps that are close together for effective, comfortable and easy to adjust back support.

Comfortable

- Thick and soft covers
- Wide width straps with narrow gaps
- Padded back cane attachments

Wipeable

Outer cover

- Water resistant Dartex®
- Easy to clean
- Durable and comfortable surface

Adjustable

Wide Velcro Straps

- Require no tools for installation
- Simple to adjust back tension
- Provides accommodating postural support



JAY®

SPECIFICATIONS

JAY GO Back	
HCPCS Code	E2611 and E2612
Weight	1.7 lbs. (18WT)
Product Width	12" to 26"
Product Height	10" – 14" or 15" – 19"
Weight Capacity	300 lbs.up to 21", 500 lbs 22" and greater
Outer Cover	Dartex® and removable