

SPECIFICATIONS

| -Item # | Overall Width (open) | Overall Width (closed) | Overall Height | Seat To-Floor Height | Back Height | Width Seat | Depth Seat | Overall Length w/Riggings | Product Weight | Carton Shipping Weight | Weight Limit |
|------------|----------------------|------------------------|----------------|----------------------|-------------|------------|------------|---------------------------|----------------|------------------------|--------------|
| 12" | 20" | 12" | 32" | 15.5" | 16"- 18" | 12" | 12" | 36" | 30 lbs. | 35 lbs. | 250 lbs. |
| 14" | 22.5" | 13" | 32" | 15.5" | 16"- 18" | 14" | 14" | 36" | 31 lbs. | 36 lbs. | 250 lbs. |
| 16" | 22.5" | 12" | 36" | 19.5" | 18"- 20" | 16" | 16" | 42" | 31.5 lbs. | 37.2 lbs. | 300 lbs. |
| 18" | 24.5" | 12" | 36" | 19.5" | 18"- 20" | 18" | 16" | 42" | 32 lbs. | 37.4 lbs. | 300 lbs. |
| 20" | 28" | 12" | 36" | 19.5" | 18"- 20" | 20" | 16" | 42" | 34 lbs. | 39 lbs. | 300 lbs. |