

## LEVANTAR Floor Lift, Manual Base 500 lbs



- Weight capacity of 600 lbs
- Base option: Manual
- Elegant, contemporary design with soft lines but still heavy duty
- Oversized sealed rear casters roll easily and lock securely
- Easily lifts resident from floor
- Redundant controls allow for use even if hand control is inoperable
- Smart & Safety Technology – LCD display
  - o Visual battery fuel indicator on Hand Pendant and Lift Control Module
  - o Auditable low battery alert fuel indicator
  - o Lift safety stop if resident or caregiver accidentally gets in the way of the boom arm
  - o Audible Alert if the lift detects weight overload
  - o Maintenance monitoring: by months or lift cycles
  - o Audible Alert that maintenance is due
- NEW Six-point cradle with spring clips offer a more secure and wider foot print reducing sling squeeze, shallow height offers greater height range
- Can be charged on, or off the lift
- Includes 1 battery and external charger stand

### PRODUCTS

### Description

FLNP500	LEVANTAR Floor Lift, Manual Base 500 lbs, 1/ea
---------	--

### SPECIFICATIONS

### Description

Safety Working Load	Levantar 500 lbs, 227 kgs
Maximum Overall Length	52"
Maximum Overall Height	76.3"
Minimum Overall Height	27.7"
Maximum height at hook	6.9"
Minimum height at hook	2.1"
Height at maximum reach	37.8" from mast axis, 38.8" from deepest mast arc
Reach at maximum height	32.4" from mast axis, 33.4" from deepest mast arc
Turn diameter	55.9"
External width when open	45.8"
Internal width when open	40.3"
External width when closed	31.6"
Internal width when closed	26.1"
Height of legs	5.3"
Ground clearance	2.9"
Front caster	3"
Rear brake caster	5"

## SPECIFICATIONS

DESCRIPTION	DESCRIPTION
Net weight	132.2 lbs
Gross weight	156.5 lbs
Battery Pack weight	7.9 lbs (5.0A), 4.9 lbs (2.9A)
Lift Frame	Lifetime Limited
Electronics	5 yr
Battery/Hand Control	1 yr

## ACCESSORIES

DESCRIPTION	DESCRIPTION
13012	Sling, 1/bx
13019	Chains for Patient Lft, 1/bx
13023	atient Lift New Style, 1/cs
13025	Sling, 1/bx
13026	Sling w/Commode Opening, 1/bx
13046	Digital Scale Patient Lift, 1/bx
13057	Level Cradle 13244, 13240,13023, 1/bx
13060	Heavy-Duty Sling, 1/bx
13061	HD Sling w/ Commode Opening, 1/bx
13232	Straps for Patient Slings, 1 set/bx
13220L	U-Sling w/ Hd Sup.-Padded, 1/bx
13220M	U-Sling w/ Hd Sup.-Padded, 1/bx
13220S	U-Sling w/ Hd Sup.-Padded, 1/bx
13220XL	U-Sling w/ Hd Sup.-Padded, 1/bx
13221L	Full Body w/Com. Opening, 1/bx
13221M	Full Body w/Com. Opening, 1/bx
13221XL	Full Body w/Com. Opening, 1/bx
13221XXL	Full Body - Bariatric, 1/bx
13222L	Full Body, 1/bx
13222M	Full Body, 1/bx
13223L	Full Body, 1/bx
13223M	Full Body, 1/bx
13224XL	Full Body, 1/bx
13227D	Sling w/Hd. Sup. & Insert Pocket, 1/bx
13233D	Com. Sling w/Head Support, 1/bx
13235D	Sling w/o Head Support, 1/bx
13238D	Com. Sling w/o Head Support, 1/bx
13251D	Com. Sling w/Head Support, 1/bx
13254L	One Pce. Com. Sling w/Pos. Strap, 1/bx
13254M	One Pce. Com. Sling w/Pos. Strap, 1/bx
13254S	One Pce. Com. Sling w/Pos. Strap, 1/bx
13262L	Divided Leg Sling, 1/bx
13262M	Divided Leg Sling, 1/bx
13263A	Stand-Up, 1/bx
13263B	Stand-Up, 1/bx
13263E	Sit to Stand Padded Sling, 1/bx
13263F	Sit to Stand Padded Sling, 1/bx
13263H	Sit to Stand Buttock Strap, 1/bx
13264A	Transfer, 1/bx

## PARTS

DESCRIPTION	DESCRIPTION
SP07-2AAX0134	leg caster 3 "FLNP/FLP500 & FLNP/FLP600, 1/ea
SP07-2AAX0133	base caster 5"FLNP/FLP500 & FLNP/FLP600, 1/ea
SP07-TA37-2139-001	Boom Actuator FLNP/FLP500, 1/ea
SP07-TA12-2139-001	Boom Actuator FLNP/FLP600, 1/ea
SP07-TA16-2139-001	Leg spreader actuator FLP500/600, 1/ea
SP07-TC12-2139-001	Control Box FLNP/FLP500 & FLNP/FLP600, 1/ea