

NEW!

LC-570 S E R I E S





The new LC-570 Series 3-position chaise lounger features a unique spring design and a one-piece soft foam seat for outstanding comfort. Removable filling in back pillows for custom lumbar support adds to the LC-570 Series' comfort, making it a superior lift chair.

Features and Benefits:

- Removable filling in back pillows for custom lumbar support
- Available in 6 sizes Petite, Petite Wide, Small, Medium, Wide and Large
- Industry exclusive dual pockets for convenience
- Available in all fabrics including Ultrafabrics® Ultraleather™, Crypton® Super Fabrics and Outlast® Technology
- Heat and massage option
- · Footrest extension option available on all models



	LC-5/0P	LC-5/0PW	LC-5/0S	LC-5/0M	LC-5/0W	LC-5/0L
	3-Position Chaise Lounger 5'2" & Below	3-Position Chaise Lounger 5'3" & Below	3-Position Chaise Lounger 5'3" & Below	3-Position Chaise Lounger 5'4" - 6'0"	3-Position Chaise Lounger 5'4" - 6'0"	3-Position Chaise Lounger 5'9" - 6'2"
Weight Capacity	375 lbs.	375 lbs.	375 lbs.	375 lbs.	375 lbs.	375 lbs.
Seat-to-Floor	18"	19"	19"	20"	20"	21"
Seat Depth	18"	19"	19"	20"	20"	22"
Seat Width	18"	24"	19"	20"	24"	22"
Top of Back to Seat	22.5"	24"	24.5"	26.5"	26.5"	27"
Back Style	Blown Fiber	Blown Fiber	Blown Fiber	Blown Fiber	Blown Fiber	Blown Fiber
	Pillow Back	Pillow Back	Pillow Back	Pillow Back	Pillow Back	Pillow Back
Fully Padded Chaise	Yes	Yes	Yes	Yes	Yes	Yes
Heat & Massage Option	Yes	Yes	Yes	Yes	Yes	Yes
Distance From Wall	15.5"	16"	16.5"	19"	19"	24"
Overall Chair Width	35"	40"	35"	36"	40"	41.5"
Chair Weight	106 lbs.	114 lbs.	110 lbs.	113 lbs.	117 lbs.	114 lbs.
	7-Year Prorated	7-Year Prorated	7-Year Prorated	7-Year Prorated	7-Year Prorated	7-Year Prorated

Note: Measurements may vary +/- 1". Specifications subject to change without prior notice.

The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice.



