



LT-980

High strength Aluminum ultra-lightweight wheelchair

- Medicare code: K0004
- Ultra-light weight/Durable aluminum frame
- Light weight only: 24 lbs (w/o footrests)
- Folding backrest for travel or storage
- Fixed armrests with molded arm pads
- Durable molded side guards
- Swing-away/Detachable footrest with heel loop
- Height adjustable footplate
- Adjustable push to lock wheel brakes
- 6" x 1 " Flat free front polyurethane casters with durable forks
- 22" x 1" Spoke wheels with high profile/polyurethane tires
- Easy-Grip push-rims
- Carry pocket on backrest
- Breathable nylon upholstery
- Weight capacity: 250 lbs
- Frame color: Black

SPECIFICATIONS (subject to change without notice) :

Model No.	seat width (in.)	seat depth (in.)	armrest height (in.)	seat height (in.)	back height (in.)	overall height (in.)	overall open width (in.)	folded width (in.)	overall length (in.)	weight w/o riggings (lb.)	weight capacity (lb.)
LT-980	18	15	9	19	19	41	26	13	28	24	250

Model No.	Description
LT-980	18", aluminum, high strength lightweight, swing-away footrests
LT-980-E	18", aluminum, high strength lightweight, elevating footrests

Optional Accessories:

- Elevating Footrests - Anti-Tippers - Brake Extension - IV Pod Carrier - Oxygen Tank Holder - Seat Belt -
- Leg Straps - Carry Pouch - Headrest - Backrest Extension



KARMAN HEALTHCARE Inc.

19255 San Jose Avenue . City of Industry . CA 91748-1418

Order.1-800-805-2762 Tel.626-581-2235 Fax.626-581-2335