



Quickie LXI/LX January 2017

EILXI/EILX

Mark For:

Submitting for:

☐ Quote

☐ Order

Date: _____

PO#: _____

Dealer Acct #: _____

Dealer: _____

Dealer Contact: _____

Dealer Address: _____

Dealer City: _____

ST: _____

ZIP: _____

Dealer Phone: () _____

Fax: () _____

Confirmation Email: _____

Confirm Via: ☐ Fax ☐ Email

ADDITIONAL SHIPPING INFORMATION

Ship To: _____

Attention: _____

Address: _____

Address: _____

Ship To City: _____

ST: _____

ZIP: _____

Ship To Phone: () _____

Fax: () _____

The **HCPCS CODES** provided are based on code verification through PDAC or our interpretation of the code definitions.

The use of **HCPCS CODES** does not ensure coverage or payment. For coverage information, verify the appropriate payer's coverage policy.

✓	Option #	HCPCS	Definition	Price	✓	Option #	HCPCS	Definition	Price
---	----------	-------	------------	-------	---	----------	-------	------------	-------

1. FRAME STYLE

☐	065WM31	K0005	Quickie LXI (Adjustable Axle Plate) 265 lbs Weight Capacity	\$1,995.00	☐	065WM30	K0004	Quickie LX (Non Adjustable Axle Plate) 250 lbs Weight Capacity	\$1,895.00
--------------------------	---------	-------	--	------------	--------------------------	---------	-------	---	------------

2. FRAME SET UP

Transit Option

Includes stabilizer bar. Only on LXI. NA w/ padded swing away armrest

☐	065TO1	Transit System	\$315.00
--------------------------	--------	----------------	----------

Frame Width

☐	065W3	12"	NC
☐	065W4	13"	NC
☐	065W5	14"	NC
☐	065W6	15"	NC
☐	065W7	16"	NC
☐	065W8	17"	NC
☐	065W9	18"	NC
☐	065W41	19"	\$125.00
☐	065W11	E2201 20"	\$500.00

Seat Depth

12"-14" NA w/ 8" casters or 24" wheels

☐	065SD5	12"	NC
☐	065SD6	13"	NC
☐	065SD7	14"	NC
☐	065SD8	15"	NC
☐	065SD9	16"	NC
☐	065SD10	17"	NC
☐	065SD11	18"	NC
☐	065SD74	19"	\$125.00
☐	065SD13	E2203 20"	\$550.00

3. FRONT SEAT HEIGHT SET UP

Front Seat Height

Rear STFH ranges from 0" to 1 1/2" lower than front STFH

☐	065SH12	15"	NC
☐	065SH13	15.5"	NC
☐	065SH14	16"	NC
☐	065SH15	16.5"	NC
☐	065SH16	17"	NC
☐	065SH17	17.5"	NC
☐	065SH18	18"	NC
☐	065SH20	18.5"	NC
☐	065SH21	19"	NC
☐	065SH23	19.5"	NC
☐	065SH24	20"	NC
☐	065SH26	20.5" (Only available on swing away)	NC

FRONT SEAT HEIGHT ORDERING GUIDE

Wheel	Caster				
	4"	5"	6"	8" POLY	8"x2"PNU
20"	15.0	15.0	16.0		18.0
	15.5	16.0	16.5		
	16.0	16.5	17.0		
	16.5	17.0	17.5		
	17.0	17.5	18.0		
	17.5	18.0			
22"	18.0				
	16.0	16.0	16.0	18.5	19.0
	16.5	16.5	16.5	19.0	
	17.0	17.0	17.0		
	17.5	17.5	17.5		
	18.0	18.0	18.0		
24"	18.5	18.5	18.5		
	19.0	19.0	19.0		
	17.00	17.0	17.0	18.5	18.0
	17.50	17.5	17.5	19.0	19.0
	18.00	18.0	18.0	19.5	19.5
	18.50	18.5	18.5	20.0	
25"	19.00	19.0	19.0		
	19.50	19.5	19.5		
	20.0	20.0	20.0		
	17.00	17.0	17.0	18.5	19.0
	17.50	17.5	17.5	19.0	19.5
	18.00	18.0	18.0	19.5	
26"	18.50	18.5	18.5	20.0	
	19.00	19.0	19.0		
	19.50	19.5	19.5		
	20.0	20.0	20.0		
	18.00	18.0	18.0	18.5	18.0
	18.50	18.5	18.5	19.0	19.5

3. FRONT SEAT HEIGHT SET UP (CONT)

Caster Wheels			
☐	065CW3	4" Low Profile Polyurethane	NC
☐	065CW31	4" Micro Lighted	\$70.00
☐	065CW72	4" x 1" Aluminum Polyurethane	\$100.00
☐	065CW44	4" x 1.25" Semi Pneumatic	\$55.00
☐	065CW64	4" x 1.5" Aluminum Soft Roll	\$175.00
☐	065CW5	5" Low Profile Polyurethane	NC
☐	065CW51	5" Micro Lighted	\$85.00
☐	065CW65	5" x 1" Aluminum Polyurethane	\$110.00
☐	065CW7	5" x 1.5" Semi-Pneumatic	\$85.00
☐	065CW54	5" Soft Roll (Grey) NA w/ Frog Legs	\$85.00
☐	065CW66	5" x 1.5" Aluminum Soft Roll	\$185.00
☐	065CW34	6" Polyurethane	NC
☐	065CW33	E2214 E2215 6" Pneumatic	\$75.00
☐	065CW35	6" x 1.5" Semi Pneumatic	\$85.00
☐	065CW68	6" x 1.5" Aluminum Soft Roll	\$195.00
☐	065CW16	7" Polyurethane	NC
☐	065CW32	8" Polyurethane	NC
☐	065CW21	E2214 E2215 8" Pneumatic	\$75.00
☐	065CW39	8" x 1.5" Semi Pneumatic	\$85.00
☐	065CW36	E2214 E2215 8" x 2" Pneumatic	\$145.00
☐	065CW37	E2214 K0108 8" x 2" Pneumatic w/Airless Insert	\$225.00
☐	065CW38	Omit Caster Wheel	NC

Caster Forks			
☐	065CF62	3"	NC
☐	065CF9	4 3/4"	NC
☐	065CF12	5 1/4"	NC
☐	065CF17	6"	NC
☐	065CF23	7"	NC
☐	065CF29	8" x 2"	NC
☐	065CF30	Omit Caster Fork and Wheel	NC

Fork Stems			
<i>Caster pin locks N/A w/ 8"x2" casters</i>			
☐	065CS1	Standard Stem	NC
☐	065CS5	3/4" Stem	NC
☐	065CS9	1 1/2" Stem	NC
☐	065CS7	K0073 Standard Stem w/Pin Locks	\$195.00
☐	065CS11	K0073 3/4" Stem w/Pin Locks	\$195.00
☐	065CS4	K0073 1 1/2" Stem w/Pin Locks	\$195.00

4. HANGER & FOOTPLATE OPTIONS

Hangers			
☐	065H145	60° Swing In/Out Hanger	NC
☐	065H139	70° Swing In/Out Hanger	NC
☐	065H158	90° Swing In/Out Front Mount <i>Only available Ang adj adult & kids footplates</i>	NC
☐	065H159	K0053 Articulating Hanger NA w/ Platform Flip Up and Alum Angle Adj Front Mount footplates.	\$350.00
☐	065H160	E0990 Elevating Hanger - Ext Mount NA w/ Platform Flip Up and Alum Angle Adj Front Mount footplates.	\$295.00
☐	065H161	E0990 Elevating Hanger - Front Mount Only available w/ front mount alum angle adj footplates and 14"-22" width	\$295.00
☐	065H32	Omit Hanger	-\$15.00

Footplate & Hanger Options (Ext Tube & Direct Mount)			
☐	065HL1	E0951 Heel Loop	\$30.00
☐	065HL5	K0038 Leg Strap	\$35.00
☐	065HL4	Omit Heel Loop	NC
☐	065HL6	Omit Leg Strap	NC
☐	065IG	Impact Guard Neoprene	\$85.00

Footplates			
<i>Heel loop or leg strap not included, See Footplate Options Section</i>			
☐	065F1	Composite (14"-20"w) Only available w/ 60°, 70°, and ALR	NC
☐	065F48	K0040 Platform Flip Up Available w/ 60°, 70°, or hemi hangers	\$190.00
☐	065F5	Aluminum Footplate	NC
☐	065F3	K0040 Aluminum Angle Adj Adult	\$185.00
☐	065F4	K0040 Aluminum Angle Adj Kids NA w/ ELR. Only on 12"-16" wide	\$185.00
☐	065F7	Omit Footplate & Ext Tube	-\$15.00

Extension Tube			
☐	065E1	Short (See Matrix for Lengths) NA w/ 90° hangers and ELR	NC
☐	065E2	Medium (See Matrix for Length)	NC
☐	065E5	Long (See Matrix for Length) NA w/ 90° hangers and ELR	NC
☐	065E14	Regular ELR	NC
☐	065E15	Long ELR	NC
☐	065E16	K0037 High Mount Only available with angle adj footplates or 60° and 70° hangers	\$55.00

Seat Sling to Top of Footrest

Hanger	Extension	Seat Sling to Top of Footrest in inches		*Minimum Allowable FSH
		Min	Max	
60° Swing In/Out	Short	11.5	15	15
	Medium	11.5	17	15
	Long	14	18.5	17
70° Swing In/Out	Short	13	18	15
	Medium	13	20	15.5
	Long	16	20.5	17.5
60° Swing In/Out w/Platform Flip Up	Short	11.5	14.5	14
	Long	14	18.5	16
70° Swing In/Out w/Platform Flip Up	Short	13	15.5	15
	Long	15.5	20	17.5
90° Front Mount	No Ext	5.5	11.5	14
	Short	16.5	20	15
	Medium	17.5	20.5	16
ALR	Long	19	20.5	17
	Short	15.5	19.5	15
	Medium	15.5	20.5	16
ELR Ext. Mount	Long	17	20.5	17
	No Ext	5.5	9.5	13

5. REAR SEAT HEIGHT SET UP

Axle			
☐	065AX19	Threaded	NC
☐	065AX3	Quick Release Steel	NC
☐	065AX13	Quad Quick Release	\$185.00

One Arm Drive			
<i>Only available w/ Mag Wheels. Select Standard Rim and either right or left mount</i>			
☐	065AL16	E0958 One Arm Drive Adj. Axle Plate	\$850.00
☐	065OA2	Standard Handrim	NC
☐	065HM1	Right Hand Mount	NC
☐	065HM2	Left Hand Mount	NC

Rear Wheel Type			
<i>Contact Built 4 Me for Spinergy Color Options</i>			
☐	065RW3	20" Mag Wheel	NC
☐	065RW8	22" Mag Wheel	NC
☐	065RW1	24" Mag Wheel	NC
☐	065RW4	20" Spoke Wheel	\$200.00
☐	065RW9	22" Spoke Wheel	\$200.00
☐	065RW2	24" Spoke Wheel	NC
☐	065RW11	26" Spoke Wheel	\$65.00
☐	065RW95	24" Spinergy SPOX Black	\$795.00
☐	065RW96	25" Spinergy SPOX Black	\$795.00
☐	065RW97	26" Spinergy SPOX Black	\$795.00
☐	065RW98	24" Spinergy LX Black	\$995.00
☐	065RW99	25" Spinergy LX Black	\$995.00
☐	065RW50	24" Mountain	\$250.00
☐	065RW25	Omit Rear Wheel and Axle (Intended Size for Set Up _____)	-\$55.00

Rear Wheel Tire			
☐	065RT1	E2211 E2212	Pneumatic NC
☐	065RT2	E2211 E2213	Pneumatic Airless Insert \$80.00
☐	065RT3	E2218	Full Poly NC
☐	065RT6	E2211 E2212	High-Pressure Clincher \$95.00
☐	065RT21	E2211 E2212	Schwalbe Marathon Plus Evolution \$125.00
☐	065RT38	E2218	SHOX Performance Poly \$95.00
☐	065RT12	E2211 E2212	Mountain Tire <i>Only available on 24" Mountain Rear Wheel</i> \$95.00

Tire	Wheel	Mag	Lite Spoke	Ultra Lightweight	Spinergy
Pneumatic		20,22,24	20,22,24,25, 26		24,25,26
Pneumatic w/Airless Insert		20,22,24	20,22,24,25,26		24,25,26
Full Poly Tire		20,22,24	20,22,24		24
High Pressure Clincher			22,24,25,26	22,24	24,25,26
Schwalbe Marathon			24,25,26	22,24	24,25,26
SHOX Performance Poly		24	24,25		24,25
22" = 489 mm 24" = 540 mm 25" = 559 mm 26" = 590 mm					

Handrim			
<i>Natural Fits available in 22", 24", and 25".</i>			
☐	065HR1	Aluminum Anodized	NC
☐	065HR4	Plastic Coated -NA on all rims	\$95.00
☐	065HR22	Natural Fit - Standard Grip	\$375.00
☐	065HR23	Natural Fit - Super Grip	\$410.00
☐	065HR7	Omit Hand Rims	NC
☐	065HR8	8 Oblique <i>Only available w/ aluminum handrims. NA w/ Kik</i>	\$325.00

6. BACKREST

Back Type			
<i>Angle Adj std w/ push handles and NA w/ all flip back armrests</i>			
☐	065BU1	Standard	NC
☐	065BT20	Angle Adj 3-5" Pivot (regular)	\$185.00
☐	065BT22	Angle Adj 5-7" Pivot (tall)	\$185.00

Angle Adjustable Backrest Options							
Angle Adjustable Pivot Height							
☐	065AH1	3"	NC	☐	065AH4	6"	NC
☐	065AH2	4"	NC	☐	065AH5	7"	NC
☐	065AH3	5"	NC				
Angle Adjustable Pre-Set Angle							
☐	065BD5	Back Angle - 15°					NC
☐	065BD4	Back Angle - 12°					NC
☐	065BD3	Back Angle - 6°					NC
☐	065BD2	Back Angle - 0°					NC
☐	065BD1	Back Angle +3°					NC
Angle Adjustable Back Tube Height (See Chart)							
☐	065B1	Low					NC
☐	065B5	Medium					NC
☐	065B9	Tall					NC

Angle Adjustable Back Height (Optional)							
☐	065UX10	12"	NC	☐	065UX20	17"	NC
☐	065UX12	13"	NC	☐	065UX22	18"	NC
☐	065UX14	14"	NC	☐	065UX24	19"	NC
☐	065UX16	15"	NC	☐	065UX26	20"	NC
☐	065UX18	16"	NC	☐	065UX28	21"	NC

JAY Back (Optional)			
☐	065JB01	J3 Back (See J3 Order Form)	Price List

Angle Adjustable Backrest Set Up

REFER TO CHART IN DETERMINING PUSH HANDLE LENGTH

PPS

SEATING SURFACE CUSHION THICKNESS

TO DETERMINE THE PROPER PIVOT-POINT HEIGHT, MEASURE BETWEEN THE PPS AND THE TOP OF THE SEATING SURFACE (A), THEN ADD THE THICKNESS OF THE CUSHION (B) IF ONE IS USED. IF A DROP-SEAT FALLS BELOW THE SEAT RAILS, SUBTRACT THIS DISTANCE FROM THE OVERALL VALUE.

REGULAR 3"-5"

TALL 5"-7"

Angle Adjustable Backrest Chart

Back Tube Height	Low	Medium	Tall
Pivot regular set at 3"	12" 13" 14"	15" 16" 17"	17" 18" 19"
Pivot regular set at 4"	13" 14" 15"	16" 17" 18"	18" 19" 20"
Pivot regular set at 5"	14" 15" 16"	17" 18" 19"	19" 20" 21"
Pivot tall set at 5"	14" 15" 16"	17" 18" 19"	19" 20" 21"
Pivot tall set at 6"	15" 16" 17"	18" 19" 20"	N/A
Pivot tall set at 7"	16" 17" 18"	19" 20" 21"	N/A

Adjustable pivot point present measurements from seat sling. If back height option is not chosen, middle option will be present.

6. BACKREST (CONT)**Backrest Upholstery**

☐	065BU1	Standard Upholstery	NC
☐	065BU60	Tension Adjustable Upholstery	\$175.00
☐	065BU14	Omit Back Upholstery	-\$15.00

Backrest Options

☐	065BU16	Reinforced Back <i>Only available w/ Std Upholstery</i>	\$30.00
☐	065SZ1	Stabilizer Bar	\$105.00

Push Handles

☐	065BP1	Push Handles Straight	NC
☐	065BP3	Push Handles 8" (Std)	\$25.00

Back Height

☐	065B92	16"-19"	NC
☐	065B86	14"-17"	\$50.00

7. ARMREST**Armrest**

☐	065AA16	E0973	Padded Swing-Away Adult	\$210.00
☐	065AA17	E0973	Padded Swing-Away Kids	\$210.00
☐	065AA5		Flip-Back	NC
☐	065AA8	E0973	Single Post Ht. Adjustable Adult	\$285.00
☐	065AA9	E0973	Single Post Ht. Adjustable Kids	\$285.00
☐	065AA6	E0973	Flip-Back Height Adjustable	\$295.00
☐	065AA14		Omit Armrest	NC

Armpad

NA w/ padded swing away armrests

☐	065AM1	Desk Length - (10")	NC
☐	065AM2	Full Length - (14")	NC

8. KEY OPTIONS & ACCESSORIES**Anti-Tip**

☐	065AT4	E0971	Rear Anti-Tip	\$140.00
--------------------------	--------	-------	---------------	----------

Wheel Lock

Additional no gos may exist depending on wheel size, seat depth and axle position

☐	065WL1		High Mount Push	NC
☐	065WL2		High Mount Pull	NC
☐	065WL30		Compact	\$125.00
☐	065WL10		Omit	NC
☐	065WL12		Do Not Mount	NC

Grade Aid

☐	065GA1	E0974	Grade Aid NA w/ poly, high-pressure clincher tires	\$150.00
--------------------------	--------	-------	--	----------

Extension Handle

☐	065EH1	E0961	6" Pair	\$75.00
--------------------------	--------	-------	---------	---------

JAY Cushions (Optional)

To order any JAY cushions on-chair, please see the JAY price list.

Cushion name or part number must be filled out with desired width and depth

☐	065JC20	JAY Bundled Cushion		
		Cushion Name or Part #		Price List
		Width: _____	Depth: _____	

Cushion Options

☐	065CU1		2" Cushion	\$120.00
☐	065CU4		Omit Seat Sling	-\$15.00

Positioning Belts

☐	065A142	E0978	Belt 1.5" Auto Buckle	\$70.00
☐	065A12	E0978	Belt 2" Aircraft Buckle	\$110.00
☐	065A13	E0978	Belt 2" Aircraft Buckle Pad	\$120.00

Accessories

☐	065A1		Quickie Backpack Black	\$55.00
☐	065A2		Quickie Kids Backpack Black	\$50.00
☐	065A8		Quickie Seat Pouch Black	\$45.00
☐	065A29	E0950	Tray Table	\$250.00
☐	065A23	K0065	Black Spoke Guards	\$115.00
☐	065A24	K0065	Quickie Clear Spoke Guards	\$115.00
☐	065A54	K0038	Extra Leg Strap Black	\$30.00
☐	065SG4	K0108	Side Guard Plastic Kids	\$125.00
☐	065SG3	K0108	Side Guard Plastic Regular	\$125.00
☐	065A82		Caddy	\$45.00
☐	065A49	E2207	Crutch Holder	\$195.00
☐	065A83		Tool Kit	\$45.00
☐	065A500		Velcro Gel Cover for Hanger/Frame	\$105.00

Military Service Patches (Optional)

Only available on tension adjustable upholstery and J3 Back

☐	065VU1		Marines Patch	\$90.00
☐	065VU2		Army Patch	\$90.00
☐	065VU3		Navy Patch	\$90.00
☐	065VU4		Air Force Patch	\$90.00
☐	065VU12		Coast Guard Patch	\$90.00
☐	065VU5		Purple Heart Patch	\$90.00
☐	065VU8		Purple Heart with Marines Patch	\$145.00
☐	065VU9		Purple Heart with Army Patch	\$145.00
☐	065VU10		Purple Heart with Navy Patch	\$145.00
☐	065VU11		Purple Heart with Air Force Patch	\$145.00
☐	065VU13		Purple Heart with Coast Guard Patch	\$145.00

8. KEY OPTIONS & ACCESSORIES (CONT)

Standard Frame Colors			
☐	065C140	Aztec Gold	NC
☐	065C2	Black <i>(Not avail as base w/Zebra)</i>	NC
☐	065C141	Black Cherry	NC
☐	065C10	Black Opal	NC
☐	065C38	Blue Opal	NC
☐	065C36	Candy Blue	NC
☐	065C14	Candy Purple	NC
☐	065C35	Candy Red	NC
☐	065C142	Electric Blue	NC
☐	065C143	Evergreen	NC
☐	065C126	Glow	NC
☐	065C144	Green Apple	NC
☐	065C147	Hot Sparkle Pink	NC
☐	065C148	Mauve Pink	NC
☐	065C103	Rootbeer	NC
☐	065C114	Sparkle Silver	NC
☐	065C145	Sunrise Orange	NC
☐	065C113	Titanium Color Paint	NC
☐	065C8	Yellow	NC

Matte Frame Colors			
☐	065C116	Matte Black	NC
☐	065C151	Matte Black Cherry	NC
☐	065C150	Matte Electric Blue	NC
☐	065C149	Matte Evergreen	NC
☐	065C152	Matte Purple	NC

Custom Frame Color			
☐	065C112	Desert Camouflage	\$295.00
☐	065C43	Mossy Oak Camouflage	\$295.00
☐	065C42	Stars & Stripes	\$295.00
☐	065C108	Zebra <i>(Choose above base color) Std with White unless above base paint is chosen NA w/ Black</i>	\$295.00

Built for Me

To inquire about any unique customization for Quickie products, contact: Phone: (800) 333 - 4000, Fax: (800) 333 9011

