



Invacare® Top End® Excelerator™ Series Handcycles

If you want a great way to exercise, cross-train or just have fun, one of the **Invacare® Top End® Excelerator™ Series Handcycles** is just what you need. This upright and extremely maneuverable three-wheel handcycle is capable of speeds up to 15 mph. The smooth Shimano® hub with its reversing drum brake delivers hands-on braking and allows the rider to tackle hills or level terrain with equal ease. The oval tube carbon steel frame construction with adjustable sliding seat makes transfers in and out easier too. There are two great models to choose from: the stocked model can be shipped right away and the custom model allows you to choose your favorite color and seat width.


The Excelerator Handcycle offers some great options too. You can order the popular mountain drive option which adds 7 extra low gears for climbing hills or the V crankset option with its wider cranks to allow more clearance for those who are tall or require additional space. No excuses...Get Cranking! Check out our childrens' models too!

Invacare Top End Excelerator Features

Carbon steel frame and fork
Adjustable seat: forwards and backwards
Adjustable back height
Fold over back upholstery
Adjustable crank height
Chain tension idler
Reflectors and safety flag
Adjustable height/angle footrest system
Shimano® internal hub with reversing drum brake

Parking brake
Single chain ring
Ovalized ergonomic hand pedals
Performance wheels with cruiser tires
Quick release axles
3° camber
Cushion and hook and loop seat restraint
Weight capacity 350 lb.

Model #: XCL/XCLSTK



1. Twist Shifter for 7 speed hub
2. Quad Twist shifter option
3. Mountain drive transmission option
4. Frictionless ovalized handpedals



Options

Crank width adapter set
Mountain drive transmission adds 7 low gears
Quad twist shifter helps those with hand impairments
Tri-pin quad handpedals
QuadGrip handpedals with 3" width extensions
Water bottle and cage
Tow bar for wheelchair
Cordless computer

Front safety light
Rear flashing light
Helmet
Seat restraint: Auto style fastener
Mirror
Crutch holder and strap
Bike rack for trailer hitch
Top End t-shirt

The Invacare® Top End® Li'l Excelerator™ and Invacare® Top End® Li'l Excelerator™-2 are sized just right for children and are available in 7 speeds and 3 speeds, respectively.

Invacare Top End Li'l Excelerator-2 Features

Aluminum frame and transport ready fork
Adjustable seat: forwards and backwards
Adjustable back height
Adjustable tension back upholstery
Adjustable crank height
Chain tension idler
Chain guard, reflectors and safety flag
Adjustable height/angle footrest system
Shimano® 3 speed hub with reversing drum brake

Locking parking brake
153 mm crankset with one chain ring
Ergonomic handpedals
20" Cruiser wheels and tires
Quick release axles
9° camber
Cushion and hook and loop seat restraint
Candy Apple Red paint
Weight capacity 175 lb.



Li'l Excelerator-2 - Handcycle



Model #: LXC2H
Candy Apple Red

1. Parking brake
2. Ergonomic handpedals
3. Hook and loop safety seat restraint

Invacare Top End Li'l Excelerator Features

Carbon steel frame and fork
Adjustable seat: forwards and backwards
Adjustable back height
Fold over back upholstery
Adjustable crank height
Chain tension idler
Full chain guard, reflectors and safety flag
Adjustable height/angle footrest system
Shimano® 7 speed hub with reversing drum brake
Parking brake
155 mm crankset with one chain ring
Ergonomic handpedals

20" Cruiser wheels and tires
Quick release axles
9° camber
Cushion and hook and loop seat restraint
Choice of 11 colors
Weight capacity 250 lb.



Model #: LXC
Royal Blue

Options for LXC2H and LXC

Foam covered horizontal handpedals
Quad twist shifter option
Mirror

Water bottle and cage
Helmet
Front safety light
Rear flashing light


Crutch holder and strap
Bike rack for trailer hitch
Seat restraint: Auto style fastener

Optional upgrades:



1. V Crankset - Only available on custom XCL

The Invacare® Top End® V crankset option accommodates consumers who might have “clearance” issues. Only available on the Invacare Top End Custom Excelerator Handcycle.



2. Mountain Drive Transmission (Only available on custom)

Available on the Excelerator handcycle
Adds 7 low gears (7 low + 7 standard = 14 total)
Easy shifting, especially for quadriplegics – just push in the center button
Lightweight planetary transmission
A must-have for hilly terrain



3. QuadGrip Handpedals

QuadGrip handpedals with 3” width extensions and QuadWraps are a great option for those with limited hand function. Add the Quad twist shifter option and the mountain drive for optimal function and fun.

Invacare Corporation
www.invacare.com

USA

One Invacare Way
Elyria, Ohio
44035-4190
(800) 333-6900

Canada

570 Matheson Blvd. E.,
Unit 8
Mississauga, Ontario
L4Z 4G4 Canada
(800) 668-5324

Top End

www.topendwheelchair.com
Sports and Recreation
Products
4501 63rd Circle North
Pinellas Park, Florida
33781-5914
(800) 532-8677
Fax: (727) 522-1007

©2018
Invacare Corporation.
All rights reserved.
Trademarks are identified
by the symbols ™ and ®.
All trademarks are owned
by or licensed to Invacare
Corporation unless
otherwise noted.
Shimano is a registered
trademark of Shimano
Inc. Specifications are
subject to change without
notification.
Form No. 18-175
180733

TOP END®