



LIBERTY FT Basic

Phone (800) 981-1540 | Fax (715) 254-0996
 sales@kimobility.com | www.kimobility.com

Company: _____

Ship To: _____

Account Number: _____ Date: _____

Address: _____

Quote PO Number: _____

City: _____ State: _____

RTS: _____

Zip: _____ Phone: _____

Marked For: _____

Email/Fax: _____

Select if chair quote has: Additional Ki Products with Different Quote/Order Number **Quote/Order Number:** _____

Liberty FT Basic E1161 **\$2795**
 250 lb weight capacity. Tilt range: 0° - 20°. Standard seat upholstery, Push to Lock wheel locks.

Effective 08/19/2020

HCPCS codes provided are based on code verification by the PDAC or our interpretation of the code definitions. Proper use of HCPCS codes does not ensure coverage or payment. For coverage information, verify the appropriate payer's coverage policy.

Frame

Seat Width
 14"-16" Seat Widths ship with 300 N Assistive Gas Strut, 18"-22" Seat Widths ship with 450 N Assistive Gas Strut

14" _____ N/C
 16" _____ N/C
 18" _____ N/C
 20" E2201 _____ \$470
 22" E2201 _____ \$470

Seat Depth

14" _____ N/C
 16" _____ N/C
 18" _____ N/C
 20" E2203 _____ \$605

Transit

Transit System _____ \$320
 Omit Transit _____ N/C

Backrest

| Back Post / Height Range / Back Height | | | | | | | |
|---|-------|---------------------------------|-----------------------------|-----------------------------|-----------------------------|-----------------------------|-----------------------------|
| Back Post | Price | Height Range | Back Height | | | | |
| <input type="checkbox"/> Height Adjustable 8° Bend with Push Handle | \$25 | <input type="checkbox"/> Short | <input type="checkbox"/> 14 | <input type="checkbox"/> 15 | <input type="checkbox"/> 16 | <input type="checkbox"/> 17 | <input type="checkbox"/> 18 |
| | | <input type="checkbox"/> Medium | <input type="checkbox"/> 16 | <input type="checkbox"/> 17 | <input type="checkbox"/> 18 | <input type="checkbox"/> 19 | <input type="checkbox"/> 20 |
| | | <input type="checkbox"/> Tall | <input type="checkbox"/> 18 | <input type="checkbox"/> 19 | <input type="checkbox"/> 20 | <input type="checkbox"/> 21 | <input type="checkbox"/> 22 |

Backrest Options

Fixed (Non-Folding) _____ Std
 Angle Adjustable (Folding) _____ \$95

Back Upholstery

If selecting both upholstery and Axiom Back, the upholstery will be installed

Standard _____ Std
 Omit Back Upholstery _____ - \$15

Footrest

| Hanger / Front Frame Bend | | | | |
|---|--------------------------|--------------------------|--------------------------|--------|
| Hanger | 60° | 70° | | Retail |
| Extension Tube Mount | <input type="checkbox"/> | <input type="checkbox"/> | | Std |
| Pro Elevating Leg Rests (PRO ELR) E0990, Rental K0195 | | | <input type="checkbox"/> | \$295 |
| Omit Hangers | | | <input type="checkbox"/> | N/C |

Footplate

Composite _____ Std
 Aluminum Angle Adjustable K0040 _____ \$190

Footrest continued

Footrest / Extension / Seat to Footrest Length

Select Extension and Seat to Footrest Length in row corresponding to desired Footrest. Footrest Length of at least 2.5" shorter than Seat Height recommended. Extension Tube Hanger measurements based on 70° Hanger Bend.

| Footrest | Extension | Seat to Footrest Length | | | | | | | | |
|--------------------------------------|---------------------------------|-------------------------------|-------------------------------|-------------------------------|-------------------------------|-------------------------------|-------------------------------|-------------------------------|-------------------------------|-------------------------------|
| Extension Tube Hangers (60°, 70°) | <input type="checkbox"/> Short | <input type="checkbox"/> 10 | <input type="checkbox"/> 10.5 | <input type="checkbox"/> 11 | <input type="checkbox"/> 11.5 | <input type="checkbox"/> 12 | <input type="checkbox"/> 12.5 | <input type="checkbox"/> 13 | <input type="checkbox"/> 13.5 | |
| | <input type="checkbox"/> Medium | <input type="checkbox"/> 11.5 | <input type="checkbox"/> 12 | <input type="checkbox"/> 12.5 | <input type="checkbox"/> 13 | <input type="checkbox"/> 13.5 | <input type="checkbox"/> 14 | <input type="checkbox"/> 14.5 | <input type="checkbox"/> 15 | <input type="checkbox"/> 15.5 |
| | <input type="checkbox"/> Long | <input type="checkbox"/> 14 | <input type="checkbox"/> 14.5 | <input type="checkbox"/> 15 | <input type="checkbox"/> 15.5 | <input type="checkbox"/> 16 | <input type="checkbox"/> 16.5 | <input type="checkbox"/> 17 | <input type="checkbox"/> 17.5 | <input type="checkbox"/> 18 |
| PRO ELR Low Pivot Setting | <input type="checkbox"/> Short | 10.5 - 15.5 | | | | | | | | |
| | <input type="checkbox"/> Long | 13 - 19 | | | | | | | | |
| PRO ELR High Pivot Setting | <input type="checkbox"/> Short | 8.5 - 17 | | | | | | | | |
| | <input type="checkbox"/> Long | 11.5 - 17.5 | | | | | | | | |

Heel Loop

- Heel Loops E0951 _____ \$30
- Omit Heel Loops _____ N/C

Calf Strap

- Velcro Adjustable K0038 _____ \$35

Hanger Release

- Option Required*
- 4-Way - Out _____ N/C

Seat Height, Casters and Rear Wheels

Seat Height / Caster Size / Rear Wheel Size

Select Seat Height corresponding to desired Rear Wheel and Caster Size

| Rear Wheel Size | Caster Size | | |
|-----------------|--|-------------------------------|-------------------------------|
| | 4" | 5" | 6" |
| 16" | Seat Height | | |
| | <input type="checkbox"/> 13.5 ¹ | <input type="checkbox"/> 14 | <input type="checkbox"/> 15 |
| | <input type="checkbox"/> 14 ¹ | <input type="checkbox"/> 14.5 | <input type="checkbox"/> 15.5 |
| | <input type="checkbox"/> 14.5 | <input type="checkbox"/> 15 | <input type="checkbox"/> 16 |
| | <input type="checkbox"/> 15 | <input type="checkbox"/> 15.5 | <input type="checkbox"/> 16.5 |
| | <input type="checkbox"/> 15.5 | <input type="checkbox"/> 16 | <input type="checkbox"/> 17 |
| | <input type="checkbox"/> 16 | <input type="checkbox"/> 16.5 | <input type="checkbox"/> 17.5 |
| | | <input type="checkbox"/> 17 | <input type="checkbox"/> 18 |
| | | <input type="checkbox"/> 17.5 | <input type="checkbox"/> 18.5 |
| | | <input type="checkbox"/> 18 | <input type="checkbox"/> 19 |
| | | <input type="checkbox"/> 18.5 | |
| | | <input type="checkbox"/> 19 | |
| 20" | <input type="checkbox"/> 13 | <input type="checkbox"/> 14 | <input type="checkbox"/> 15 |
| | <input type="checkbox"/> 13.5 | <input type="checkbox"/> 14.5 | <input type="checkbox"/> 15.5 |
| | <input type="checkbox"/> 14 | <input type="checkbox"/> 15 | <input type="checkbox"/> 16 |
| | <input type="checkbox"/> 14.5 | <input type="checkbox"/> 15.5 | <input type="checkbox"/> 16.5 |
| | <input type="checkbox"/> 15 | <input type="checkbox"/> 16 | <input type="checkbox"/> 17 |
| | <input type="checkbox"/> 15.5 | <input type="checkbox"/> 16.5 | <input type="checkbox"/> 17.5 |
| | <input type="checkbox"/> 16 | <input type="checkbox"/> 17 | <input type="checkbox"/> 18 |
| | | <input type="checkbox"/> 17.5 | <input type="checkbox"/> 18.5 |
| | | <input type="checkbox"/> 18 | <input type="checkbox"/> 19 |
| | | <input type="checkbox"/> 18.5 | |
| | | <input type="checkbox"/> 19 | |

| Rear Wheel Size | Caster Size | | |
|-----------------|-------------------------------|-------------------------------|-------------------------------|
| | 4" | 5" | 6" |
| 22" | Seat Height | | |
| | <input type="checkbox"/> 14 | <input type="checkbox"/> 14 | <input type="checkbox"/> 15 |
| | <input type="checkbox"/> 14.5 | <input type="checkbox"/> 14.5 | <input type="checkbox"/> 15.5 |
| | <input type="checkbox"/> 15 | <input type="checkbox"/> 15 | <input type="checkbox"/> 16 |
| | <input type="checkbox"/> 15.5 | <input type="checkbox"/> 15.5 | <input type="checkbox"/> 16.5 |
| | <input type="checkbox"/> 16 | <input type="checkbox"/> 16 | <input type="checkbox"/> 17 |
| | | <input type="checkbox"/> 16.5 | <input type="checkbox"/> 17.5 |
| | | <input type="checkbox"/> 17 | <input type="checkbox"/> 18 |
| | | <input type="checkbox"/> 17.5 | <input type="checkbox"/> 18.5 |
| | | <input type="checkbox"/> 18 | <input type="checkbox"/> 19 |
| | | <input type="checkbox"/> 18.5 | |
| | | <input type="checkbox"/> 19 | |
| 24" | <input type="checkbox"/> 15 | <input type="checkbox"/> 14.5 | <input type="checkbox"/> 15 |
| | <input type="checkbox"/> 15.5 | <input type="checkbox"/> 15 | <input type="checkbox"/> 15.5 |
| | <input type="checkbox"/> 16 | <input type="checkbox"/> 15.5 | <input type="checkbox"/> 16 |
| | | <input type="checkbox"/> 16 | <input type="checkbox"/> 16.5 |
| | | <input type="checkbox"/> 16.5 | <input type="checkbox"/> 17 |
| | | <input type="checkbox"/> 17 | <input type="checkbox"/> 17.5 |
| | | <input type="checkbox"/> 17.5 | <input type="checkbox"/> 18 |
| | | <input type="checkbox"/> 18 | <input type="checkbox"/> 18.5 |
| | | <input type="checkbox"/> 18.5 | <input type="checkbox"/> 19 |
| | | <input type="checkbox"/> 19 | |

¹ Available only with Low Profile tires

Seat Height, Casters and Rear Wheels continued

Caster

| | | Caster Size | | | |
|-------------|-----------|-------------|--------------------------|--------------------------|--------------------------|
| | | Price | 4" | 5" | 6" |
| Caster Type | 1" Poly | Std | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> |
| | 1.5" Poly | \$50 | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> |

Rear Wheel

| | | Rear Wheel Size | | | | |
|-----------------|------------------------|-----------------|--------------------------|--------------------------|--------------------------|--------------------------|
| | | Price | 16" | 20" | 22" | 24" |
| Rear Wheel Type | Mag | Std | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> |
| | MAXX Performance Spoke | N/C | | | | <input type="checkbox"/> |
| | MAXX Performance Spoke | \$180 | | <input type="checkbox"/> | <input type="checkbox"/> | |

Wheel Locks

Wheel Lock

Option Required

Push to Lock _____ N/C

Wheel Lock Options

Extension Handles E0961 _____ \$67

Accessories

Anti-Tip

Rear Anti-Tips E0971 _____ \$150

Omit Rear Anti-Tips _____ N/C

Other Accessories

1.5" Auto Buckle E0978 _____ \$70

Seat Belt Mounting Hardware _____ \$30

O₂ Holder E2208 _____ \$295

Mounted to the left side. Not available with Fold Down Back.

Cane and Crutch Holder _____ \$195

Rear Tire Type

Full Profile

Street Pneumatic E2211, E2212 _____ Std

Full Poly E2218 _____ N/C

Pneumatic with Airless Insert E2211, E2213 _____ \$80

Low Profile

Low Poly E2218 _____ N/C

Available only with 16" Mag Wheels

Handrim

Not Available with 16" Rear Wheels

Aluminum Anodized _____ Std

Plastic Coated _____ \$115

Axle

Quick Release _____ Std

Armrests

Armrest

Selection Required

Height Adjustable T-Arm (8.5" - 13.5") E0973

Desk _____ \$265

Full _____ \$265

Arm Pad

Standard Arm Pad _____ N/C

Frame Color

Gloss Colors

Black _____ N/C

Candy Blue _____ N/C

Candy Red _____ N/C

Matte Colors

Charcoal Gray _____ N/C

Shadow Black _____ N/C