



SAFE TRANSFERS AND MOVEMENT™ PROGRAM

Hill-Rom® Liko® Sling Sizing Guide

Liko Sling	Size	Recommended Weight Range
Original HighBack/ Soft Original HighBack	Medium	88-176 lb
	Large	154-264 lb
	X-Large	242-550 lb
	XX-Large	440-1100 lb
Solo HighBack (Disposable)	Medium	88-176 lb
	Large	154-265 lb
	X-Large	243-551 lb
Universal	Medium	132-198 lb
	Large	176-264 lb
	X-Large	242-660 lb
Universal <u>High Level</u>	Medium	132-198 lb
	Large	176-264 lb
Original	Medium	88-176 lb
	Large	154-264 lb
	X-Large	240-660 lb
Amputee/HighBack Amputee	Medium	88-132 lb
	Large	132-176 lb
	X-Large	154-440 lb
Ultra Sling	Short/Regular/Long	Up to 1100 lb
Comfort Plus	Medium	88-176 lb
	Large	154-243 lb
	X-Large	220-660 lb
Silhouette	Medium	99-165 lb
	Large	143-220 lb
Repo Sheet®/Solo Repo Sheet®	Regular	Up to 440 lb
Repo Sheet®	Ultra	Up to 1100 lb
Solo Repo Sheet®	Short	Up to 440 lb
Solo Repo Sheet®	Generous	Up to 660 lb
Sabina® Lift with SupportVest™ sling	Medium	Vest Width: 35 in, with straps 55 in
	Large	Vest Width: 41 in, with straps 61 in
	X-Large	Vest Width: 47 in, with straps 63 in



Hill-Rom®

Liko® Sling Sizing Guide

Liko Sling	Size	Recommended Weight Range
Sabina® Lift SafetyVest™ sling	Medium	<176 lb Vest Width: 35 in
	Large	>154 lb Vest Width: 41 in
Sabina® Lift ComfortVest™ sling	Small	<130 lb
	Medium	110-175 lb
	Large	155-240 lb
	X-Large	>220 lb
HygieneSling	Medium	88-198 lb
	Large	176-264 lb
	X-Large	220-660 lb
HygieneVest	Medium	154-220 lb
	Large	198-286 lb
	X-Large	286-660 lb
	XX-Large	330-660 lb
MasterVest™ sling	Medium	132-176 lb
	Large	176-220 lb
	X-Large	220-286 lb
	XX-Large	264-660 lb
LiftPants™	Medium	132-176 lb
	Large	176-220 lb
	X-Large	220-286 lb
	XX-Large	264-660 lb
Ultra LiftPants™	XX-Large	220-1100 lb Patient Height: 67-83 in
	XXX-Large	220-1100 lb Patient Height: 67-83 in

Color Guide:

■ Small = Orange;
 ■ Medium = Yellow;
 ■ Large = Blue;
 ■ X-Large = Black;
 XX-Large = White

See appropriate Liko Sling Instruction Guidelines for complete information.