

Invacare® Top End® Little Excelerator™ Handcycle



PROVIDER / CUSTOMER INFORMATION:			
Provider Information:		End User Information:	
Date:	Name:	Height:	
IVC Acct #:	Ship to Address:	Weight (lbs.):	
Contact:	Address Line 2:	Leg Length (in.):	
Phone #:	City:	Hip Width (in.):	
Fax/Email:	State:	Zip:	Phone #:
<input type="radio"/> Ship to Provider	<input type="radio"/> Ship to End User		Email:

ORDER REQUEST TYPE:		
<input type="radio"/> Quote	<input type="radio"/> Order	PO # <input type="text"/>

BASE MODEL:		
<input checked="" type="radio"/> 10008	Top End Little Excelerator Handcycle	\$4,225

Assembly: <i>(see definitions)</i>		
<input checked="" type="radio"/> 10048	Assembled	\$0
Shipping:		
Shipping Fee <i>(Continental US Only - See definitions)</i>		\$150

SEAT WIDTH:		
<input checked="" type="radio"/> 15013	14" Seat Width	\$0

CAMBER ANGLE:		
<input checked="" type="radio"/> 45004	9 Degree Camber	\$0

WHEELBASE: <i>(Front to Rear Axle)</i>		
<input checked="" type="radio"/> 75074	56" Wheelbase	\$0

WHEEL SETS:		
<input checked="" type="radio"/> 50120	20" Front Wheel & 20" Rear Wheels with Quick Release Axles	\$0

TIRES:		
<input checked="" type="radio"/> 50212	20" Treaded Tires	\$0

FOOTREST:		
<input checked="" type="radio"/> 60002	Dual, Height/Angle Adjustable	\$0



DRIVE CONTROL:		
CRANK TYPE:		
<input checked="" type="radio"/> 65074	Crank Arms: 155mm Long by 12" Wide	\$0
HAND PEDAL TYPE: <i>(see definitions)</i>		
<input type="radio"/> 65137	Medium 3D Printed with Power Plates	\$0
<input type="radio"/> 65128	Horizontal Hand Pedals	\$127
<input type="radio"/> 65134	Tri-Pin Hand Pedals	\$361
<input type="radio"/> 65135	Quad Cuff Hand Pedals	\$481
<input type="radio"/> 65133	Quad-Grips Hand Pedals with Wrist Wraps & 3" Crank Width Extensions	\$962
SHIFTING:		
<input type="radio"/> 65060	Twist Activated Shifter with Internal Gearing	\$0
<input type="radio"/> 65061	Quad Twist Shifting Adaptor <i>(see definitions)</i>	\$181
BRAKING:		
<input checked="" type="radio"/> 65080	Reverse Activated Drum Brake	\$0
<input checked="" type="radio"/> 65090	Locking Parking Brake	\$0
CHAINRING:		
<input checked="" type="radio"/> 65100	Single 38 Tooth Chainring	\$0

FRAME COLOR:				
	55018	Black	Wet Black	\$0
	55019	Black	Black Prism	\$0
	55017	Black	Stealth Black	\$0
	55010	Blue	Midnight Blue	\$0
	55011	Blue	Royal Blue	\$0
	55012	Blue	Electric Blueberry	\$0
	55023	Grey	Stealth Grey	\$0
	55053	Green	Agriculture Green	\$0
	55008	Green	Grasshopper Green	\$0
	55054	Green	Neon Green	\$0
	55006	Orange	Sunkissed Orange	\$0
	55005	Orange	Simple Orange	\$0
	55051	Orange	Neon Orange	\$0
	55015	Pink	Bubblegum Pink	\$0
	55057	Pink	Neon Pink	\$0
	55055	Purple	Lilac Purple	\$0
	55001	Red	Simple Red	\$0
	55002	Red	Tangerine Red	\$0
	55003	Red	Raspberry Red	\$0
	55050	Red	Neon Red	\$0
	55020	Silver	Silver Metallic	\$0
	55016	White	White Snow	\$0
	55007	Yellow	Safety Yellow	\$0
	55052	Yellow	Neon Yellow	\$0

Invacare® Top End® Little Excelerator™ Handcycle



HANDCYCLE SPECIFICATIONS:	
Weight:	45 pounds
Turning Radius:	10 feet for 360 degree circle, 7 feet on a 3-point turn
Transfer Height:	17.5" from ground with 2" cushion
Frame Style:	Single main tube, upright position, 1-3/4" X .065" ovalized carbon steel tubing
Seat:	Seat adjust forward and rearward with sliding mechanism
Fork:	Carbon Steel
Propulsion / Shifting:	Internal front hub, twist activated shifter mounted on crank support
Cranks:	155mm adjustable height
Back:	Adjustable height
Footrest:	Dual adjustable height and angle, includes straps
Hand Pedals:	Ovalized ergonomic with frictionless bearings
Brake:	Reverse activated drum and Locking Parking Brake mounted on fork
Wheelbase:	56" axle to axle
Camber:	9 degree
Wheels:	20" metal spoked wheels
Axles:	Quick Release
Tires / Tubes:	20" X 1.75" tires (40 PSI)
Positioning Straps:	Hook and loop strap included
Cushions:	Seat cushion included
Safety / Visibility:	Flag and wheel reflectors included
Chain Guard:	Full length, shock resistant composite material

DEFINITIONS:	
Assembly and Shipping:	Arrives factory assembled, rear wheels and some accessories will need to be installed. Shipping rate is \$150 for all shipments in the continental US. Freight quotes for shipments outside the continental US are available upon request.
Quad Twist Shifter Adaptor:	Three pronged adaptor which allows user to tap for shifting.
Tri-pin Hand Pedals:	3-Point contact hand pedal similar to hand controls used for vehicles. Allows for the hands to be kept in place for people with limited grip strength.
Quad Grips Hand Pedals:	Combination of Tri-Pin hand pedal and Quad Cuff hand pedal. Allows for the hands to be kept in place for people with limited grip strength.
Quad Cuff Hand Pedals:	Includes Power Plates on the top and the bottom of the pedal to assist through both up and down stroke. Includes an adjustable metal cuff with a hook and look strap to secure hands.

Warranty Info:

Invacare warrants the frames when purchased new and unused to be free from defects in materials and workmanship for a period of three (3) years from the date of purchase from Invacare or a dealer, with a copy of the seller's invoice required for coverage under this warranty. Invacare warrants the upholstered materials (seat and back) and remaining components of this product when purchased new and unused to be free from defects in materials and workmanship for a period of thirteen (13) months from date of purchase from Invacare or a dealer, with a copy of the seller's invoice required for coverage under this warranty