



## **EW-M41**

Features include a maximum speed of 5mph, a weight capacity of 350 lbs., a range of 15 miles, a 22 AH battery pack for longer travel range, LED headlight, easy take apart mechanism, on board battery charger, and delta handlebar.

### EW-M41 Specifications:

Overall Length  
44"

Overall Width  
23"

Overall Height  
37"

Wheels:  
Front 9"  
Rear 9"

Weight w/ Batteries  
121 lbs.

Max Speed  
5mph

Weight Capacity  
350 lbs.

Ground Clearance  
2"  
Grade Climbable  
8 degree

Turning Radius  
55"  
Brake  
Electromagnetic

Seat Type  
Swivel Light Weight Foldable

Seat Width  
18"

Motor Size  
250 W

Battery Size  
(2)12V 22AH

Battery Weight  
34lbs.

Travel Range  
15 Miles

Battery Charger  
On Board

Electronics  
On/Off Key Switch, Battery Indicator, Speed Control Knob