



EW-M43 Folding Power Wheelchair

Key Specs:

Top Speed: 3.7 mph

Turning Radius: 24.25"

Battery Charge Distance: 9.7 miles

Weight Capacity: 300 lbs.

Flat Free Tires

Folding Armrest

Large Aluminum Footplate

General Specifications:

Swing-Away joystick: Yes

Folded Size: 31"H x 15.5"D x 23.5"W

Lap Belt: Yes

Back Seat Pocket: Yes

Front Suspension: Yes

Under Seat Removable Storage: Yes

Optional Left or Right Hand Joystick: Yes

Frame Material: Aluminum

Top Speed: 3.7 mph

Heaviest Piece: 56 lbs.

Battery Charge Distance: 9.7 mph

Drive Type: Rear Wheel Drive

Foldable: Yes

Footplate: Large Aluminum

Seat Width: 17.25"

Seat Depth: 17.5"

Seat to Floor Height: 22"

Overall Width: 24.25"

Overall Length: 42"

Overall Height: 39"

Joystick Options: Right or Left

Frame Type: Aluminum Transportable

Battery: 24V 10AH Lithium Battery

Battery Included: Yes

Battery Charger Type: Off Board

Rear Wheel Size: 12"

Front Wheel Size: 8"

Flat Free Tires: Yes

Base Product Weight: 60 lbs.

Battery Pack Weight: 4 lbs.

Weight Without Batteries: 56 lbs.

Adjustable Anti-Tip Wheels