

EW-M48 Power Wheelchair

Key Specs:

Seat Width: 20" deluxe captain seat

Ultra-Comfort Front, Rear & Center

Suspension

Top Speed: 4.2 mph

Turning Radius: 24"

Weight Capacity: 300 lbs.

General Specifications:

Swing-Away joystick: Yes

Lap Belt: Yes

Turning Radius: 24"

Optional Left or Right Hand Joystick: Yes

Brake: Electro-Mechanical

Frame Material: Aluminum

Top Speed: 4.2 mph

Adjustable Variable Speed: Yes

Self-Diagnostic Warning Lights: Yes

Battery Charge Distance: 12.4 miles

Drive Type: 2 Motor Mid-Wheel Drive

Grade Climbable: 6 degree

Foldable: Yes

Footplate: Large flip-up Footplate

Adjustable Footplate: Yes, 3-Angle
Adjustment

Ground Clearance: 2.4"

Seat Type: 20" Swivel Captain Seat with
Adjustable Headrest & Seat Height

Adjustable Armrest: Yes, Width and
Height

Removable Seat: Yes

Seat Width: 20"

Seat Depth: 19"

Seat to Floor Height: 21"

Overall Width: 24"

Overall Length: 39"

Overall Height: 42.5"

Electric Horn: Yes

Battery Gauge: Yes

Free Wheel Levers: Yes

Joystick Options: Right or Left

Extendable Joystick Adjustment: Yes

Frame Type: Aluminum

Battery: (2) 12V 36AH
Battery Included: Yes
Battery Charger Type: Off Board 5A
Middle Wheels: 10"
Rear Wheels Size: 6"
Front Wheels Size: 6"
Solid Flat Free Wheels: Yes
Base Product Weight: 181 lbs.
Battery Pack Weight: 49.6 lbs.
Weight Without Batteries: 131.4 lbs.
Anti-Tip Wheels: Yes
Rear & Side Reflectors: Yes
Weight Capacity: 300 lbs.