

EVALUATION FORM **MANTIS FRAME RUNNER**

The **MANTIS** by Freedom Concepts gives you the chance to walk, run, stride and glide wherever you go. This mobility device provides support in a non-restrictive way, allowing the user to take the lead and have the freedom needed to feel independent. Its design allows for indoor and outdoor use, so whether it's in the community or at school, you can keep up with the action.

Three-Step Order Process:

1. Select your size.
2. Select accessories.
3. Call our team if you have any questions.



 **Call 1-800-661-9915 or visit**
FREEDOMCONCEPTS.COM

User Information:

User Name: _____

Date of birth (YYYY-MM-DD): _____

Diagnosis: _____

Inseam: Inches: _____ Centimeters: _____

Height: Inches: _____ Centimeters: _____

Weight: LBS: _____ Kilograms: _____

If our current options do not meet your individual needs please contact us about custom work. Our team of engineers are committed to helping everyone find success with our products.

SIZING AND DIMENSIONS: Recommended User Size Ranges / Frame Dimensions

Mantis Sizing	Weight capacity	Recommended Inseam range	Device Width	Device Length	Min. seat height	Max. seat height	Product weight
Small	140 lbs	18"-24"	29.5"	60"	16"	24.5"	39 lbs
Medium	175 lbs	22"-29"	33.5"	64.5"	20.5"	29"	45 lbs



Sizing and dimensions can be reduced if needed. **The MANTIS** comes standard with quick-release, removable wheels for increased portability and easier storage.

MODEL SIZE - see front page chart for sizing guide

MANTIS - Size Small



SMALL (MANTISSM)

OR

MANTIS - Size Medium



MEDIUM (MANTISMED)

FRAME COLOR - select one

Standard Color (STD)



- | | |
|---|--|
| ☐ Electric Green
(GREEN) | ☐ Candy Blue
(CANDY BLUE) |
| ☐ Cherry Red
(RED) | ☐ Hot Rod Purple
(HOT ROD PURPLE) |

Premium Color



- | | |
|-------------------------|------------------------------|
| Snow (WHITE) | Arctic Teal (TEAL) |
| Fossil (FOSSIL) | Blaze Orange (ORANGE) |
| Midnight (BLACK) | Radiant Pink (PINK) |
| Chrome (CHROME) | |

Custom Color (CUSTOM COLOR)

Indicate custom color note here:

DECAL COLOR - select one



BLACK (6401802) (STD)



WHITE (6401803) (STD)

NAMEPLATE - optional



CUSTOM NAMEPLATE (NAMEPLATE CO)

FRONT BRAKE LOCATION select one

- | | |
|----------------------------------|--------------|
| Left Handlebar (2113219) | (STD) |
| Right Handlebar (2113220) | (STD) |

TIRE TYPE - select one

SMOOTH SOLID TIRES (STD) Small: (2112811) Medium: (2112814) <i>Ideal for indoor and outdoor use, with limited speeds and/or with caregiver present</i>	 PNEUMATIC Small: (2112809) Medium: (2112812) <i>Ideal for more independent riders, wanting to ride outdoors or on uneven surfaces.</i>
	

SEAT - select one

	
UNICYCLE STYLE SEAT (STD) (2114800)	MOUNTAIN BIKE STYLE SEAT (2114801)

CHEST PAD FABRIC - select one

	
WIPEABLE VINYL FABRIC (STD) (2115800)	BREATHABLE REMOVABLE FABRIC (2115801)

MANTIS™ order notes

Customer Name: _____

Email: _____

Phone Number: _____

Shipping Address: _____

City: _____ Prov/State: _____

Postal/Zip Code: _____ Federal ID/SS#: _____

Credit Card Number: _____

Name On Card: _____

Exp Date: _____ / _____ CVV / CVC: _____

Sub-Total	
Tax	
Freight	
TOTAL	

Terms & Conditions Prices and frame colors are subject to change without notice. Prepaid or net 30 days O.A.C. All taxes extra. Freight extra. Manufacturer reserves the right to change prices, colors, features, dates, locations and products without notice. Residential delivery charge may apply. Duty and brokerage included where applicable.

Privacy Policy Freedom Concepts Inc. collects, uses and discloses personal information of the Customer and/or Rider to properly assemble and manufacture the mobility device and may contact the Rider's therapist for additional personal information to assist with this process. The collection, use and disclosure of the personal information of the Customer and/or the Rider is governed by any consent which may, from time to time, be provided, and by the Freedom Concepts Inc. Privacy Policy which is available on our website www.freedomconcepts.com or by contacting our office. By disclosing your information here, you hereby warrant that you have the authority to bind the User and hereby authorize the collection, use, and disclosure of personal information as herein provided and in accordance with the Freedom Concepts Inc. Privacy Policy.



INCLUDED STANDARD FEATURES



1. Removable Wheels:

Tool-free, unidirectional removable wheels with marks for easy placement, improving portability and storage.

4. Hoop Handlebars:

These rounded handlebars allow for multiple different grips and holding styles to accommodate mobility needs. The foam grip provides a comfortable and easier steering experience.

2. Dropper Seat:

The selected seat is mounted on a dropper seat post that is telescopic and can be lowered by a lever throttle for transfers on, and released to raise up to support the user. The quick-release lever makes seat height adjustments easy, while the dropper post can be raised or lowered hands-free with a press of the remote lever.

5. Steering Damper:

Controls and absorbs impact to steering on uneven surfaces and helps to bring the handlebars back to centre.

3. Adjustable Chest Pad:

The post and angle of the chest pad are adjustable to help the user get into a supportive, comfortable position. Additionally, it can be locked at fixed angles or set to move dynamically while in motion.

6. Abductor Pad:

This padding wraps around the lower frame posts to protect knees and legs during use if adduction occurs.