

# The Cozzia MC-520

---



MC-520 CRA




MC-520 BRI



MC-520 SAD



**Cranberry**



**Saddle**



**Briarwood**

The MC-520 is the only lift chair in the world that helps you get out, *but also helps you get in.*

The MC-520 features a lift function, retractable arms, an adjustable headrest, air massage, and heat.

**cozzia®**

# The Cozzia MC-520

Model #	MC-520
Chair Sitting Dimensions (LxWxH)	35.9" x 30" x 44"
Chair Upright/Lifted Dimensions (LxWxH)	37.8" x 30" x 59.5"
Chair Full Recline Dimensions (LxWxH)	73.3" x 30" x 22.9"
Shipping Dimensions (LxWxH)	54.3" x 30.3" x 33.7"
Weight Capacity	300 LBS
Headrest Adjustable Height	3"
Seat Width Between Arms	21"
Seat Height	20"
Seat Depth	20"
Back Height from Seat	28"
Distance from Wall	15"
Recline Position #	Infinite Recline Positions
Fabric Type	Bonded Fabric
Trendelenburg Position	✓
Memory Foam Seating	✓
Heat Therapy	✓
Air Massage	✓
Air Lumbar Support	✓
Retractable Arm	✓
Adjustable Headrest	✓
Footrest Extension	✓
Casters for Easy Moving	✓
Battery Backup	✓

## Easy in. Easy out.

The MC-520 has an electric lift function to give you easy access to the chair. Both armrests of the MC-520 are power adjustable to make getting in to and out of the chair seamless.



## Relax with heat and massage.

The MC-520 has many features built into it that make it our most advanced lift chair. These include: Stretch air massage, Relax air massage, and heating.



## Adjustable headrest.

The headrest of the MC-520 is adjustable in two ways. You can alter the headrest's height, and you're able to move the wings of the headrest in and out.



WWW.COZZIAUSA.COM

**cozzia**®

TEL: 887-977-0656 / FAX: 800-521-4712

861 S. OAK PARK ROAD, COVINA, CA 91724